

## Product datasheet for **SC307951**

### SLC22A25 (NM\_199352) Human Untagged Clone

#### Product data:

|                           |                        |
|---------------------------|------------------------|
| Product Type:             | Expression Plasmids    |
| Tag:                      | Tag Free               |
| Symbol:                   | SLC22A25               |
| Synonyms:                 | HIMTP; UST6            |
| Mammalian Cell Selection: | Neomycin               |
| Vector:                   | pCMV6-Entry (PS100001) |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_199352, the custom clone sequence may differ by one or more nucleotides

```
ATGGCCTTTCAGGACCTCCTAGATCAAGTTGGAGGCCTGGGGAGATTCCAGATCCTTCAGATGGTTTTC
TTATAATGTTCAACGTCATAGTATACCATCAAACCTCAGCTGGAGAATTTCGAGCATTTCATCTGATCA
TCGCTGCTGGGTTTATATACTGGACAATGACACTATCCCTGACAATGACCCTGGGACCCTCAGCCAGGAT
GCCCTCCTGAGAATCTCCATCCCATTCGACTCAAATCTGAGGCCAGAGAAGTGTGTCGCTTGTCCATC
CCCAGTGGAAGCTCATTTCATCTGAATGGGACCTTCCCCAACACGAGTGAGCCAGATACAGAGCCCTGTGT
GGATGGCTGGGTATATGACCAAAGCTCCTTCCCTTCCACCATTGTGACTAAGTGGGATCTGGTATGCGAA
TCTCAACCACTGAATTCAGTAGCTAAATTTCTATTTCATGGCTGGAATGATGGTGGGAGGCAACCTATATG
GCCATTTGTCAGACAGGTTTGGGAGAAAGTTTCGTGCTCAGATGGTCTTACCTCCAGCTCGCCATTGTAGG
CACCTGTGCGGCCTTTGCTCCACCATCCTCGTATACTGCTCCCTGCGCTTCTTGGCTGGGGCTGCTACA
TTTAGCATCATTGAAATACTGTTTGTAAATTTAGAGTGGATAACTCACCATTCTGTGCCATGGCAT
TGACATTGACACTTTGTGCTGCTAGTATTGGACATATAACCCTGGGAAGCCTGGCTTTTGTATTGAGCA
CCAGTGCATCCTCCAGTTGGTGATGTCTGCACCATGCTTTGTCTTTCTGTTCTCAAGGTGGCTGGCA
GAGTCTGCTCGGTGGCTCATTATCAACAACAAACCAGAAAGAGGGCTTAAAGGAACCTAGAAAAGCTGCAC
ACAGGAATGGAATGAAGAATGCTGAAGACATCCTAACCATGGAGGTTTTGAAATCCACCATGAAGCAAGA
ACTGGAGGCAGCAGAAAAAGCATTCTCTTTGTAATTGCTCCGCATACCCACATATGTAAGAAGATC
TGTTCCTGTCTTTGTGAGATTTGCAAGTACCATCCCTTTTGGGGCCTTACTTTGCACCTCCAGCATC
TGGGAAACAATGTTTCTGTTGCACTCTCTTGGTGCAGTCACCTCCTGGCCAATTGTGTTGCACC
TTGGGCACTGAATCAGATGAGCCGTCGACTAAGCCAGATGCTTCTCATGTTCTACTGGCAACCTGCCTT
CTGGCCATCATATTTGTCCTCAAGAAATGCAGACCCTGCGTGTGGTTTTGGCAACCCTGGGTGTGGAG
CTGCTTCTCTTGGCATTACCTGTTCTACTGCCCAAGAAAATGAACTAATTCCTTCCATAATCAGGGGAAG
AGCTACTGGAATCACTGGAACTTTGCTAATATTGGGGGAGCCCTGGCTTCCCTCATGATGATCCTAAGC
ATATATTCTCGACCCTGCCCTGGATCATCTATGGAGTCTTTGCCATCCTCTGTCCTTGTGCTCTCC
TCCTTCTGAAACCAGGAACCGCTCTTCTTGACAGCATCCAGGATGTGGAATGAGGGAGTAATAG
CCTAGCTGCCCTCAGAGGAGCTCTGTGCTATAG
```

**Restriction Sites:** Please inquire

**ACCN:** NM\_199352

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <a href="#">NM_199352.2</a> , <a href="#">NP_955384.2</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>RefSeq Size:</b>           | 1644 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 1644 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 387601                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>            | <a href="#">Q6T423</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 11q12.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Families:</b>      | Transmembrane                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |