

Product datasheet for **SC307957**

TPD52L2 (NM_199360) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TPD52L2 (NM_199360) Human Untagged Clone
Tag:	Tag Free
Symbol:	TPD52L2
Synonyms:	D54; TPD54
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_199360, the custom clone sequence may differ by one or more nucleotides ATGGACTCCGCCGCAAGATATCAACCTGAATTCTCCTAACAAAGGTCTGCTGTCTGAC TCCATGACGGATGTTCTGTGACACAGGTGTGGCTGCCCGACTCCTGTGTTGAGGGT CTGACAGAGGCTGAGGAGGAGGCTCAGGGCTGAGCTTACCAAGGTGGAAGAGGAAATT GTCACTCTGCGCCAGGTCTGGCAGCCAAGGAGAGGCACTGTGGAGAGCTCAAGAGGAGG CTGGCCCTCTCCACCCTGGGGGAGCTGAAACAGAACCTGTCCAGGAGCTGGCATGACGTG CAGGTCTCTAGCGCCTATGTGAAAACCTTCTGAGAACTTGGAGAGTGGAATGAGAAAGT ACCCAGTCAGACCTCTACAAGAAGACTCAGGAACTCTTTCACAGGCAGGACAGAAGACT TCAGCTGCCCTGTCCACAGTGGGCTCTGCCATCAGCAGGAAGCTTGGAGACATGAGGGCT CATCCATTTTCACACTCCTTTAGCAGCTACTCCATCCGCACTCAATAAGTATGCCAGCC ATGAGGAACTCTGCGACCTTCAAGTCGTTTGGAGACCGAGTTGGGACCATAAAGTCTAAG GTTGTGGGTGACAGAGAGAACGGCAGTGACAACCTCCCTTCTCAGCGGGGAGTGGTGAC AAGCCCTGTCCGATCCCGCACCTTTCTAA
Restriction Sites:	Please inquire
ACCN:	NM_199360
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199360.1](#), [NP_955392.1](#)

RefSeq Size: 2426 bp

RefSeq ORF: 690 bp

Locus ID: 7165

UniProt ID: [O43399](#)

Cytogenetics: 20q13.33

Gene Summary: This gene encodes a member of the tumor protein D52-like family. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. Expression of this gene may be a marker for breast cancer and acute lymphoblastic leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 12. [provided by RefSeq, Aug 2011]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.