

Product datasheet for **SC308048**

Galectin 8 (LGALS8) (NM_201544) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Galectin 8 (LGALS8) (NM_201544) Human Untagged Clone
Tag:	Tag Free
Symbol:	Galectin 8
Synonyms:	Gal-8; PCTA-1; PCTA1; Po66-CBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC308048 representing NM_201544. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATGTTGTCCTTAAACAACCTACAGAATATCATCTATAACCCGGTAATCCCGTTTGTGGCACCATT
CCTGATCAGCTGGATCCTGGAACTTTGATTGTGATACGTGGGCATGTTCCCTAGTGACGCAGACAGATTC
CAGGTGGATCTGCAGAATGGCAGCAGCATGAAACCTCGAGCCGATGTGGCCTTTTCATTTCAATCCTCGT
TTCAAAGGGCCGGCTGCATTGTTTGAATACTTTGATAAATGAAAAATGGGACGGGAAGAGATCACC
TATGACACGCCTTTCAAAGAGAAAAGTCTTTGAGATCGTGATTATGGTGCTGAAGGACAAATTCAG
GTGGCTGTAATGAAAACATACTCTGCTCTATGGCCACAGGATCGGCCAGAGAAAATAGACACTCTG
GGCATTATGGCAAAGTGAATATCACTCAATTGGTTTTAGCTTCAGCTCGGACTTACAAAGTACCCAA
GCATCTAGTCTGGAAGTACAGAGATAAGTAGAGAAAATGTTCCAAAGTCTGGCACGCCCCAGCTTAGG
CTGCCATTCGCTGCAAGGTTGAACACCCCATGGGCCCTGGACGAACTGTGCTCGTTAAAGGAGAAGTG
AATGCAAATGCCAAAAGCTTTAATGTTGACCTACTAGCAGGAAAATCAAAGGATATTGCTCTACACTTG
AACCCACGCCTGAATATTAAGCATTGTAAGAAATCTTTTCTCAGGAGTCTGGGGAGAAGAAGAG
AGAAATATTACCTCTTCCCATTTAGTCTGGGATGTACTTTGAGATGATAATTTATTGTGATGTTAGA
GAATTCAGGTTGCAGTAAATGGCGTACACAGCCTGGAGTACAAACACAGATTTAAAGAGCTCAGCAGT
ATTGACACGCTGGAAATTAATGGAGACATCCACTTACTGGAAGTAAGGAGCTGGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_201544
Insert Size:	954 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_201544.2</u>
RefSeq Size:	6212 bp
RefSeq ORF:	954 bp
Locus ID:	3964
UniProt ID:	<u>O00214</u>
Cytogenetics:	1q43
MW:	35.8 kDa

Gene Summary:

This gene encodes a member of the galectin family. Galectins are beta-galactoside-binding animal lectins with conserved carbohydrate recognition domains. The galectins have been implicated in many essential functions including development, differentiation, cell-cell adhesion, cell-matrix interaction, growth regulation, apoptosis, and RNA splicing. This gene is widely expressed in tumoral tissues and seems to be involved in integrin-like cell interactions. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (3) differs in the 5' UTR and lacks an in-frame segment of the coding region, compared to variant 1. Variants 2 and 3 encode a shorter isoform (b) that is missing an internal segment compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments. CCDS Note: This CCDS representation uses the 5'-most in-frame start codon, which is conserved in higher primates. An alternative downstream start codon, which is more widely conserved, also exists. It is possible that leaky scanning by ribosomes would allow the downstream start codon to be used, at least some of the time. The use of the downstream start codon would result in a protein that is 1 aa shorter at the N-terminus. There is no experimental evidence showing which start codon is preferentially used in vivo.