

Product datasheet for **SC308070**

CACNB2 (NM_201597) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CACNB2 (NM_201597) Human Untagged Clone
Tag:	Tag Free
Symbol:	CACNB2
Synonyms:	CAB2; CACNLB2; CAVB2; MYSB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC308070 representing NM_201597.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGTCCAAGGGACATGTCCAAGTCGCCTCCCACAGCGGGCGGCGGTGGCGCAGGAGATCCAGATG
GAACCTGTAGAGAACGTGGCTCCCGGGGGCGCTCGGAGCCCGCCGACAGTCATATGAAAAGGAGCC
AGAAGGAAAAACAGATTTAAAGGATCTGATGGAAGCACGTCATCTGATACTACCTCAAATAGTTTGT
CGCCAGGGTTCGGCAGACTCCTACACTAGCCGTCCATCCGATTCCGATGTATCTCTGGAGGAGGACCGG
GAGGCAGTGCGCAGAGAAGCGGAGCGGCAGGCCAGGCACAGTTGAAAAAGCAAAGACAAAGCCCGTT
GCATTTGCGGTTTCGGACAAATGTCAGCTACAGTGCGGCCATGAAGATGATGTTCCAGTGCCTGGCATG
GCCATCTCATTGAAGCAAAGATTTTCTGCATGTTAAGGAAAAATTAACAATGACTGGTGGATAGGG
CGATTGGTAAAAGAAGGCTGTGAAATCGGATTCATTCCAAGCCAGTCAAAGTAAAACATGAGGCTG
CAGCATGAACAGAGACCAAGCAAGGAAATTCTACTCCAGTAAATCAGGAGGAAATTCATCATCCAGT
TTGGGTGACATAGTACCTAGTTCAGAAAAATCAACACCTCCATCATCTGGTGAAAAATCTGCAGATGAA
CAAGACCAGTGGAAAACAGGCTTGTTCGCGGTTTACTACAGAGCACACTCCTCCGTATGATGTG
GTACCTTCCATGCGACCAGTGGTCTAGTGGGCCCTTCTCTGAAGGGCTACGAGGTCACAGATATGATG
CAAAAAGCGCTGTTTGAATTTTTAAAACACAGATTTGAAGGGCGGATATCCATCACAAGGGTACCAGCT
GACATCTCGCTTGCCAAACGCTCGGATTAACAATCCCAGTAAGCACGCAATAATAGAAAGATCCAAC
ACAAGGTCAAGCTTAGCGGAAGTTCAGAGTGAATCGAAAGGATTTTTGAACTTGCAAGAACATTGCAG
TTGGTGGTCCCTGACGCGGATACAATTAATCATCCAGCTCAACTCAGTAAAACCTCCTGGCCCCATT
ATAGTATATGTAAGATTTCTTCTCCTAAGTTTTACAAAGTTAATAAAATCTCGAGGGAAATCTCAA
GCTAAACACCTCAACGTCCAGATGGTAGCAGTGATAAACTGGCTCAGTGCCTCCAGAGTGTTCGAT
GTGATCTTGGATGAGAACCAGCTTGAGGATGCCTGTGAGCACCTTGCCGACTATCTGGAGGCTACTGG
AAGGCCACCCATCTCCAGCAGTAGCCTCCCAACCCTCCTTAGCCGTACATTAGCCACTTCAAGT
CTGCCTCTTAGCCCCACCTAGCCTCTAATTCACAGGGTCTCAAGGTGATCAGAGGACTGATCGCTCC
GCTCCTATCCGTTCTGCTTCCCAAGCTGAAGAAGAACCCTAGTGTGAACCAGTCAAGAAATCCCAGCAC
CGCTTCTCCTCCTCAGCCCCACACCACAACCATCGCAGTGGGACAAGTCCGCGCTCCTCCAGGCAAGAG
ACATTTGACTCGGAAACCAGGAGAGTCGAGACTCTGCCTACGTAGAGCCAAAGGAAGATTATCCCAT
GACCACGTGGACCACTATGCCTCACACCGTGACCACAACCACAGAGACGAGACCCACGGGAGCAGTGAC
CACAGACACAGGGAGTCCCGCACCGTTCGCGGACGTGGATCGAGAGCAGGACCAACGAGTGAAC
AAGCAGCGCAGCCGTATAAATCCAAGGATCGCTACTGTGAAAAGGATGGAGAAGTGATATCAAAAAA
CGGAATGAGGCTGGGGAGTGAACAGGGATGTTTACATCCGCCAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:

ACCN: NM_201597

Insert Size: 1911 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_201597.2
RefSeq Size:	4006 bp
RefSeq ORF:	1911 bp
Locus ID:	783
UniProt ID:	Q08289
Cytogenetics:	10p12.33-p12.31
Protein Families:	Druggable Genome, Ion Channels: Other
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM), MAPK signaling pathway
MW:	71.1 kDa
Gene Summary:	<p>This gene encodes a subunit of a voltage-dependent calcium channel protein that is a member of the voltage-gated calcium channel superfamily. The gene product was originally identified as an antigen target in Lambert-Eaton myasthenic syndrome, an autoimmune disorder. Mutations in this gene are associated with Brugada syndrome. Alternatively spliced variants encoding different isoforms have been described. [provided by RefSeq, Feb 2013]</p> <p>Transcript Variant: This variant (4) lacks an alternate in-frame exon and contains an alternate in-frame exon, in the 5' coding region, compared to variant 2. The encoded isoform (4) is shorter than isoform 2.</p>