

Product datasheet for **SC308143**

TMEM81 (NM_203376) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM81 (NM_203376) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM81
Synonyms:	HC3107; KVL2788; UNQ2788
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC308143 representing NM_203376. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**GCGATCGCC**
 ATGAAGGTTTTAGCCACTAGTTTTGTCCTTGGGAGCCTGGGGTTGGCCTTCTACCTGCCTTTGGTGGTG
 ACTACACCTAAACACTGGCCATCCCTGAGAAGCTGCAAGAAGCTGTGGGAAAGTTATCATCAATGCC
 ACAACCTGTACTGTCACCTGTGGCCTTGCTATAAGGAGGAGACCGTCTGTGAGGTGGGCCCTGATGGA
 GTGAGAAGGAAATGTCAGACTCGGCGCTTAGAATGTCTGACCAACTGGATCTGTGGGATGCTCCATTTT
 ACCATTCTCATTGGCAAGGAATTTGAGCTTAGCTGTCTGAGTTCAGACATCTGGAGTTTGGACAGGAA
 GCTTTCCGGTTCACCTGGAGACTTGCTCGAGGTGTCATCTCCACTGACGATGAGGTCTTCAAACCTTT
 CAAGCCAACTCCCACTTTGTGAAGTTTAAATATGCTCAGGAGTATGACTCTGGGACATATCGCTGTGAT
 GTGCAGCTGGTAAAAAATTGAGACTCGTCAAGAGGCTCTATTTTGGGTTGAGGGTCCTTCTCCTAAC
 TTGGTGAATCTGAATTTCCATCAGTCACCTACTGAGGATCAGAAGTTAATAGATGAGGGATTGGAAGTT
 AATCTGGACAGCTACTCCAAGCCTCACCACCAAGTGGAAAAAGAAGGTGGCGTCAGCCTTGGGAATA
 GGAATTGCCATTGGAGTGGTGGTGGCGTGTGGTGGAGATTGTCCTCTGTGCGCTAAGGGGGGGCCTG
 CAGCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_203376
Insert Size:	768 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_203376.1</u>
RefSeq Size:	1350 bp
RefSeq ORF:	768 bp
Locus ID:	388730
UniProt ID:	<u>Q6P7N7</u>
Cytogenetics:	1q32.1
Protein Families:	Transmembrane
MW:	28.5 kDa