

Product datasheet for **SC308176**

ASCL4 (NM_203436) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ASCL4 (NM_203436) Human Untagged Clone
Tag:	Tag Free
Symbol:	ASCL4
Synonyms:	ASH-4; bHLHa44; HASH4
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	<p>>OriGene sequence for NM_203436 edited</p> <p>TGCACCAGCACCTCTGCTTCTAAGATTCTGTTTCGTCTTCTTCTATTGAGAGATTGACCT CTTGAGTGATTTGTGTGCTTTCCGGCAATTGATGGAGACGCGTAAACCGCGGAACGGCT GGCTTGGCATACTCGCTGCGCACCAGCGCCCTGGGCGTTCCGGGGACCTGCCCGACT CCCGCGGAGGGACCCCTCAGGGTCGCCCTGCGTCTGGACGCCGCGTGCTGGGAGTGGGC GCGCAGCGGCTGCGCACGGGGATGGCAGTACTTGCCCGTGCCGCTGGACAGCGCCTTCGA GCCCGCTTCTCCGCAAGCGCAACGAGCGCGAGCGGCGAGCGGGTGCGTGCCTGAACGA GGGCTATGCGCGCTCCGAGACCACCTGCCCCGGGAGCTGGCAGACAAGCGCCTCAGCAA AGTGGAGACGCTCCGCGCTGCCATCGACTACATCAAGCACCTGCAGGAGCTGCTGGAGCG CCAGGCCTGGGGGCTCGAGGGCGCGCCGGCGCCGTCCCCAGCGCAGGGCGGAATGCAA CAGCGACGGGGAGTCCAAGGCCTCTTCGCGCCTTCGCCAGCAGCGAGCCCGAGGAGGG GGGCAGCTAGCGAGCGCCCGAACTGGCCAGGACCCCGCGCCCGCACAGCGCGCAGC CGGGCGCTCAACCTAAGTCTCTTCTGAAGGTGGTTTGATTCTTAATCTGGCATCTTCT CCAGGCCG</p>
Restriction Sites:	Please inquire
ACCN:	NM_203436
Insert Size:	700 bp


[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_203436.1 , NP_982260.1
RefSeq Size:	2254 bp
RefSeq ORF:	519 bp
Locus ID:	121549
UniProt ID:	Q6XD76
Cytogenetics:	12q23.3
Gene Summary:	Basic helix-loop-helix transcription factors, such as ASCL4, are essential for the determination of cell fate and the development and differentiation of numerous tissues (Jonsson et al., 2004 [PubMed 15475265]).[supplied by OMIM, Mar 2008]