

Product datasheet for **SC308252**

CAGE1 (NM_205864) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CAGE1 (NM_205864) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAGE1
Synonyms:	bA69L16.7; CT3; CT95; CTAG3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC308252 representing NM_205864.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGACAGAGAAGCCAGAATTTCAAAGTCAAGTGTATAATTATGCAAAAGACAACAATATAAGCAAGAC
TCATTTAAGGAAGAAAATCCAATGGAACTAGTGTCTTGCAAAATACAGACCAACTTGGTAACGAGTAT
TTTAGACAGCCTCCTCCTAGAAGCCCGCTCTAATCCACTGCAGTGGAGAGATGCTGAAATTTACAGAA
AAGTCACTGGCTAAAAGTATTGCCAAGGAATCTGCCCTCAACCTAGCCAACCTCCAAGCTTCTTATGT
AAGACTGCAGTTCCTTCAAAGAAATACAGAATTATGGGAGATTCTGAGATGTCAGTCAGTTATGAA
AAGGAAGTCACAGCAGAGGGTGTGGAGAGGCCAGAAATGTCTCAACTTGGTCTTCGGCAGGCATTTCC
TGGAGGAGTGAAGCATGTCGGGAGAACTGTGAGATGCCTGACTGGGAGCAAAGTCTGAAAGCTTACAA
CCTGTTCAAGAGGACATGGCTTTAAATGAAGTCTTGCAAGAAATTAACATACTAACAGAAAGCAGGAA
GTGCGAATCCAAGAACTCCAGTGTAGTAACCTGTATTTAGAGAAGAGGGTTAAAGAACTACAGATGAAG
ATTACCAACAGCAAGTGTTCATTGATGTCATCAATAAGCTAAAGGAGAATGTTGAAGAATTAATTGAA
GACAAATACAAAATAATCCTAGAGAAGAATGATACTAAAAAGACATTGCAGAATTTGGAAGAGGTTTGA
GCTAACACGCAAAAACATCTTCAGGAATCCAGGAATGACAAGGAAATGTTACAGCTTCAATTTAAGAAG
ATCAAGGCTAATTATGTGTGTTTACAGGAAAGGTACATGACTGAAATGCAACAAAAAATAAATCTGTA
AGTCAGTATTTAGAGATGGACAAAACCTTAAGCAAGAAAGAAGAGGTTAGAGAGACTACAACAACCTC
AAAAAGAAGTGGAAAAGGCCACAGCTTCTGCTTTGGACTTGTGAAACGGGAAAAAGAGGCCCAAGAA
CAAGAGTCTTGTCTTTACAGGAGGAATCCAGAACTTGAAAAGGAAAACCTGGAAGAAAGACAGAAA
CTGAAATCTAGACTTGAGAAATTGCTCACTCAAGTTAGAAATTTGCAATTTATGTCTGAAATGAAAGA
ACGAAGAATATAAACTTCAGCAGCAAAATCAACGAAGTAAAAATGAGAATGCAAACTTAACAGCAG
GTTGCAAGGAGTGAAGAGCAAAATTTATGCCCTAAATTTGAGACAGCTCAGTTAAAGGATCAATTAGAG
GAAGTCTTGAAGTCAGATATTACCAAGGATACAAAACGACACATTCTAATCTGCTCCCGGATTGCTCA
CCTTGTGAAGAGAGGCTGAATCCTGCAGATATAAAAGAGCTTCTCAGCTGGCCTCCAAATGCACAGT
CTTCTGGCTCTGATGGTGGGACTTCTCACATGCCAGGACATCATCAATTCTGATGCTGAACATTTCAA
GAGAGTGAGAAGGTTAGTGATATAATGTGCAAAAAGTGAAGAGCCTCCATCTTAAAAAGAAAATTTA
GATAAAGAGGTCATAGACTGTGATTGATGAAGCAAGAGTATCAGAGATGTACCTACCCTTCTGGGA
GCCAACTGGATAAGTACCACAGTCTAATGAGGAGCTTGATTTCTTGGTTACATCATATGAAGAAATC
ATTGAATGTGCTGACCAAGGCTTGCAATATCCCACTCCAGATTGCACATTTGGAAGAGAGAAATAAA
CATTTAGAGGATTTAATTAGAAAGCCCAGAGAAAAGCCAGAAAACCAAGATCAAAAAGCTTAGAAAAT
CATCCGAGTCCATGACCATGATGCCAGCTCTTTTAAAGAAAATAGAAATGATTAGATTAA
ACGGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_205864

Insert Size: 1926 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_205864.2](#)

RefSeq Size: 2851 bp

RefSeq ORF: 1926 bp

Locus ID: 285782

UniProt ID: [Q8TC20](#)

Cytogenetics: 6p24.3

MW: 74.7 kDa