

## Product datasheet for **SC308347**

### LTK (NM\_206961) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	LTK (NM_206961) Human Untagged Clone
Tag:	Tag Free
Symbol:	LTK
Synonyms:	TYK1
Vector:	<u>pCMV6 series</u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_206961, the custom clone sequence may differ by one or more nucleotides

```

ATGGGGCTGCTGGGGACAGCTGCTGGTGTGGTTCGGAGCCGCGGGCGCCATTCTCTGCTCT
AGCCCCGGGTCCCAGGAGACTTTTCTGCGGTCTCGCCCTGCCGCTGGCAAGTCCCAGC
CCCCGGGACCCGAAAGTCAGCGCCCCGCTAGTATCTTGGAGCCAGCCTCCCCGCTGAAT
TCTCCGGGCACCGAGGGGTCTTGGCTGTTTTCTACCTGCGGGGCCAGCGCCGGCATGGG
CCCACACAGACACAATGTGACGGGGCGTACGCGGGGACCAGCGTGGTGGTACCCTGGGG
GCCGCCCCGGCAGCTGAGAGGCGTGCAGCTGTGGCGCGTGCCGGGCCCTGGCCAGTATCTG
ATCTCAGCCTACGGAGCCGCGGGCGCAAGGCGCCAAGAACCACCTGTGCGGGGCGCAT
GGCGTCTTCGTCTCAGCAATCTTCTCCCTCGGTCTCGGGGAGTCGCTGTACATCCTGGTG
GGGCAGCAGGGAGAGGACGCCTGTCCCGGAGGTAGCCCGAGAGCCAGCTCGTCTGCCTC
GGGGAGTCTCGAGCCGTTGAAGAGCACGCGGCGATGGATGGGAGCGAAGGGGTCCCAGGG
TCGCGGCGCTGGGCGGGAGGTGGCGGGGTGGCGGGGGCCACCTACGTTTTCCGGGTG
CGCGCTGGCGAGCTGGAACCGTTGCTGGTGGCGGCGGAGGCGGGGTGCGGCCTACCTG
AGGCCGCGGGACCGAGGCCGACTCAGGCCTCCCCGAGAACTGGAGAACCCTCGGAG
GCGCCCCGGGAGCGCGGGAGAGGCGGGCGGCAGGGGGCGACGTTTCAGAGACTGACAAC
CTCTGGGCTGATGGGGAAGATGGAGTATCCTTCATACACCCAGCAGCGAGCTCTTCTCTG
CAGCCTCTGGCAGTCACCGAGAACCACGGAGAGGTAGAGATCCGAAGGCACCTCAACTGC
AGTCACTGCCCTTGGAGACTGCCAATGGCAGGCAGAGCTCCAGCTGGCTGAATGCCTG
TGCCAGAAAGGCATGGAGCTAGCTGTGGATAACGTCACCTGCATGGACCTGCACAAGCCC
CCAGGCCCTCTGGTCTGATGGTGGCTGTGGTGGCAACCTCAACTGAGCCTCCTTATG
GTGTGTGGGTCTGATTCTGGTGAAGCAGAAGAAGTGGCAGGGCCTGCAGGAGATGAGG
CTGCCAGCCCTGAGCTTGGCTGAGCAAGCTTGAACCTTGCATCAGGACAGCCCCC
AATCCCTATTATTGCCAGGTGGGCTTGGCCCGGCCAGTCTGGCCTCTGCCACCAAGT
GTCACCGAGGTTTCCCCAGCAATGTTACTCTGCTCAGAGCCCTGGGCCATGGTGCTTT
GGGGAGGTGTATGAGGGACTGGTAATTGGCCTTCTGGGACTCCAGTCCCCTGCAGGTA
GCTATCAAGACCCTGCCAGAACTCTGCTCGCCTCAGGATGAGCTGGATTTCTCATGGAG
GCCCTCATCATCAGCAAGTTTCCCATCAGAACATTGTGCGGTGTGTGGGGCTCAGCCTC
AGGGCCACCCCTCGCCTATTCTGCTGGAAGTGTCTGGAGGGGACATGAAGAGTTTC
CTGAGGCACAGTCGGCCACACCTGGGCCAGCCATCACCTCTGGTCATGCGGGACCTGCTG
CAACTGGCCCAGGACATAGCCCAGGGCTGCCACTACCTGGAGGAAAATCACTTCATCCAC
AGGGATATTGCCCGCCGAACTGCCTGCTGAGCTGCGCTGGACCCAGCCAGTGGCCAAG
ATTGGGGACTTTGGGATGGCACGAGATATCTACCGGGCCAGTTATTACCGCAGGGGGGAC
CGGGCCTTGCTCCAGTCAAGTGGATGCCCCAGAGGCCTTCTGGAGGGCATCTTCACA
TCCAAGACAGATTCTGGTCTTTTGGGGTGTGCTCTGGGAGATCTTCTCACTGGGCTAC
ATGCCCTATCCTGGGCGCACCAACCAGGAGGTGCTGGACTTCGTCTGGAGGAGGCCGG
ATGGACCCTCCTAGGGGCTGCCAGGCCTGTGTACCGCATCATGACCCAGTGTGGCAG
CAGGACCCGGATGTGCTGAATTCACCTCTGCCAATGGAGCTGGGGCCACCCAGAGGAG
GAAGGGACTTCTGGGCTGGGGAACAGATCTTTGGAGTGCCTAAGACCCCCACAGCCCCAG
GAACTGAGTCCAGAGAAGTTGAAAAGCTGGGGAGGTAGCCCTCTTGGCCCTGGCTGTCC
TCTGGCCTCAAGCCCTCAAATCCAGGGCCTCAAACCTCAGAACCTTTGGAATCCCACT
TATCGCTCCTGA
    
```

**Restriction Sites:** Please inquire

**ACCN:** NM\_206961

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_206961.1</a></u> , <u><a href="#">NP_996844.1</a></u>
<b>RefSeq Size:</b>	2952 bp
<b>RefSeq ORF:</b>	2412 bp
<b>Locus ID:</b>	4058
<b>UniProt ID:</b>	<u><a href="#">P29376</a></u>
<b>Cytogenetics:</b>	15q15.1
<b>Protein Families:</b>	Druggable Genome, Protein Kinase, Transmembrane
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a member of the ros/insulin receptor family of tyrosine kinases. Tyrosine-specific phosphorylation of proteins is a key to the control of diverse pathways leading to cell growth and differentiation. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]</p> <p>Transcript Variant: This variant (2) is lacking one of the internal coding exons compared to variant 1. This results in an isoform (2) that is missing a 61 aa segment compared to isoform 1.</p>