

Product datasheet for **SC308411**

ADAM15 (NM_207194) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ADAM15 (NM_207194) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADAM15
Synonyms:	MDC15
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_207194, the custom clone sequence may differ by one or more nucleotides

```

ATGCGGCTGGCGTGCTCTGGGCCCTGGGGCTCCTGGGCGGGGAGCCCTCTGCCTTCC
TGGCCGCTCCCAAATATAGGTGGCACTGAGGAGCAGCAGGCAGAGTCAGAGAAGGCCCG
AGGGAGCCCTTGGAGCCCAGGTCTTCAGGACGATCTCCAATTAGCCTCAAAAAGGTG
CTTCAGACCAGTCTGCCTGAGCCCCTGAGGATCAAGTTGGAGCTGGACGGTGACAGTCAT
ATCCTGGAGCTGCTACAGAATAGGGAGTTGGTCCCAGGCCGCCAACCCCTGGTGTGGTAC
CAGCCCGATGGCACTCGGGTGGTCAGTGAGGGACACACTTTGGAGAACTGCTGTACCAG
GGAAGAGTGGGGGATATGCAGGCTCCTGGGTGTCCATCTGCACCTGCTCTGGGCTCAGA
GGCTTGGTGGTCTGACCCAGAGAGAAGCTATACCCTGGAGCAGGGGCTGGGGACCTT
CAGGGTCTCCCATTTTCGGAATCCAAGATCTCCACCTGCCAGGCCACACCTGTGCC
CTGAGCTGGCGGAATCTGTACACACTCAGAAGCCACCAGAGCACCCCTGGGACAGCGC
CACATTCGCCGGAGGCGGGATGTGGTAACAGAGACCAAGACTGTGGAGTTGGTATTGTG
GCTGATCACTCGGAGGCCAGAAATACCGGACTTCCAGCACCTGCTAAACCGCACACTG
GAAGTGGCCCTCTTGTGGACACATTCTCCGGCCCTGAATGTACGAGTGGCACTAGTG
GGCCTGGAGGCCTGGACCCAGCGTGACCTGGTGGAGATCAGCCAAACCCAGCTGTACC
CTCGAAAACCTTCCACTGGCGCAGGGCACATTTGCTGCCTCGATTGCCCATGACAGT
GCCAGCTGGTACTGTTTCTCTGGGCTACGGTGGGCATGGCCATTCAGAAC
TCCATCTGTTCTCCTGACTTCTCAGGAGGTGTGAACATGGACCACTCCACCAGCATCCTG
GGAGTCGCTCTCCATAGCCATGAGTTGGGCCACAGCTGGGCTGGACCATGATTTG
CCTGGGAATAGCTGCCCTGTCCAGGTCCAGCCCAGCCAAGACTGCATCATGGAGGCC
TCCACAGACTTCTACCAGGCCTGAACCTCAGCAACTGCAGCCGACGGGCCCTGGAGAAA
GCCCTCTGGATGGAATGGGAGCTGCCTCTTGAACGGCTGCCTAGCCTACCCCTCCTG
GCTGCTTTCTGCGGAAATATGTTTGTGGAGCCGGGCGAGCAGTGTGACTGTGGCTTCTG
GATGACTGCGTCGATCCCTGCTGTGATTCTTTGACCTGCCAGCTGAGGCCAGGTGCACAG
TGTGCATCTGACGGACCTGTTGTCAAATGCCAGCTGCGCCCGTCTGGCTGGCAGTGT
CGTCTACCAGAGGGGATTGTGACTTGCCTGAATTCTGCCAGGAGACAGCTCCAGTGT
CCCCCTGATGTCAGCCTAGGGGATGGCGAGCCCTGCGCTGGCGGGCAAGCTGTGTGCATG
CACGGGCGTTGTGCTCCTATGCCAGCAGTGCCAGTCACTTTGGGGACCTGGAGCCAG
CCCCTGCGCCACTTTGCCTCCAGACAGCTAATACTCGGGGAAATGCTTTTGGGAGCTGT
GGGCGCAACCCAGTGGCAGTTATGTGCTCCTGCACCCCTAGAGATGCCATTTGTGGGCG
CTCCAGTGCCAGACAGGTAGGACCCAGCCTCTGCTGGGCTCCATCCGGGATCTACTCTGG
GAGACAATAGATGTGAATGGGACTGAGCTGAACTGCAGCTGGGTGCACCTGGACCTGGGC
AGTGATGTGGCCAGCCCTCCTGACTCTGCCTGGCACAGCCTGTGGCCCTGGCCTGGTGT
TGTATAGACCATCGATGCCAGCGTGTGGATCTCCTGGGGGCACAGGAATGTGAAGCAA
TGCCATGGACATGGGGTCTGTGACAGCAACAGGCACTGCTACTGTGAGGAGGGCTGGGCA
CCCCCTGACTGCACCACTCAGCTCAAAGCAACCAGCTCCCTGACCACAGGGCTGCTCCTC
AGCCTCCTGGTCTTATTGGTCTGGTGTGCTTGGTGGCAGCTACTGGTACCGTGCCCGC
CTGCACCAAGCAGCTCTGCCAGCTCAAGGGACCCACCTGCCAGTACAGGGCAGCCCAATCT
GGTCCCTCTGAACGGCCAGGACCTCCGAGAGGGCCCTGCTGGCACGAGGCACTAAGCAG
GCTAGTGTCTCAGCTTCCCGGCCCCCTTCCAGGCCGCTGCCGCTGACCCCTGTGTCC
AAGAGACTCCAGTCTCAGGGGCCAGCCAAGCCCCACCCCAAGGAAGCCACTGCCTGCC
GACCCCAAGGGCCGTTGCCATCGGGTACCTGCCCGGCCAGGGGCTGGAATCCCGCC
CTAGTGGTACCCTCCAGACCAGCGCCACCCTCCGACAGTGTCTCGCTCTACCTCTGA

```

Restriction Sites: Please inquire
ACCN: NM_207194

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207194.1, NP_997077.1</u>
RefSeq Size:	2898 bp
RefSeq ORF:	2520 bp
Locus ID:	8751
UniProt ID:	<u>Q13444</u>
Cytogenetics:	1q21.3
Protein Families:	Druggable Genome, Protease, Transmembrane
Gene Summary:	<p>The protein encoded by this gene is a member of the ADAM (a disintegrin and metalloproteinase) protein family. ADAM family members are type I transmembrane glycoproteins known to be involved in cell adhesion and proteolytic ectodomain processing of cytokines and adhesion molecules. This protein contains multiple functional domains including a zinc-binding metalloprotease domain, a disintegrin-like domain, as well as a EGF-like domain. Through its disintegrin-like domain, this protein specifically interacts with the integrin beta chain, beta 3. It also interacts with Src family protein-tyrosine kinases in a phosphorylation-dependent manner, suggesting that this protein may function in cell-cell adhesion as well as in cellular signaling. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3) lacks an in-frame coding exon compared to variant 6. The resulting isoform (3) lacks an internal region, as compared to the longest isoform 6.</p>