

Product datasheet for SC308425

MBNL1 (NM_207293) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MBNL1 (NM_207293) Human Untagged Clone
Tag:	Tag Free
Symbol:	MBNL1
Synonyms:	EXP; MBNL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC308425 representing NM_207293. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGCATCGCC
ATGGCTGT TAGTGT CACACCAATTCGGGACACAAAATGGCTAACACTGGAAGTATGTAGAGAGTTCCAG
AGGGGGACTTGCTCACGGCCAGACACGGAATGTAATTTGCACATCCTTCGAAAAGCTGCCAAGTTGAA
AATGGACGAGTAATCGCCTGCTTTGATT CATTGAAAGGCCGTTGCTCCAGGGAGAAGCTGCAAAATATCTT
CATCCACCCCCACATTTAAAAACGCAGTTGGAGATAAATGGACGCAATAACTTGATT CAGCAGAAGAAC
ATGGCCATGTTGGCCAGCAAATGCAACTAGCCAATGCCATGATGCCTGGTGCCCAATTACAACCCGTG
CCAATGTTTT CAGTTGCACCAAGCTTAGCCACCAATGCATCAGCAGCCGCTTTAATCCCTATCTGGGA
CCTGTTTCTCCAAGCCTGGTCCCAGGAGATCTTGCCGACTGCACCAATGTTGTTACAGGGAATCCG
GGTGTCCCTGTACCTGCAGCTGCTGCAGCTGCTGCACAGAAATTAATGCGAACAGACAGACTTGAGGTA
TGTCGAGAGTACCAACGTGGCAATTGCAACCGAGGAGAAAATGATTGTCGGTTTGCTCATCTGCTGAC
AGCACAATGATTGACACCAATGACAACACAGTCACTGTGTGTATGGATTACATCAAAGGGAGATGCTCT
CGGGAAAAGTGCAAACTTTTATCCCCCTGCACATTTGCAAGCCAAGATCAAGGCTGCCAATACCAG
GTCAACCAGGCTGCAGCTGCACAGGCTGCAGCCACCGCAGCTGCCATGACTCAGTCGGCTGTCAAATCA
CTGAAGCGACCCCTCGAGGCAACCTTTGACCTGGGAATTCCTCAAGCTGTAATTTCCCAATACCAACAG
AGGCCTGCTCTTGAAAAACCAACGGTGCCACCGCAGTCTTTAACACTGGTATTTTCCAATACCAACAG
GCTCTAGCCAACATGCAGTTACAACAGCATACAGCATTTCTCCCACAGTTCCCATGGTGACGGTGCT
ACGCCAGCCACTGTGTCGCAGCAACAACATCTGCCACAAGTGTTCCTTCGCTGCAACAGCCACAGCC
AACCAGATACCCATAATATCTGCCGAACATCTGACTAGCCACAAGTATGTTACCCAGAT GTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

ACCN:	NM_207293
Insert Size:	1167 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207293.1</u>
RefSeq Size:	5390 bp
RefSeq ORF:	1167 bp
Locus ID:	4154
UniProt ID:	<u>Q9NR56</u>
Cytogenetics:	3q25.1-q25.2
MW:	41.8 kDa

Gene Summary:

This gene encodes a member of the muscleblind protein family which was initially described in *Drosophila melanogaster*. The encoded protein is a C3H-type zinc finger protein that modulates alternative splicing of pre-mRNAs. Muscleblind proteins bind specifically to expanded dsCUG RNA but not to normal size CUG repeats and may thereby play a role in the pathophysiology of myotonic dystrophy. Mice lacking this gene exhibited muscle abnormalities and cataracts. Several alternatively spliced transcript variants have been described but the full-length natures of only some have been determined. The different isoforms are thought to have different binding specificities and/or splicing activities. [provided by RefSeq, Sep 2015]

Transcript Variant: This variant (3) contains one alternate in-frame exon and lacks an alternate in-frame exon in the 3' coding region, compared to variant 1. The encoded isoform (c) is longer than isoform a. Other names for this isoform are EXP42 and MBNL1subscript42.

Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.