

Product datasheet for **SC308478**

GADL1 (NM_207359) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GADL1 (NM_207359) Human Untagged Clone
Tag:	Tag Free
Symbol:	GADL1
Synonyms:	ADC; CSADC; HuADC; HuCSADC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC308478 representing NM_207359.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGCAGCGACTCGGACCGCCAGTGCCTGTGGACGGAGATATTGATCAACAAGAGATGATTCCAAGT
AAGAAGAATGCTGTTCTGTGGATGGGGTTGTGCTGAATGGTCCACAACAGATGCAAAAGCTGGAGAA
AAATTTGTTGAAGAGGCCTGTAGGCTAATAATGGAAGAGGTGGTTTTGAAAGCTACAGATGTCAATGAG
AAGGTGTGTGAATGGAGGCCTCCTGAACAACGAAACAGCTTCTTGATTTGGAGATGAGAGACTCAGGC
GAGCCACCCATAAACTATTGGAACCTGTGCGGGATGTCATACACTACAGTGTCAAACTAACCACCCA
AGATTTTTCAACCAATTGTATGCTGGACTTGATTACTCCTTGGTGGCCCGATTTATGACCGAAGCA
TTGAATCCAAGTGTTTATACGTATGAGGTGTCCCAGTGTCTGTTAGTGAAGAAGCGGTTCTGAAG
AAAATGATTGAATTTATTGGCTGAAAAGAAGGGGATGGAATTTAACCCAGGTGGCTCAGTGTCCAAT
ATGTATGCAATGAATTTAGCTAGATACAAATATTGTCCTGATATTAAGGAAAAGGGGCTGTCTGGTTCG
CCAAGATTAATCCTTTTACATCTGCAGAGTGCATTACTCTATGAAGAAGGCAGCCTCTTTTCTGGG
ATTGGCACTGAGAATGTTTGCTTTGTGGAAACAGATGGAAGAGGTAAAATGATACCTGAGGAACTGGAG
AAGCAAGTCTGGCAAGCCAGAAAAGAGGGGGCAGCACCGTTTTCTGTCTGTGCCACTTCTGGTACAAC
GTGTTGGGAGCTTTTGACCCTCTGGATGAAATAGCAGACATCTGCGAGAGGCACAGCCTCTGGCTTCAT
GTAGATGCTTCTTGGGGTGGCTCAGCTTGTATGTCGAGGAAGCACCGCAAGCTTCTGCATGGCATCCAC
AGGGCTGACTCTGTGGCCTGGAACCCACACAAGATGCTGATGGCTGGGATCCAGTGTGTCTCTCCTT
GTGAAAGACAAATCTGATCTTCTTAAAAATGCTACTCTGCCAAGGCATTTACCTTCCAGCAGGAT
AAATCTATGATGTGAGCTATGACACAGGAGACAAGTCTATCCAGTGTAGCAGAAGACCAGATGCATTC
AAGTTCTGGATGACCTGGAAGGCCCTGGGTACATTAGGCCTTGAAGAAAGAGTTAATCGTCTCTTGCT
TTATCTAGGTACCTAGTAGATGAAATCAAGAAAAGAGAAGGATTCAAGTTACTGATGGAACCTGAATAT
GCCAATATTTGCTTTTGGTACATTCCACCAGCCTCAGAGAGATGGAAGAAGACCAGGTTCTGGGCA
AACTTAATTTGGTGGCCCGCCAGCCATTAAGGAGAGGATGATGAAGAAGGGAAAGCTTGATGCTGGGCTAC
CAGCCGCACCGGGAAAGGTCAACTTCTCCGCCAGGTGGTGTATCAGCCCTCAAGTGAGCCGGGAGGAC
ATGGACTTCTCCTGGATGAGATAGACTTACTGGGTAAGACATGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_207359

Insert Size: 1566 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_207359.2
RefSeq Size:	3655 bp
RefSeq ORF:	1566 bp
Locus ID:	339896
UniProt ID:	Q6ZQY3
Cytogenetics:	3p24.1-p23
Protein Families:	Druggable Genome
MW:	59.2 kDa
Gene Summary:	May catalyze the decarboxylation of L-aspartate, 3-sulfino-L-alanine (cysteine sulfinic acid), and L-cysteate to beta-alanine, hypotaurine and taurine, respectively. Does not exhibit any decarboxylation activity toward glutamate.[UniProtKB/Swiss-Prot Function]