

Product datasheet for SC308651

PATE2 (NM_212555) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PATE2 (NM_212555) Human Untagged Clone
Tag:	Tag Free
Symbol:	PATE2
Synonyms:	C11orf38; PATE-M; UNQ3112
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC308651 representing NM_212555. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGCTTGTCTCTTTCTCCTGGGCACAGTCTTTCTGCTCTGCCCATATTGGGGTGAACCTTCATGACCCCT
 ATAAAAGCGACTGAAATAATGTGTATGAATGTAAAAATATCATCTTGGGTTATGCTATGGTGTCATG
 ACATCCTGCTCCCTGAAGCATAAACAGTCCTGTGCAGTTGAGAACTTTACATCCTTACAAGGAAAGGG
 CAGAGCATGTATCATTATTCAAACGTGCTGTATGACCAGCTGTGAGGACATCAACTTCCTGGGGTTC
 ACGAAGAGGGTAGAGCTCATCTGTTGTGATCATAGTAAGTACTGCAACCTCCCTGAGGGAGTTTAG
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_212555
Insert Size:	342 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_212555.2</u>
RefSeq Size:	1640 bp
RefSeq ORF:	342 bp
Locus ID:	399967
UniProt ID:	<u>Q6UY27</u>
Cytogenetics:	11q24.2
Protein Families:	Secreted Protein
MW:	13 kDa