

Product datasheet for **SC308667**

ZNF543 (NM_213598) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF543 (NM_213598) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF543
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC308667 representing NM_213598.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGCCTCGCGCAGGTGTCTGTGACCTTTGAGGATGTGGCTGTGACATTCACCCAGGAGGAGTGG
GGACAGTTGGATGCAGCCAGAGAACCTTGTATCAGGAGGTGATGCTGGAGACCTGCCGACTTCTCATG
TCTCTGGGCTGTCTTTGTTCAAACCAGAGCTGATCTACCAGTTGGATCAGAGACAGGAGCTATGGATG
GCTACAAAAGACCTCTCCCAAAGCTCCTATCCAGGTGACAACACAAAACCAAGACCACAGAGCCTACC
TTTTCTCACCTGGCCTTGCCGTGAGGAAGTCTTACTCCAGGAACAACCTGACACAAGGAGCCTCAAAGAAC
TCCCAATTAGGGCAATCCAAGGATCAGGATGGGCCATCTGAAATGCAAGAAGTCCACTTGAAAAATAGGG
ATAGGCCCCAGCGGGGAAGCTGCTGGAGAAAATGAGTTCTGAACGTGATGGTTTGGGGTCAGATGAT
GGTGTATGTACAAAGATTACACAGAAACAAGTTTCAACAGAAGGTGATCTCTATGAATGTGATTCACAT
GGACCAGTTACAGATGCCTTGATTCGCGAAGAGAAAAATTCCTATAAATGTGAGGAATGCCGGAAGTG
TTTAAAAAGAATGCCCTCCTTGTTTCAGCATGAACGGATTACACTCAAGTGAAGCCCTATGAATGCACA
GAGTGTGGGAAAACCTTTAGCAAGAGCACTCATCTTCTCAGCACCTCATCATCCACACTGGGGAGAAG
CCCTATAAGTGCATGGAGTGTGGGAAGGCTTTTAAACCGCAGGTCACACCTCACACGGCACCAGCGGATT
CACAGTGGAGAGAAGCCTTATAAGTGCAGTGAATGTGGAAGGCCTTCACCCACCCTCCACTTTTGTC
TTGCATCAGAGGACCACTGGAGAAAAACCTTTGTGTGCAAGAGTGTGGCAAAGCCTTTTCGAGAT
AGGCCAGGTTTCATTGACACTACATCATCCACACGGGAGAGAAGCCCTATGAGTGCATTGAGTGTGGG
AAGGCCTTCAACCGCGGTCATACCTCACGTGGCACCAACAGATTCACACTGGAGTGAAACCTTTGAA
TGCAACGAGTGTGAAAAGCTTTTTCGAGAGTGCAGACCTCATTCAACACTACATTATCCACACTGGG
GAGAAGCCCTATAAGTGCATGGAGTGTGGGAAGGCTTCAACCGTAGGTACACCTCAAGCAGCATCAA
CGGATTCACACTGGGAGAAAGCCTTATGAATGCAGTGAATGTGGAAAGGCCTTCACCCACTGCTCCACT
TTGTCTTGCAATAAAGGACCCACACAGGAGAAAAACCTATGAATGCAAGAATGTGAAAAGCCTTT
AGTGATAGGGCAGACCTCATTGCGCACTTCAGCATCCACACTGGAGAGAAACCTATGAGTGCCTGGAG
TGTGGAAGGCCTTCAACCGCAGCTCACACCTCACGAGGCACCAACAGATTCACACTGGAGAGAAACCC
TATGAATGCATCCAGTGTGGAAAGCCTTTGCCGAGCGCAACCTTATTCGACACTCCATCATTAC
ACTGGAGAGAAGCCGTATGAATGCAGTGAATGTGGAAAGGCTTTTAAATCGCGGCTCATCCCTCACACAT
CATCAAAGGATTCATACTGGGAGAAACCTACCATTGTAAACAGATGTGGGAAGACCTTTTATGACTGCA
CAGACTTCAGTCAACATCCAGGAACCTTTATTAGGGAAGAGTTTTTGAATATCACCCTGAAGAAAAAT
CTGTGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_213598

Insert Size: 1803 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_213598.3</u>
RefSeq Size:	3671 bp
RefSeq ORF:	1803 bp
Locus ID:	125919
UniProt ID:	<u>Q08ER8</u>
Cytogenetics:	19q13.43
MW:	68.4 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]