

Product datasheet for **SC308671**

ZNF789 (NM_213603) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF789 (NM_213603) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF789
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC308671 representing NM_213603. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGTTCCACACGAGGGGGAAGGAGCTGCTCTCGTTTGAGGATGTGGCGATGTAATTCACAGAGAG
GAGTGGGGCCACCTCAACTGGGGTCAGAAGGACCTCTACCGAGATGTGATGTTGGAGAATACAGGAAC
ATGGTCTTGCTGGGATTTCACTTTCCAAACCTGAGATGATCTGTGAGCTGGAGAATGGGACGAGCAG
TGGATCCTGGATCTACCGAGAATGGGAATAGGAAGGCTCCGGTAGTGCTTGCCAGGTTCTGAAGCC
AGACACAAGATGAAAAAGCTAACTCCAAACAGAAATTTCTGAAGATTTAGAGTCATATAAGATATCA
GTGGTAATGCAGGAATCAGCTGAGAACTTTGAGAAAAGTTACATAAGTGTAAGAATTTGTGGACAGT
TGCAGGCTTACTTTCCCTACTAGTGGTGATGAATACAGCAGGGGCTTCCTTCAAAACCTTAACCTTATT
CAAGATCAGAATGCGCAACAAGGTGGAAGCAGGGCAGATATGATGAGGATGGCAAAACCTTCAATCAA
AGATCTTTGCTTTTGGGCATGAGCGAATTCTCACAAGAGCAAAGTCTTATGAATGCAGTGAATGTGGA
AAAGTCATTAGGCGTAAGGCATGGTTTGATCAACATCAAAGAATTCACCTTTTAGAGAATCCTTTTGAG
TGTAAGGTCTGTGGCAAGCCTTCAGACAGCGGTGAGCTCTTACGGTCCATAAACAGTGTACCTGCAA
AACAGCCATACAGATGTCATGACTGTGGAAGTGTTTTCGGCAGCTCGCGTATCTTGTGAACATAAG
AGGATTCACACCAAGAAAAACCTTATAAATGTAGCAAATGTAAAAACGTTTAGTCAGAATTCACCC
CTTATTCGACATCAGGTGATCCATAGTGGAGAAAAACGCCATAAATGCCTTGAGTGTGAAAAAGCCTTT
GGCCGGCATTCAACCTTCTATGTCATCAACAGATTACAGTAAACCGAACCCATAAATGCAGTGAA
TGTGGACAGTCCTTTGGTAGGAATGTGGATCTCATTGAGCATCAAAGAATCCATACAAAGGAGGAATTC
TTTCAATGTGGAGAATGTGGGAAACGTTTAGTTTTAAGAGGAATCTTTTCGACATCAGGTCAATCAC
ACTGGAAGCCAACCTACCAATGTGTCATATGTGGAAAATCTTTCAAGTGGCACACAAGCTTTATTAAG
CACCAGGGCACTCACAAGGACAGATATCCACA**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_213603
Insert Size:	1278 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_213603.2</u>
RefSeq Size:	1654 bp
RefSeq ORF:	1278 bp
Locus ID:	285989
UniProt ID:	<u>Q5FWF6</u>
Cytogenetics:	7q22.1
MW:	50 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).