

## Product datasheet for **SC308833**

### RHCE (NM\_020485) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	RHCE (NM_020485) Human Untagged Clone
Tag:	Tag Free
Symbol:	RHCE
Synonyms:	CD240CE; RH; Rh4; RH30A; RHC; RHCE(152N); RHE; RhIVb(J); RHIXB; RHNA; RHPI; RhVI; RhVIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC308833 representing NM_020485. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGCATCGCC
ATGAGCTCTAAGTACCCGCGGTCTGTCCGGCGCTGCCTGCCCTCTGGGCCCTAACACTGGAAGCAGCT
CTCATTCTCCTCTCTATTTTTTTACCCACTATGACGCTTCCTTAGAGGATCAAAAGGGGCTCGTGGCA
TCCTATCAAGTCGGCCAAGATCTGACCGTGATGGCGGCCCTTGGCTTGGGCTTCTCACCTCAAATTTT
CGGAGACACAGCTGGAGCAGTGTGGCCTCAACCTCTTCATGCTGGCGCTTGGTGTGCAGTGGGCAATC
CTGCTGGACGGCTTCTGAGCCAGTTCCTCCTGGGAAGGTGGTCATCACACTGTTCAAGTATTGGCTG
GCCACCATGAGTGCTATGTCGGTCTGATCTCAGCGGGTCTGTCTTGGGAAGGTCAACTTGGCGCAG
TTGGTGGTGTGGTGGTGGAGGTGACAGCTTAGGCACCCTGAGGATGGTCATCAGTAATATCTTC
AACACAGACTACCACATGAACCTGAGGCACTTCTACGTGTTTCGACGCTATTTTGGGCTGACTGTGGCC
TGGTGCCTGCCAAAGCCTCTACCCAAGGGAACGGAGGATAATGATCAGAGAGCAACGATACCCAGTTTG
TCTGCCATGCTGGGCGCCCTCTTCTGTGGATGTTCTGGCCAAGTGTCAACTCTGCTCTGCTGAGAAGT
CCAATCCAAGGAAGAATGCCATGTTCAACACCTACTATGCTCTAGCAGTCAGTGGTGCAGCCATC
TCAGGGTCATCCTTGGCTCACCCCAAAGGAAGATCAGCATGACTATGTGCACAGTGGGCTTGGCA
GGAGGCGTGGCTGTGGTACCTCGTGTACCTGATCCCTTCTCCGTGGCTTGGCATGGTGTGGGATT
GTGGCTGGGCTGATCTCCATCGGGGAGCCAAGTGCCTGCCGGTGTGTTGTAACCGAGTCTGGGATT
CACCACATCTCCGTGATGCACTCCATCTTCAGCTTGTGGTCTGCTTGGAGAGATCACCTACATTGTG
CTGCTGGTCTTCACTGTCTGGAACGCAATGGCATGATTGGCTTCCAGGTCCTCCTCAGCATTGGG
GAACTCAGCTTGGCCATCGTGATAGCTCTCACGCTGGTCTCCTGACAGGTTTGCTCCTAAATCTCAA
ATATGGAAGCACCTCATGTGGCTAAATATTTTGATGACCAAGTTTTCTGGAAGTTTCTCATTTGGCT
GTTGGATTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>Plasmid Map:</b>	□
<b>ACCN:</b>	NM_020485
<b>Insert Size:</b>	1254 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_020485.4</a>
<b>RefSeq Size:</b>	1618 bp
<b>RefSeq ORF:</b>	1254 bp
<b>Locus ID:</b>	6006
<b>UniProt ID:</b>	<a href="#">P18577</a>
<b>Cytogenetics:</b>	1p36.11
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	45.5 kDa

**Gene Summary:**

The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene which encodes both the RhC and RhE antigens on a single polypeptide and a second gene which encodes the RhD protein. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes. A mutation in this gene results in amorph-type Rh-null disease. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms. [provided by RefSeq, Aug 2016]

Transcript Variant: This variant (1) encodes the longest isoform of this protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.