

Product datasheet for SC308867

IQCC (NM_018134) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	IQCC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC308867 representing NM_018134. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGCCAGAGCTGCTGGTTCGGAAGGTGCTGCATTGCAGGCCTGCGTCCGGGGCTTCTTGGTCCGA
CGCCAGTTCAGAGCCTGCGAGCTGAGTATGAGGCGATTGTACGAGAGGTGAGGGCGACCTGGGCACG
CTTCAGTGGACCGAGGGCCGATTCCAGGCCGCGATTCTCCAGAGAAGGCAAAATCCCATCAGACC
TGGAAAGCAGGAGACAGGGTAGCAATCCAGAGCAGGGGCTGTGGAACCACTTCCCATGTGAAGAGTCT
GAGGGAGAGGCCACCTGGGAGGAGATGGTGTCTGAAGAAGTCAGGAGAGAGCTCAGCAAAATCAAGGAAGC
CTCTGCAGAGATCACAGCTCCTGGCTTCAGATGAAGCAGAACAGGAACCCAGCCAAGAGAAGACCAGA
GACACGACAAGGATGGAGAATCCAGAAGCCACAGATCAAAGACTGCCCCACAGCCAACCTCAGCTTCAA
GAGCTTCAGTACCACCGCAGCCACTTGGCCATGGAATTGCTGTGGCTGCAACAGGCCATCAATAGCCGT
AAGGAGTACCTACTTCTTAAACAAACACTGAGATCCCCAGAGGCGGGCCGATCAGAGAGGAACCCCGC
GTGTTCTAGAACATGGGAACAGGCCTGTGAGAGGGACAGTCACAACCAAGCGCACCACTGGAGGAC
CAGTCTACAGAGACAGGACCACTGGAGAGCTAGAACAGGAGGATGACTCCTGTACAGGGTCAAATCA
CCCCACAGATCCCCAGGAAGTTTGGCCACTACACAAAAAACATTGCTGGGGCTAAGTGCAGAGAACCA
TGCTACAGCAAGTCTGGACCACCGTCGTCTATACCATCAAACAGCCAGGCCTTGGGGGACAGGCTCACC
AAAGGGCCAGACGATGGAAGACAGACCTTTGGAGGGACCTGCCTGCTGCAGATGAAATCCTGGAGGAC
CAGACCCCCAGAGGTTTAAACCTAGGAACCATGTCCCAGGAAGTCCAGGACACAGCTGTCTGCACTC
TATGAGGACTCAAATATTAAGGAGATGTCTCCAGAAAACCTAGACCACAAAGAGCCTGACTGCCGAACA
GTCAGGACACAAGATTGGGCTCTCAGAGGACCACATCATCTGGGATGGTACCTTGGGGGGCCAGAG
CATAGTGTCTCGATCTCTGGAGGACTAAACCCCAAGGCCAGGCCCCCACTGATAGAAGCTCCAGA
GATGGAACCTCCAATGAGCCTAGTCATGAAGGACAGAAAAAGCAGAGGACTATACCATGGAGATCAAAG
TCACCTGAGATTCTGTCTTCTACAAAGGCAGGCTGTACAGGAGAGGAACAGTGGAGGGGCAGGCCATGG
AAACAGAACCACTGGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

```



Restriction Sites:	Sgfl-MluI
ACCN:	NM_018134
Insert Size:	1401 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_018134.2</u>
RefSeq Size:	2053 bp
RefSeq ORF:	1401 bp
Locus ID:	55721
UniProt ID:	<u>Q4KMZ1</u>
Cytogenetics:	1p35.2
MW:	53 kDa