

Product datasheet for **SC308963**

DYRK1B (NM_006484) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DYRK1B (NM_006484) Human Untagged Clone
Tag:	Tag Free
Symbol:	DYRK1B
Synonyms:	AOMS3; MIRK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC308963 representing NM_006484.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCCGTCCCACCGGGCCATGGTCCCTTCTCTGGCTTCCAGGGCCCCAGGAGCACACGCAGGTATTG
CCTGATGTGCGGCTACTGCCTCGGAGGCTGCCCTGGCTTCCGGGATGCAACCTCAGCCCCGCTGCGT
AAGCTCTCTGTGGACCTCATCAAGACCTACAAGCACATCAATGAGGTATACTATGCGAAGAAGAAGCGG
CGGGCCCAGCAGGCGCCACCCAGGATTCGAGCAACAAGAAGGAGAAGAAGGTCCTGAACCATGTTAT
GATGACGACAACCATGACTACATCGTGCGCAGTGGCGAGCGCTGGCTGGAGCGCTACGAAATTGACTCG
CTCATTGGCAAAGGCTCCTTTGGCCAGGTGGTAAAGCCTATGATCATCAGACCCAGGAGCTTGTGGCC
ATCAAGATCATCAAGAACAAAAGGCTTTCTGAACAGGCCAGATTGAGCTGCGGCTGCTGGAGCTG
ATGAACCAGCATGACACGGAGATGAAGTACTATATAGTACACCTGAAGCGGCACTTCATGTTCCGGAAC
CACCTGTGCTGGTATTTGAGCTGCTGTCTACAACCTGTACGACCTCCTGCGCAACACCCACTTCCGC
GGCGTCTCGTGAACCTGACCCGGAAGCTGGCGCAGCAGCTCTGCACGGCACTGCTTTTCTGGCCAGC
CCTGAGCTCAGCATCATCTACTGCGACCTCAAGCCCGAAAACATCTTGCTGTGCAACCCCAAGCGCAGC
GCCATCAAGATTGTGGACTTCGGCAGCTCCTGCCAGCTTGGCCAGAGGATCTACCAGTATATCCAGAGC
CGCTTCTACCGCTCACCTGAGGTGCTCCTGGGCACACCCTACGACCTGGCCATTGACATGTGGTCCCTG
GGTGCATCCTTGTGGAGATGCACACCGGAGAGCCCTTTCAGTGGCTCCAATGAGGTCGACCAGATG
AACCGCATTGTGGAGGTGCTGGGCATCCCACCGGCCGCCATGCTGGACCAGGCGCCCAAGGCTCGCAAG
TACTTTGAACGGCTGCCTGGGGTGGCTGGACCCTACGAAGGACGAAAGAACTCAGGAAGGATTACCAG
GGCCCCGGGACACGGCGGCTGCAGGAGGACCTGGTGTGCGCATGCTGGAGTATGAGCCCGCCCGCCGC
ATCAGCCCCCTGGGGCTCTGCAGCACGGCTTCTCCGCCGACGCGCCGACGAGGCCACCAACACGGGC
CCGGCAGGCAGCAGTGCCTCCACCTCGCCCGGCCCTCGACACCTGCCCTCTTCCAGCACCGCCAGC
TCCATCTCCAGTTCTGGAGGCTCCAGTGGCTCCTCCAGTGACAACCGGACCTACCCTACAGCAACCGA
TATTGTGGGGCCCTGGGCCCTTATCACAGACTGTGAGATGAACAGCCCCAGTCCCACCCTCCAG
CCGCTGCGGCCCTGGGCAGGGGTGATGTGCCCAAGACACATCAAGCCCTGCCTCTGCCTCGTCA
CTGCCTGGGACCGGGGCCAGTTACCCCCAGCCCCGATACCTTGGTGTCCCCATCACCAACCTCA
CCACCACCCCGGAGCTGATGGATGTGAGCCTGGTGGCGGCCCTGCTGACTGCTCCCCACCTACCCCA
GCGCCTGCCCCCAGCACCCGGCTGCCTCAGCCCTCCGGACTCGGATGACTGGAGGTCGTCCACCCCTC
CCGCTCTGATGACCCTGCCACTCTGGGGCTCACCTGGGCCTCCGTGGTGTACCCAGAGCACAGCA
GCCAGCTCGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_006484

Insert Size: 1806 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_006484.2
RefSeq Size:	2504 bp
RefSeq ORF:	1806 bp
Locus ID:	9149
UniProt ID:	Q9Y463
Cytogenetics:	19q13.2
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
MW:	66.3 kDa
Gene Summary:	<p>This gene encodes a member of a family of nuclear-localized protein kinases. The encoded protein participates in the regulation of the cell cycle. Expression of this gene may be altered in tumor cells, and mutations in this gene were found to cause abdominal obesity-metabolic syndrome 3. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014]</p> <p>Transcript Variant: This variant (3, also known as 'c') lacks an alternate in-frame segment in the 3' coding region, compared to variant 1. The encoded isoform (p66, also known as isoform c) is shorter than isoform p69.</p>