

## Product datasheet for SC309240

### DIPK2B (NM\_176819) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DIPK2B (NM_176819) Human Untagged Clone
Tag:	Tag Free
Symbol:	DIPK2B
Synonyms:	4930578C19Rik; bA435K1.1; CXorf36; DIA1R; EPQL1862; PRO3743
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC309240 representing NM_176819. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGCCCCAGCTGGGGCTGAGGCTGCCGCCCTCCGCCCTGGCTGGCTGGCCCTGCTGCTGTGGGTC
TCAGCCCTGAGCTGTTCTTTCTCCTTGCCAGCTTCTTCCCTTTCTTCTGCTGGTCCCAAGTCAGAAC
AGCTACAATTTTGAAGGACTTTCTCGGTCTTGATAAATGCAATGCCTGCATCGGGACATCTATTTGC
AAGAAGTTCTTTAAAGAAGAAATAAGATCTGACAACCTGGCTGGCTTCCACCTTGGACTGCCTCCCGAT
TCCTTGCTTTCTATCCTGCAAATTAAGATGATTCCAAAATCTGGCGCCCTGTGGAGATCTTTAGA
CTGGTCAGCAAATATCAAACGAGATCTCAGACAGGAGAAATCTGTGCCTCTGCATCAGCCCCAAAGACC
TGCAGCATTGAGCGTGTCTGCGGAAAACAGAGAGGTTCCAGAAATGGCTGCAGGCCAAGCGCCTCAG
CCGGACCTGGTGCAGGGCCTGGCCAGCCGCTCCTGCGCTGCCCTTCGCAGCGACTCCTGGATCGCGTG
GTCAGGGCCTATGCAGAGGTGGCCGACGCCGAGCATCTTCATGGACCACTTACCAGCCGTGACAAG
CTGCGCCTGCTCTACACGCTGGCTGTCACTCGCACCCCATCCTCCTACAGATCTTCCCTGGGGCTGAG
GGATGGCCGCTGCCAAGTACCTGGGCTCCTGTGGCAGATTCTCGTCAGCACCAGCACCAGACCGCTG
CAGGAATTCTATGATGCACCCCAAGATCAGGCAGCCGACCTTGCCTACCAGCTCCTGGGTGTCCTGGAG
TCTTTGAGGAGCAACGATTTGAACTATTTCTTCTACTTCAACCCACATTGATGCAGGCATGTTCCGGCTC
TTTAAACAACGGGCATCTGTTTCATCCGGGATGCCAGTGCAGTGGGGCTCATCGACAAGCAGGAAGGCAGC
CAAGAAGCCAACAGGGCAGGAGAGAATAAAGACATTTTTAGCTGCCTGGTTTCCGGCTGCCAGGCCAG
CTGCCCTCCTGCCAAAGCATCTCTGAGAAGCAGAGCCTGGTGTGCTGAGAAAGTTGCTGCCTCGA
CTTCTCCAGGGGAGGTTCCCTCCCAAGTCAAGACGACATAGACTCCATCCTTGTTCAGTGTGGGGAC
AGCATCCGCCAGACCCAGAAGTCTTGGGGCCGCCAGCCAGCTGAAAGACATCTTGGGCCCTGAGA
ACGTGTGACTCCAGATTTGCCTATCGTTACCCAGATTGCAAATATAACGATAAGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_176819
<b>Insert Size:</b>	1302 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_176819.3</a></u>
<b>RefSeq Size:</b>	4735 bp
<b>RefSeq ORF:</b>	1302 bp
<b>Locus ID:</b>	79742
<b>UniProt ID:</b>	<u><a href="#">Q9H7Y0</a></u>
<b>Cytogenetics:</b>	Xp11.3
<b>Protein Families:</b>	Secreted Protein
<b>MW:</b>	48.6 kDa