

## Product datasheet for **SC309276**

### **EXOC3L1 (NM\_178516) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	EXOC3L1 (NM_178516) Human Untagged Clone
Tag:	Tag Free
Symbol:	EXOC3L1
Synonyms:	EXOC3L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC309276 representing NM\_178516.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGACTCAGCAGCCAAGGATGAGATGCAGCCGGCGTTGTCCCCTGGACCTGAGTGGCCAGAGCAGGAG
CGGGCAGAGCAGCTGGCCCGGGTGCAGCGCTCAAGTGGGCTCAGGCATCTTCTACCGCCGGAGCAG
CTGGCCAGGCTAGGCCAGTACCGCAGCCGCGAGGTGCAGCGTACCTGCTCCCTGGAATCGCGCCTAAG
TCAGTGATGCAGTCATACCTGGAAGCGTGCAGACTGGTGTGTGGCAGCTGGCCAGGCCATTGAGGTG
GTGCAGGGAACCCGGGAGGCCCTGAGCCAGGCCGTGGTTGCTCCAGGGCATGTCCAGGCCTTACAG
ACTCTAGAGCCCCTACGGGAGCGGGTTGCCAGACAAGCAACTGCAGGCCCTGTCTCACCTGCTGCCT
CGGCTGCGGGCAGTGGCCGTGCAGTGTCCACACACAGACTCTGATTGATGGCCAACAGTCTTGGAG
GCATATGTGAGCCTTCGGGAGCTGGAGCAGCTGCGAGAGGATACGTGGGCACCCCTGGGGGGCTGGAG
TTGCCAGTCTTCCAGGGCTGGACCTTCTGTTGAGGCACTGGGCCAGGCTGTGGAAGCAGCTGCAGGG
GCCGCAGGAAGCTGGCAGGGAGACCAGCCCTGTTGGTGGCTGCTGTGCGTGTGGCGGAGGTGGAG
ACTGGACGAACAACCCCTGGGCCAGTCCCCGGGACTGGCGTCAGCGCTGTCTGAGGGCACTACAG
GAGGGCCTGGAGCAGGCCACTTTGGGTCACTCTGCTGCCTGCACCAGGGGCCCTACCAGGGTGGCTG
GAGGCTCTGCGAGTGGCCCTGCCAGTTCGAGTTGGCCACAGCTGAGGCACTAGTAGCGCCTTGTGCCCG
CCACAGTACAACGTGGTCCAGCTATGGGCCACACGCTGCATAGTGGTTTGCGCCGAGCCTGCAGAAC
CTCCTTGACAGGGCCTGAGCTAGAAGTGCAGTGCCTTCGCCCTTGTGCACTGGGCACTGCATGTGTAC
CTGGGGCAGGAAATGATGGGGAGCCTGGAGTTGGGGCTGAGGCTGATGTGTCCAGCTGGAGCCCTT
CTGACCTTGGAGAATGAGCAGCTGGAGGCAACATTTGTGGCCAACATCCAGGCAAGTGTGTCTCAG
TGGCTGCAGAAATGCACTGGATGGGAGGTAGCTGAGTGGGGCCGGGAGCATGGGCCAACACAGACCCG
TCTGGCTCCTATTACTACCAATGCCAGCATTGTGCTGCAGATCCTGGAAGAGAACATTGCTGTGGCC
AGCCTGGTCAGTGAAGTCACTGCAACAGCGAGTGCATGGCATGGCACTGTGCAACTGGGCACATTCTTG
AGGAGCTTCAAGTGTCTGATCCGATTCTCCCGAGACCACTTCAAGGGGAAATCAATGGCCCTCAT
TACGTGCCCTACCTACTGGCCGCCCTCAACCACAAGTCAAGCACTCAGCTCCTCAGTGTCTGTCTGCAG
CTGGACGGGGCGCCTTCAAGGGCCTTGGCTCCAGTGGAGCTGCCTGGACGAGTGCAGAGGAGGATC
TACCGCTTGGTGTGGAGGCGCTGCAGGCGGAGCTCCAGCCCTGTTCCGGGATCTGCCCTCGGCCAA
TGGCTGTGAGCCCTGAGCTCCTGCAAAGTGTGTGTAACGGACGGGGCGCTTCTGCCGGGACTTCTGG
CGCGTGCAGAACCCACGGTTCAGCTGTGCTGGCTGAGGCCGAGCGTGCCTGGTGTCTCCAGTACCTG
AGCGCGCTGATGCAAGGCCGCTGGTGTGCCGCGAGCCGACGAGAGGCCAGGCGGCCGAGCGCCTG
CGGCACGATGCTGCCAGCTTCAAGCAGCTTTTCTCAGTTTGGGCTGGAGGAGAACCGCACTGCAGC
CCGGTGTGCTCGCCCTGAGGGAGCTGTAACCTCCGCGACCCCGCGCTGCTGGGCTGGAGGTGGCT
GGCCTGCGGCAACAATTTCCGACGTGAGCGAGGACCAGTCTCCGCCCTTGGGCTGCGCGGGGAC
TTGTCCCGGGAGCAGCCTGGCCGCGCTCAGTCCCTGCAGGCTGCGCTGCCGCCGTGCCCCGCGCG
AGCCGCGCGTCTTTCAGCTAGTGCCCGCCCGCGCTGCGCCGGCCTCTGCTGCCCTCGGGG
TCCTGCGCCCGAGCCCTGCTGCTCGAGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:** □

**ACCN:** NM\_178516

**Insert Size:** 2241 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_178516.3</a>
<b>RefSeq Size:</b>	2505 bp
<b>RefSeq ORF:</b>	2241 bp
<b>Locus ID:</b>	283849
<b>UniProt ID:</b>	<a href="#">Q86V11</a>
<b>Cytogenetics:</b>	16q22.1
<b>MW:</b>	81.7 kDa
<b>Gene Summary:</b>	As part of the exocyst, may play a role in regulated exocytosis of insulin granules. [UniProtKB/Swiss-Prot Function]