

Product datasheet for **SC309314**

BAGE2 (NM_182482) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BAGE2 (NM_182482) Human Untagged Clone
Tag:	Tag Free
Symbol:	BAGE2
Synonyms:	CT2.2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_182482 edited
 ATGGTGGTGGCAACAGAGATGGCAGCGCAGCTGGAGTGTTAGGAGGGCGGCTGAGCGGT
 AGGAGTGGGGCTGGAGCAGTAAGATGGCGGCCGGAGCAGTAAGATGGCGGCTGGAGTGGT
 TTTTCTGGCATTGTCTGCCAGCTGCTCCAAGCCAGGCTGATGAAGGAGGAGTCCCCTGT
 GGTGAGCTGGAGTTGGAGCCTGAAGATGGCACAGCTCTCGATGTGCATTTTGTAAAGCAC
 TTTGGAGCCACTATCAAATGCTGTGAAGAGAAATGTACCCAGATGTATCATTATCCTTGT
 GCTGCAGGAGCCGACACCTTTCAGGATTTCAGTCACATCTTCTGCTTTGTCCAGAACAC
 ATTGACCAAGCTCCTGAAAGATCGAAGGAAGATGCAAACTGTGCAGTGTGCGACAGCCCG
 GGAGACCTCTTAGATCAGTTCTTTTGTACTACTTGTGGTCAGCACTATCATGGAATGTGC
 CTGGATATAGCGTTACTCCATTAACCGTGCAGGTTGGCAATGTCCTGAGTGCAAAGTG
 TGCCAGAACTGCAAACAATCGGGAGAAGATAGCAAGATGCTAGTGTGTGATACGTGTGAC
 AAAT

Restriction Sites: Please inquire

ACCN: NM_182482

Insert Size: 600 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: The open reading frame of this clone has been fully sequenced and found to be a perfect match to the protein associated with this reference, NM_182482.1.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_182482.1, NP_872288.1</u>
RefSeq Size:	1891 bp
RefSeq ORF:	330 bp
Locus ID:	85319
Cytogenetics:	21p11.2
Protein Families:	Druggable Genome