

Product datasheet for **SC309355**

ZFP200 (ZNF200) (NM_198088) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | ZFP200 (ZNF200) (NM_198088) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | ZFP200 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC309355 representing NM_198088. Blue=Insert sequence Red=Cloning site Green=Tag(s) |

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATGGCTGCAAAAGTGGTTCCTATGCCCAAGCCAAAGCAGTCCTTTATACTGAGAGTTCGCCA
GACTCCAAGCTGGGCCAAGACCTACTTCGAGATGCCACTAACGGGCCAAGACCTCCACCAGCTAGTG
CTGGAGCACTTCCCTCACCTTCTTGCCCAAGCCAAGCCTGGTCCAGCCCAGTCAGAAAGTCAAGGAGACC
TTGGTTATTATGAAAGATGTGAGCTCAAGCCTTCAGAACAGAGTGCATCCTCGTCCCTTGGTGAAGCTT
CTGCCCAAGGAGTCCAAAAGGAACAAGAGACAGTGTCTCTGTATTTGAAAGCTAACCCCTGAGGAGCTG
TGTTCTTTGAGGATTTGAATGTATTTCACTGCCAGGAAGAATGTGTGAGCTTGGATCCTACTCAACAA
CTCACGTGAGAGAAGGAAGATGACAGCAGTGTCCGGGAAATGATGTTACTGGCAGTCAATGGCAGTAAT
CCTGAAGGTGAAGATCCTGAGAGGGAACCTGTAGAAAATGAAGATTATAGAGAAAAGTCTTCAGATGAT
GATGAAATGGATTCTTCTTGGTCTCTCAGCAGCCTCCCGATAACCAGGAAAAGGAACGACTAAATACA
TCCATTCCACAAAAAGGAAAAATGAGAAATCTGTTAGTTACCATTGAGAAATGATACTCCTCTAGAGGAA
CTCTCAAATATGTAGACATCAGTATTATTGCCCTTACTCGAAATCGGAGGACAAGGAGATGGTACTACT
TGTCCTACTGTGTGGAAACAGTTAATGAAAGTCTTACCTCATTTCCACCAGAGGACCCACACTGGA
GAAAAACCCTATGACTGTAATCACTGTGGGAAAAGCTTCAATCATAAAACAAACCTCAATAAACATGAG
CGAATTCATACAGGAGAGAAACCTTATTCTGTCTCAGTGTGGAAAAAATTCCGTGAGAAATTCAT
CGGAGTCGTGATGAAGGAATCCATATAAGGGAGAAGATATTTAAGTGTCCAGAATGTGGGAAAACCTTC
CCAAAGAATGAGGAGTTTGTGCTTCATCTGCAGAGTCATGAGGCTGAGAGACCATATGGTTGCAAAAAA
TGTGGGAGAAGATTTGGTCGGCTGTCAAACGTACCCGGCATGAGAAAACCCACTCAGCCTGTAAGACC
CGAAAGCAGAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

| | |
|--------------------|-----------|
| Restriction Sites: | Sgfl-MluI |
| ACCN: | NM_198088 |



[View online »](#)

| | |
|-------------------------------|---|
| Insert Size: | 1188 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_198088.2</u> |
| RefSeq Size: | 3306 bp |
| RefSeq ORF: | 1188 bp |
| Locus ID: | 7752 |
| UniProt ID: | <u>P98182</u> |
| Cytogenetics: | 16p13.3 |
| Protein Families: | Transcription Factors |
| MW: | 45.5 kDa |
| Gene Summary: | <p>Could have a role in spermatogenesis.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to transcript variant 1. Both variants 1 and 2 encode the same isoform (1).</p> |