

Product datasheet for **SC309368**

MLX (NM_198205) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MLX (NM_198205) Human Untagged Clone
Tag:	Tag Free
Symbol:	MLX
Synonyms:	bHLHd13; MAD7; MXD7; TCFL4; TF4
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF sequence for NM_198205 edited ATGACGGAGCCGGGCGCCTCTCCCGAGGACCCTTGGGTCAAGGTGGAGTATGCCTACAGC GACAACAGCCTGGACCCCGATGATGAGGACAGTGATTACCACCAGGAGGCCTACAAGGAG TCCTACAAAGACCGGCGCGCGCGCACACACTCAGGCTGAGCAGAAGAGGAGGGACGCC ATCAAGAGAGGCTATGATGACCTTCAGACCATCGTCCCCACTTGCCAGCAGCAGGACTTC TCCATTGGCTCCCAAAAGCTCAGCAAAGCCATCGTTCTACAAAAGACCATTGACTACATT CAGTTTTTGCACAAGGAGAAGAAAAAGCAGGAGGAGGAGGTGCCACGTTACGCAAGGAT GTCACCGCCCTAAAGATCATGAAAGTGAAGTATGAGCAGATTGTGAAGGCACACCCGGAC AACCCCATGAAGGGGAGGACCAGGTCTCTGACCAGGTCAAGTTCAACGTGTTTCAAGGC ATCATGGATTCCCTGTTCCAGTCTTCAATGCCTCCATCTCAGTGGCCAGCTTCCAGGAG CTGTCAGCGTGTGTCTCAGCTGGATCGAGGAGCACTGTAAGCCTCAGACCCTGCGGGAG ATTGTGATTGGCGTCCTGCACCAATTGAAAAACCAGCTTTACTGA
Restriction Sites:	Please inquire
ACCN:	NM_198205
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The ORF of this clone has been fully sequenced and found to contain one SNP compared with NM_198205.1.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_198205.1](#), [NP_937848.1](#)

RefSeq Size: 2316 bp

RefSeq ORF: 645 bp

Locus ID: 6945

UniProt ID: [Q9UH92](#)

Cytogenetics: 17q21.2

Protein Families: Druggable Genome, Transcription Factors

Gene Summary: The product of this gene belongs to the family of basic helix-loop-helix leucine zipper (bHLH-Zip) transcription factors. These factors form heterodimers with Mad proteins and play a role in proliferation, determination and differentiation. This gene product may act to diversify Mad family function by its restricted association with a subset of the Mad family of transcriptional repressors, namely, Mad1 and Mad4. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (1) encodes the shortest isoform (alpha).