

Product datasheet for **SC309489**

Dystrophin (DMD) (NM_004022) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dystrophin (DMD) (NM_004022) Human Untagged Clone
Tag: Tag Free
Symbol: DMD
Synonyms: BMD; CMD3B; DXS142; DXS164; DXS206; DXS230; DXS239; DXS268; DXS269; DXS270; DXS272; MRX85

Vector: pCMV6 series

Fully Sequenced ORF: >NCBI ORF sequence for NM_004022, the custom clone sequence may differ by one or more nucleotides

```

ATGCCATCTTCCTTGATGTTGGAGGTACCTGCTCTGCCAGATTTCAACCGGGCTTGACA
GAACTTACCGACTGGCTTCTCTGCTTGATCAAGTTATAAAATCACAGAGGGTGATGGTG
GGTGACCTTGAGGATATCAACGAGATGATCATCAAGCAGAAGGCAACAATGCAGGATTTG
GAACAGAGGCGTCCCAGTTGGAAGAACTATTACCGCTGCCAAAATTTGAAAAACAAG
ACCAGCAATCAAGAGGCTAGAACAATCATTACGGATCGAATTGAAAGAATTCAGAATCAG
TGGGATGAAGTACAAGAACCTTCAGAACCGGAGGCAACAGTTGAATGAAATGTTAAAG
GATTCAACACAATGGCTGGAAGCTAAGGAAGAAGCTGAGCAGGTCTTAGGACAGGCCAGA
GCCAAGCTTGAGTCATGGAAGGAGGTCCTATACAGTAGATGCAATCCAAAAGAAAATC
ACAGAAAACAAGCAGTTGGCCAAAGACCTCCGCCAGTGGCAGACAAATGTAGATGTGGCA
AATGACTTGGCCCTGAAACTTCTCCGGGATTATTCTGCAGATGATACCAGAAAAGTCCAC
ATGATAACAGAGAATATCAATGCCTCTTGAGAAGCATTCAAAAAGGGTGAGTGAGCGA
GAGGCTGCTTTGGAAGAACTCATAGATTACTGCAACAGTTCCCCTGGACCTGGAAAAG
TTTCTTGCTGGCTTACAGAAGCTGAAACAACCTGCCAATGTCCTACAGGATGCTACCCGT
AAGGAAAGGCTCCTAGAAGACTCCAAGGGAGTAAAAGAGCTGATGAAACAATGGCAAGAC
CTCCAAGGTGAAATTGAAGCTCACACAGATGTTTATCACAACTGGATGAAAACAGCCAA
AAAATCCTGAGATCCCTGGAAGGTTCCGATGATGCAGTCTGTTACAAAGACGTTTGGAT
AACATGAACTTCAAGTGGAGTGAACCTCGGAAAAAGTCTCTCAACATTAGGTCCCATTG
GAAGCCAGTTCTGACCAGTGAAGCGTCTGCACCTTTCTCTGCAGGAACTTCTGGTGTGG
CTACAGCTGAAAGATGATGAATTAAGCCGGCAGGCACCTATTGGAGGCGACTTTCAGCA
GTTCAGAAGCAGAACGATGTACATAGGGCCTTCAAGAGGGAATTGAAAACCTAAAGAACCT
GTAATCATGAGTACTCTTGAGACTGTACGAATATTTCTGCAGAGCAGCCTTTGGAAGGA
CTAGAGAACTCTACCAGGAGCCAGAGAGCTGCCTCCTGAGGAGAGAGCCAGAATGTC
ACTCGGCTTCTACGAAAGCAGGCTGAGGAGGTCAATACTGAGTGGGAAAAATTGAACCTG
CACTCCGCTGACTGGCAGAGAAAAATAGATGAGACCCTTGAAGACTCCAGGAACCTCAA
GAGGCCACGGATGAGCTGGACCTCAAGCTGCGCCAAGCTGAGGTGATCAAGGGATCCTGG
CAGCCCGTGGGGATCTCCTCATTGACTCTCTCAAGATCACCTCGAGAAAGTCAAGGCA
CTTCGAGGAGAAAATTGCGCCTCTGAAAGAGAAGCTGAGCCACGTCAATGACCTTGCTCGC
CAGCTTACCACCTTGGGCATTGAGTCTCACCGTATAACCTCAGCACTCTGGAAGACCTG

```



[View online »](#)

```

AACACCAGATGGAAGCTTCTGCAGGTGGCCGTCGAGGACCGAGTCAGGCAGCTGCATGAA
GCCACAGGGACTTTGGTCCAGCATCTCAGCACTTTCTTCCACGTCTGTCCAGGGTCCC
TGGGAGAGAGCCATCTCGCCAAACAAAGTGCCCTACTATATCAACCACGAGACTCAAACA
ACTTGCTGGGACCATCCCAAAATGACAGAGCTCTACCAGTCTTTAGCTGACCTGAATAAT
GTCAGATTCTCAGCTTATAGGACTGCCATGAAACTCCGAAGACTGCAGAAGGCCCTTTCG
TTGGATCTCTTGAGCCTGTCAGCTGCATGTGATGCCTTGGACCAGCACACCTCAAGCAA
AATGACCAGCCCATGGATATCCCTGCAGATTATTAATTGTTTGACCACTATTTATGACCCG
CTGGAGCAAGAGCACAAACAATTTGGTCAACGTCCCTCTCTGCGTGGATATGTGTCTGAAC
TGCTGTGAATGTTTATGATACGGGACGAACAGGGAGGATCCGTGTCCTGTCTTTTAAA
ACTGGCATCATTTCCCTGTGTAAAGCACATTTGGAAGACAAGTACAGATACCTTTTCAAG
CAAGTGGCAAGTTCAACAGGATTTTGTGACCAGCGCAGGCTGGGCCTCCTTCTGCATGAT
TCTATCCAAATCCAAGACAGTTGGGTGAAGTTGCATCCTTTGGGGCAGTAACATTGAG
CCAAGTGTCCGGAGCTGCTTCCAATTTGCTAATAATAAGCCAGAGATCGAAGCGGCCCTC
TTCCTAGACTGGATGAGACTGGAACCCAGTCCATGGTGTGGCTGCCCGTCTGCACAGA
GTGGCTGTGCAGAACTGCCAAGCATCAGGCCAAATGTAACATCTGCAAAGAGTGTCCA
ATCATTGGATTACAGTACAGGAGTCTAAAGCACTTTAATTATGACATCTGCCAAAGCTGC
TTTTTTCTGGTCGAGTTGCAAAAGGCCATAAAATGCACTATCCCATGGTGGAAATATTGC
ACTCCGACTACATCAGGAGAAGATGTTTCGAGACTTTGCCAAGGTAATAAAAAACAATTT
CGAACCAAAAGGTATTTTGCAGCATCCCCGAATGGGCTACCTGCCAGTGCAGACTGTC
TTAGAGGGGGACAACATGGAAACGCCTGCCTCGTCCCCTCAGCTTTCACACGATGATACT
CATTACGCATTGAACATTATGCTAGCAGGCTAGCAGAAATGGAAAACAGCAATGGATCT
TATCTAAATGATAGCATCTCTCCTAATGAGAGCATAGATGATGAACATTTGTTAATCCAG
CATTACTGCCAAAGTTTGAACCAGGACTCCCCCTGAGCCAGCCTCGTAGTCTGCCCCAG
ATCTTGATTTCCCTTAGAGAGTGAGGAAAAGAGGGGAGCTAGAGAGAATCCTAGCAGATCTT
GAGGAAGAAAACAGGAATCTGCAAGCAGAATATGACCGTCTAAAGCAGCAGCACGAACAT
AAAGGCCTGTCCCACTGCCGTCCCCTCCTGAAATGATGCCACCTCTCCCAGAGTCCC
CGGGATGCTGAGCTCATTGCTGAGGCCAAGCTACTGCGTCAACACAAAGGCCGCTGGAA
GCCAGGATGCAATCCTGGAAGACCACAATAAACAGCTGGAGTACAGTTACACAGGCTA
AGGCAGCTGCTGGAGCAACCCAGGCAGAGGCCAAAGTGAATGGCACAACGGTGTCTCT
CCTTCTACCTCTCTACAGAGTCCGACAGCAGTACGCCTATGCTGCTCCGAGTGGTTGGC
AGTCAAATTCGGACTCCATGGGTGAGGAAGATCTTCTCAGTCTCCCCAGGACACAAGC
ACAGGGTTAGAGGAGGTGATGGAGCAACTCAACAACCTCTCCCTAGTTC AAGAGGACAC
AATGTAGGAAGTCTTTCCACATGGCAGATGATTTGGGCAGAGCGATGGAGTCTTAGTA
TCAGTCATGACAGATGAAGAAGGAGCAGAATAA
    
```

- Restriction Sites:** Please inquire
- ACCN:** NM_004022
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004022.1](#), [NP_004013.1](#)

RefSeq Size: 7339 bp

RefSeq ORF: 3693 bp

Locus ID: 1756

UniProt ID: [P11532](#)

Cytogenetics: Xp21.2-p21.1

Protein Pathways: Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM), Viral myocarditis

Gene Summary: This gene spans a genomic range of greater than 2 Mb and encodes a large protein containing an N-terminal actin-binding domain and multiple spectrin repeats. The encoded protein forms a component of the dystrophin-glycoprotein complex (DGC), which bridges the inner cytoskeleton and the extracellular matrix. Deletions, duplications, and point mutations at this gene locus may cause Duchenne muscular dystrophy (DMD), Becker muscular dystrophy (BMD), or cardiomyopathy. Alternative promoter usage and alternative splicing result in numerous distinct transcript variants and protein isoforms for this gene. [provided by RefSeq, Dec 2016]

Transcript Variant: Dp140 transcripts use exons 45-79, starting at a promoter/exon 1 located in intron 44. Dp140 transcripts have a long (1 kb) 5' UTR since translation is initiated in exon 51 (corresponding to aa 2461 of dystrophin). In addition to the alternative promoter and exon 1, differential splicing of exons 71-74 and 78 produces at least five Dp140 isoforms. Of these, this transcript (Dp140ab) lacks exons 71 and 78 and encodes a protein with a unique C-terminus.