

Product datasheet for **SC309557**

RHCE (NM_138617) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RHCE (NM_138617) Human Untagged Clone
Tag:	Tag Free
Symbol:	RHCE
Synonyms:	CD240CE; RH; Rh4; RH30A; RHC; RHCe(152N); RHE; RhIVb(J); RHIXB; RHNA; RHPI; RhVI; RhVIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC309557 representing NM_138617. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGAGCTCTAAGTACCCGCGGTCTGTCCGGCGCTGCCTGCCCTCTGGGCCCTAACACTGGAAGCAGCT
CTCATTCTCCTCTTCTATTTTTTACCCTATGACGCTTCCTTAGAGGATCAAAGGGGCTCGTGGCA
TCCTATCAAGTCGGCCAAGATCTGACCGTGATGGCGGCCCTTGGCTTGGGCTTCTCACCTCAAATTC
CGGAGACACAGCTGGAGCAGTGTGGCCTCAACCTTTCATGCTGGCGCTTGGTGTGCAGTGGCAATC
CTGCTGGACGGCTTCTGAGCCAGTTCCTCCTGGGAAGGTGGTCATCACACTGTTCAAGTATTGGCTG
GCCACCATGAGTCTATGTCGGTCTGATCTCAGCGGGTGTCTTGGGAAGGTCAACTGGCGCAG
TTGGTGGTGTGGTGTGGTGGAGGTGACAGCTTAGGCACCCTGAGGATGGTCATCAGTAATATCTTC
AACACTTATGTGCACAGTGCAGTGTGGTGTGGCAGGAGCGTGGCTGTGGGTACCTCGTGTACCTGATCCCT
TCTCCGTGGCTTGCCATGGTGTGGTCTTGTGGCTGGGCTGATCTCCATCGGGGAGCCAAGTGCCCTG
CCGGTGTGTTGTAACCGAGTGTGGGATTACCACATCTCCGTGCATGCACTCCATCTTCAGCTTGTG
GGTCTGCTGGAGAGATCACCTACATTGTGCTGCTGGTGTTCATACTGTCTGGAACGGCAATGGCATG
TTTGCTCCTAAATCTCAAATATGGAAGCACCTCATGTGGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites:	Sgfl-MluI
Plasmid Map:	□
ACCN:	NM_138617
Insert Size:	804 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_138617.3
RefSeq Size:	1223 bp
RefSeq ORF:	804 bp
Locus ID:	6006
UniProt ID:	P18577
Cytogenetics:	1p36.11
Protein Families:	Transmembrane
MW:	28.7 kDa
Gene Summary:	<p>The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene which encodes both the RhC and RhE antigens on a single polypeptide and a second gene which encodes the RhD protein. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes. A mutation in this gene results in amorph-type Rh-null disease. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms. [provided by RefSeq, Aug 2016]</p> <p>Transcript Variant: This variant (4), also called RhVI, lacks three internal exons resulting in a frameshift and use of an upstream stop codon, as compared to variant 1. Isoform 4 lacks an internal segment of 45 aa and has a unique C-terminus as compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>