

Product datasheet for SC309562

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ST3GAL3 (NM_174967) Human Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: ST3GAL3 (NM_174967) Human Untagged Clone

Tag: Tag Free Symbol: ST3GAL3

Synonyms: DEE15; EIEE15; MRT12; SIAT6; ST3GALII; ST3Gal III; ST3GalIII; ST3N

Vector: pCMV6 series

Fully Sequenced ORF: >NCBI ORF sequence for NM_174967, the custom clone sequence may differ by one or more

nucleotides

ATGGGACTCTTGGTATTTGTGCGCAATCTGCTGCTAGCCCTCTGCCTCTTTCTGGTACTG
GGATTTTTGTATTATTCTGCGTGGAAGCTACACTTACTCCAGTGGAGGAGGACTCCAAT
TCAGTGGTTCTTTCCTTTGACTCCGCTGGACAAACACTAGGCTCAGAGTATGATCGGTTG
GGCTTCCTCGAATCTGGACTCTAAACTGCCTGCTGAATTAGCCACCAAGTACGCAAAC
TTTTCAGAGGGAGCTTGCAAGCCTGGCTATGCTTCAGCCTTGATGACGGCCATCTTCCCC
CGGTTCTCCAAGCCAGCACCCATGTTCCTGGATGACTCCTTTCGCAAGTGGGCTAGAATC
CGGGAGTTCGTGCCGCCTTTTGGGATCAAAGGTCAAGTCCTGGACGCACAATATCCAGCG

AGAGAAAGAGTTTCTGCGGAAGCTGGTGAAAGCTCGCGTCATCACTGA

Restriction Sites: Please inquire **ACCN:** NM 174967

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a

point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative

RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning

into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).





Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM 174967.1, NP 777627.1

RefSeq Size: 1630 bp RefSeq ORF: 468 bp Locus ID: 6487 **UniProt ID:** Q11203 Cytogenetics: 1p34.1

Protein Families: Secreted Protein, Transmembrane

Protein Pathways: Glycosphingolipid biosynthesis - lacto and neolacto series, Keratan sulfate biosynthesis,

Metabolic pathways

Gene Summary: The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of

sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein is normally found in the Golgi apparatus but can be proteolytically processed to a soluble form. This protein is a member of glycosyltransferase family 29. Mutations in this gene have been associated with a form of autosomal recessive nonsymdromic cognitive disability as well as infantile epileptic encephalopathy. Multiple transcript variants encoding several different

isoforms have been found for this gene. [provided by RefSeq, Jul 2017]

Transcript Variant: This variant (5) has multiple differences in the coding region, one of which results in a frameshift compared to variant 1. The resulting isoform (e, also called B7) has a distinct C-terminus and is shorter compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.