

Product datasheet for **SC309700**

RASGRP 4 (RASGRP4) (NM_170604) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RASGRP 4 (RASGRP4) (NM_170604) Human Untagged Clone
Tag:	Tag Free
Symbol:	RASGRP4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC309700 representing NM_170604.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAACAGAAAAGACAGTAAGAGGAAGTCCCACCAGGAATGCACCGAAAAATAGGAGGGCGAGGCCGG
CCCCGCCAAGTGCGCCGCCACAAGACATGCCCCAGCCCTCGGAAATCAGCAAGGTCATGGCTTCCATG
AACCTGGGCTGCTGAGTGAGGGCGGCTGCAGCGAAGATGAGCTGCTGGAGAAATGCATCCAGTCCCTC
GATTCAGCTGGCAGCCTGTGCCACGAGGACCACATGCTCAACATGGTGTGGCCATGCACAGCTGGGTG
CTGCCGTCGCCGACCTGGCTGCCCGCTGCTGACCTCATACCAGAAGGCCACAGGGGACACCCAGGAG
CTGAGACGGCTGCAGATCTGTACCTGGTCAGGTACTGGCTGATGCGACACCCTGAGGTGATGCACCAG
GATCCCCAGCTAGAAGAAGTCAAGGTCGTTTCTGGGCCACCGTGGCCCGGAGGGCAACTCAGCCAG
AGAAGACTGGGAGACTCTTCTGACCTCCTGAGCCCTGGTGGCCCTGGCCCCACTCCCAATGAGCAGC
CCAGGCCCTGGGCAAAAAGCGCAAAGTGCCTTGTCTTTTCGACCACTTGGAGACGGGGAGCTGGCTCAG
CACCTCACCTACCTGGAGTTCGGTCTCCAGGCTATCACGCCACAGGACTGCGGAGCTACGTTTTG
CAGGGCTCAGTACGAGGCTGCCCGCCCTGGAGGGCTCCGTAGGTCTCAGCAACAGCGTGTCCCGCTGG
GTGCAGGTGATGGTGTGAGCCGTCGCCGGCCCTACAGCGTGCACAGGTGCTGGACAAGTTCATTAC
GTGGCACAGAGGCTCCACCAGCTGCAGAAATTTCAACACGCTGATGGCAGTACAGGGGGCCTGTGTAC
AGTGCCATCTCCAGACTCAAGGACTCCCATGCCACCTGAGCCCTGACAGCACCAAGGCCCTCCTGGAG
CTCACTGAGTCTTGCCTCCACAACAACTACGCCCGTACCGCCGACCTGGGCTGGCTGCGCGGGT
TTCCGGCTGCCTGACTGGGCGTGACCTCAAGGACCTGGTGTCCCTGCATGAGGCACAGCCGACAGG
TTGCCGTACGGCCGCTGCACCTACCCAAGCTGAACAACCTCTACCTGCGGCTGCAGGAGCTGTGGCC
CTCCAAGGGCAGCATCCACCCTGCAGCGCAATGAGGATCTGTCACCTGCTCACGCTCTCCCTGGAC
CTTTCTACACGGAAGACGAGATCTATGAGCTTTCTTATGCCCGGAGCCGCTTGTCCCAAGAGCCTG
CCACCCTCCCCCTCAATGCACCTCTGGTGGTGGAGTGGGCCCTGGTGTGACACCAAGCCGGACAGG
GTCACACTGGGTGCGCATGTGGAGCAGCTGGTGGAGTCTGTGTTCAAGAATTATGACCCTGAAGGCCGA
GGAACAATCTCTCAGGAGGACTTTGAGCGACTCTCGGCAATTTTCCCTTCGCTGCCATGGGCTTAC
CCACCCACGCCAGGGGAGAGGATCCTTCAGCAGAGAGGAGCTGACAGGGTACCTGCTCCGGGCCAGC
GCCATCTGCTCCAAGTTGGGCTGGCTTCTGACACCTTCCATGAGGTACCTCCGAAAGCCTACC
TTCTGCGACAGCTGCAGTGGCTTCTCTGGGTGTACCAAGCAAGGCTACCGCTGTCGGGAGTGGGG
CTGTGTTGCCACAAACTGCAGAGACCAGGTGAAGGTAGAATGTAAGAAGAGGCCAGGGCCAAGGGC
GATGCAGGACCCCGGAGCTCCTGTCCATCCACACCAGCTCCCATGCCAGCTGTGGCTCCGAGGAA
AATCACTCTACACGCTATCCCTGGAGCCTGAGACTGGGTGCCAGCTTCGCCATGCCTGGACCCAGACT
GAATCCCCACACCCTTCTGGGAAACAGATACGGTCCCCTGCCCGGTGATGGACCCACCATCAACTGCA
TCCTCCAAGCTGGATTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_170604

Insert Size: 2022 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_170604.2
RefSeq Size:	3235 bp
RefSeq ORF:	2022 bp
Locus ID:	115727
UniProt ID:	Q8TDF6
Cytogenetics:	19q13.2
Protein Families:	Druggable Genome
Protein Pathways:	MAPK signaling pathway
MW:	74.9 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the Ras guanyl nucleotide-releasing protein (RasGRP) family of Ras guanine nucleotide exchange factors. It contains a Ras exchange motif, a diacylglycerol-binding domain, and two calcium-binding EF hands. This protein was shown to activate H-Ras in a cation-dependent manner in vitro. Expression of this protein in myeloid cell lines was found to be correlated with elevated level of activated RAS protein, and the RAS activation can be greatly enhanced by phorbol ester treatment, which suggested a role of this protein in diacylglycerol regulated cell signaling pathways. Studies of a mast cell leukemia cell line expressing substantial amounts of abnormal transcripts of this gene indicated that this gene may play an important role in the final stages of mast cell development. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2009]</p> <p>Transcript Variant: This variant (a) represents the longest transcript and encodes the longest isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.</p>