

Product datasheet for **SC309723**

CREB5 (NM_182898) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CREB5 (NM_182898) Human Untagged Clone
Tag:	Tag Free
Symbol:	CREB5
Synonyms:	CRE-BPA; CREB-5; CREBPA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC309723 representing NM_182898.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGATTTATGAGGAATCCAAGATGAATTTGGAGCAGGAGAGGCCGTTTGTCTGCAGTGCCCCAGGCTGC
TCCAGCGCTTCCAACAGAGGACCATCTGATGATTCATAGGCACAAACATGAAATGACTTTGAAGTTT
CCTTCAATAAAAAACAGACAATATGTTATCAGATCAAACCTCCGACCCCAACGAGATTCCTGAAGAAGTGC
GAGGAGGTGGGCCTCTTACGCGAGCTGGACTGCTCCCTGGAGCACGAGTTGAGGAAGGCTCAGGAAGAG
GAGAGCAGCAAGCGGAATATCTCGATGCATAATGCAGTTGGTGGGGCCATGACGGGGCCCGAACTCAC
CAGCTTAGCAGCGCTCGGCTGCCAACCATGACACCAACGTTGTGATTGAGCAAGCCATGCCGTCGCT
CAGTCCAGCTCTGTCATCACTCAGGCACCTTCCACCAACCGCCAGATCGGGCCTGTCCAGGCTCTCTA
TCTTCTGCTACATCTCCACAACAGACAGAGACGCCCATGCCAGCCTCCATGCCTGGGACCCTGCC
AACCTACAATGCCAGGATCTTCCGCCGCTTGTGATGCCAAATGGAGCGACAATGTCAGTGAAGTCCAGC
ATCATGGGGATGCAAGGTCCAAATCTCAGCAACCCTGTGCTTCTCCCCAGGTCCAGCCAATGCATTCA
GAAGCCAAAATGAGGTTGAAGGCTGCATTGACTACCACCTGTGCCATGTCAAATGGGAACATGAAC
ACCATGGGACACATGATGGAGATGATGGGCTCCCGGCAGGACCAGACGCCACACCATCACATGCACTCG
CACCCGCATCAGCACCAGACACTGCCACCCATCACCTTACCCACACCAGCACCAGCACCAGCACAC
CATCCTCACCTCAACCCATCACCCAGCAGAACCATCCACATCACCACTCCCATTTCCACCTTATGCA
CACCCAGCACATCACCAGACCTCGCCACATCCGCCCTGCACACCCGCAACCAAGCACAGGTTTACCA
GCAACACAACAGATGCAGCCAACCCAGACAATACAGCCACCCAGCCACAGGGGGCGCCGGCGAAGG
GTGGTAGACGAGGATCCGGACGAGAGGGCGGGAATTTCTGAAACGGAACCGGGCAGCTGCCACCCGC
TGCAGACAGAAGAGGAAGTCTGGGTGATGTCATTGAAAAGAAAGCAGAAGAAGTCAACCCAGACAAC
ATGCAGCTTCAGAATGAAGTGTCTATGTTGAAAAATGAGGTGGCCAGCTGAAACAGTTGTTGTTAACA
CATAAAGACTGCCAATAACAGCCATGCAGAAAAGAAATCAAGGATATCTAAGTCCAGAGAGTAGCCCT
CCTGCTAGTCTGTCCCAGCTTGTCCCAGCAACAAGTATCCAGCATAATACCATCACTACTTCTCTCA
TCGGTCAGCGAGGTGGTAGGAAGCTCCACCCTCAGCCAGCTCACCACTCACAGAACAGACCTGAATCCG
ATTCTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:

ACCN: NM_182898

Insert Size: 1527 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_182898.3
RefSeq Size:	8554 bp
RefSeq ORF:	1527 bp
Locus ID:	9586
UniProt ID:	Q02930
Cytogenetics:	7p15.1
Protein Families:	Transcription Factors
Protein Pathways:	Huntington's disease, Prostate cancer
MW:	56.9 kDa
Gene Summary:	<p>The product of this gene belongs to the CRE (cAMP response element)-binding protein family. Members of this family contain zinc-finger and bZIP DNA-binding domains. The encoded protein specifically binds to CRE as a homodimer or a heterodimer with c-Jun or CRE-BP1, and functions as a CRE-dependent trans-activator. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript, and encodes the longest isoform (alpha). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>