

Product datasheet for **SC309724**

DRP1 (DNM1L) (NM_005690) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DRP1 (DNM1L) (NM_005690) Human Untagged Clone
Tag:	Tag Free
Symbol:	DNM1L
Synonyms:	DLP1; DRP1; DVLP; DYMPLE; EMPF; EMPF1; HDYNIV; OPA5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC309724 representing NM_005690.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGGCGCTAATTCTGTGCATAAACAAGCTCCAGGACGTCTTCAACACGGTGGGCGCCGACATCATC
CAGCTGCCTCAAATCGTCGTAGTGGGAACGCAGAGCAGCGAAAGAGCTCAGTGCTAGAAAGCCTGGTG
GGGAGGGACCTGCTTCCAGAGGTACTGGAATTGTCACCCGGAGACCTCTCATTCTGCAACTGGTCCAT
GTTTCACAAGAAGATAAACCGAAAAACAAGGAGAAGAAAATGGGGTGAAGCAGAAGAATGGGGTAAA
TTTCTTACACCAAAAAAAGCTTTACACGGATTTTGTGAAATTCGACAAGAAATGAAAATGAAACA
GAAAGAATTTAGGAAATAATAAGGGAGTAAGCCCTGAACCAATTCATCTTAAGATTTTTTACCCCAAC
GTTGTCAATTTGACACTTGTGGATTTGCCAGGAATGACCAAGGTGCCTGTAGGTGATCAACCTAAGGAT
ATTGAGCTTCAAATCAGAGAGCTCATTCTCGGTTTCATCAGTAATCCTAATTCATTATCCTCGCTGTC
ACTGCTGCTAATACAGATATGGCAACATCAGAGGCACTTAAAATTTCAAGAGAGGTAGATCCAGATGGT
CGCAGAACCCTAGCTGTAATCACTAAACTTGATCTCATGGATGCGGGTACTGATGCCATGGATGTATTG
ATGGGAAGGGTTATTCCAGTCAAACCTTGAATAATTGGAGTAGTTAACAGGAGCCAGCTAGATATTAAC
AACAGAAGAGTGAACCTGATTCATCCGTGATGAGTATGCTTTTCTTAAAAGAAATATCCATCTCTG
GCCAATAGAAATGGAACAAAGTATCTTGCTAGGACTCTAAACAGGTTACTGATGCATCACATCAGAGAT
TGTTTACCAGAGTTGAAAACAAGAATAAATGTTCTAGCTGCTCAGTATCAGTCTCTTCTAAATAGCTAC
GGTGAACCCGTGGATGATAAAAGTGCTACTTTACTCCAATTATTACCAATTTGCCACAGAATATTGT
AACACTATTGAAGGAAGTCAAAATATATTGAACTTCGGAGCTATGCGGTGGTCTAGAATTTGTTAT
ATTTTCCATGAGACTTTTGGGCGAACCTTAGAATCTGTTGATCCACTTGGTGGCCTAACACTATTGAC
ATTTTACTGCCATTAGAAATGCTACTGGTCTCGTCTGCTTTATTTGTCCTGAGGTTTCATTTGAG
TTACTGGTGAAGCGGCAAAATCAAACGTCTAGAAGAGCCAGCCCTCGCTGTGTGGAAGTGGTTCATGAG
GAAATGCAAAGGATCATTGAGCACTGTAGCAATTACAGTACACAGGAATTGTTACGATTTCTAAACTT
CATGATGCCATAGTTGAAGTGGTACTTGTCTTCTCGTAAAAGTTGCCTGTTACAAATGAAATGGTC
CATAACTTAGTGGCAATTGAACTGGCTTATATCAACACAAAACATCCAGACTTTGCTGATGCTTGTGGG
CTAATGAACAATAATAGAGGAACAAAGGAGAAACAGGCTAGCCAGAGAATTACCTTCAGCTGTATCA
CGAGACAAGTTGCATCTGGAGGTGGTGGGTTGGAGATGGTGTCAAGAACCAACCACAGGCAACTGG
AGAGGAATGCTGAAAACCTCAAAAGCTGAAGAGTTATTAGCAGAAGAAAAATCAAAACCCATTCCAATT
ATGCCAGCCAGTCCAAAAAGGTCATGCCGTGAACCTGCTAGATGTGCCAGTTCCTGTTGCACGAAAA
CTATCTGCTCGGGAACAGCGAGATTGTGAGGTTATTGAACGACTCATTAAATCATATTTTCTCATTGTC
AGAAAGAATATTCAGACAGTGTGCCAAAGGCAGTAATGCATTTTTTGGTTAATCATGTGAAGACT
CTTCAGAGTGAGCTAGTAGGCCAGCTGTATAAATCATCTTATTGGATGATCTTCTGACAGAATCTGAG
GACATGGCACAGCGCAGGAAAGAAGCAGCTGATATGCTAAAGGCATTACAAGGAGCCAGTCAAATATT
GCTGAAATCCGGGAGACTCATCTTTGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:

ACCN: NM_005690

Insert Size: 2100 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005690.4
RefSeq Size:	4504 bp
RefSeq ORF:	2100 bp
Locus ID:	10059
UniProt ID:	O00429
Cytogenetics:	12p11.21
Domains:	dynamain_2, dynamin, GED
Protein Pathways:	Endocytosis, Fc gamma R-mediated phagocytosis
MW:	78.1 kDa
Gene Summary:	<p>This gene encodes a member of the dynamin superfamily of GTPases. The encoded protein mediates mitochondrial and peroxisomal division, and is involved in developmentally regulated apoptosis and programmed necrosis. Dysfunction of this gene is implicated in several neurological disorders, including Alzheimer's disease. Mutations in this gene are associated with the autosomal dominant disorder, encephalopathy, lethal, due to defective mitochondrial and peroxisomal fission (EMPF). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jun 2013]</p> <p>Transcript Variant: This variant (3) lacks two consecutive in-frame exons in the 3' coding region, compared to variant 1. The encoded isoform (3) is shorter, compared to isoform 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>