

Product datasheet for **SC309784**

ABHD15 (NM_198147) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ABHD15 (NM_198147) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABHD15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC309784 representing NM_198147. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCCGCTGGGGCGCCGCTCGCGCTCATCTTGGCCGTGCTCGCCCTTCTCGGCCTGCTCGGCCCG
CGGCTCCGGGGACCTGGGGCGCGCCGTGGAGAGAGGACCCTGCCGGGGGCCAAGACCGAGACGAC
GGGAGGAGGCGGACGGCGGAGGCCCGCGGACCAGTTTCAGCGACGGGCGGAGCCACTGCCGGGAGGG
TGCAGCCTTGTTCGAAGCCGTCGGCCCTGGCCAGTGCCTGCTGCGCGCCCTGCCGGCCTCAGAGGCG
CTGGAGGCCGGCCGCGCTCTGGTTCTCCGGGCCACCTGCAGACCCTCTGCCACTTCGTCCTGCC
GTAGCCCTGGCCTGAGCTGGCCGGGAGTACCTGCAGTTGGCGGACGATGGGCTAGTGCCCTGGAC
TGGGTGGTAGGACCTTGTGTTCCGGGGCCCGGATCACCAGCGCCGGGGCCTTCCTGCGGTGCTTCTG
GTGATCCCAATGCGTGGGTGCGCTCACCCGCAACGTGCTCGGCCTTTGCTTGCCTGCGCCCTGGAGCGC
GGCTACTACCCGGTCATCTTCCATCGCCGCGGCCACCACGGTTGCCCACTGGTCAGCCCCGGCTGCAG
CCTTTCGGGGACCCGTCCGACCTCAAGGAGGCGGTACATACATCCGCTTCCGACACCCGGCGGCCG
CTGTTCCGGGTGAGCGAAGGCTCGGGCTCGGCGTGTCTCTGCTACCTGGCGAGTGGGCTCCTCC
AGCTACGTGACAGGCGCCGCTGCATCTCGCCGTGCTGCGTGCAGAGTGGTTCGAGGCCGGCCTG
CCCTGGCCCTACGAGCGGGCTTTCTGCTCCACCAGAAGATCGCCCTCAGCAGGTATGCCACAGCCCTG
GAGGACTGTGGACACCAGCAGACTGTTCAAGGAGCCGTTCCCTTCGAGAGTTTGAGGAGGCTCTCTTC
TGCCACACAAAAGCTTCCCATCAGCTGGGATACCTACTGGGACCGCAACGACCCGCTCCGGGATGTC
GATGAGGACGCCGTGCCTGTGCTGTATCTGCAGTGCTGACGACCCCGTGTGTGGACCCCGAGACCAC
ACTCTGACAACTGAACCTTCCACAGCAACCCCTACTTCTTCTCCTGCTCAGTCGCCACGGAGGCCAC
TGTGGCTTCTGCGCCAGGAGCCCTTGCCAGCCTGGAGCCATGAGGTATCTTGGAGTCTTCCGGGCC
TTGACTGAGTCTTCCGAACGGAGGAGAGGATTAAGGGCTGAGCAGGCACAGAGCTTCTTCTTGGG
GGCCGTCGTGGGGGAGCCTTGACAGAGCGGGAAGTCTTCTCTTCCAACCTGGAGGAGATCTTT
AACTGGAAGCGATCATAACAAGTGA
ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
Plasmid Map:	□
ACCN:	NM_198147
Insert Size:	1407 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198147.2</u>
RefSeq Size:	3365 bp
RefSeq ORF:	1407 bp
Locus ID:	116236
UniProt ID:	<u>Q6UXT9</u>
Cytogenetics:	17q11.2
Protein Families:	Transmembrane
MW:	51.8 kDa