

## Product datasheet for **SC309785**

### **C18orf34 (CCDC178) (NM\_198995) Human Untagged Clone**

#### **Product data:**

|                           |                                                     |
|---------------------------|-----------------------------------------------------|
| Product Type:             | Expression Plasmids                                 |
| Product Name:             | C18orf34 (CCDC178) (NM_198995) Human Untagged Clone |
| Tag:                      | Tag Free                                            |
| Symbol:                   | C18orf34                                            |
| Synonyms:                 | C18orf34                                            |
| Mammalian Cell Selection: | Neomycin                                            |
| Vector:                   | pCMV6-Entry (PS100001)                              |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                |



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_198995, the custom clone sequence may differ by one or more nucleotides

```

ATGACTGAAAACAAGACAGTTTCCTCTTCTTCCACTAGAGATGATCAAACCAATATAGGT
TTAACATGTCAGGAAGTAAAGGCTCTCAGAGAGAAGGCATGGTCAAGGACAAATGAAGGC
AATGCCATGTCTCAAAGTTTGGTTCTATATGGAGCCTCTAAGGAGAACAGTGAAGGTTTT
CATGAAAGTAAATGACAAATACTGAAGGGGTGAATAAAGGCATTACTTTAGCTACCCA
TGTCGACGTACAGCTGTGCCGTAGTAAATATTCCAGCACCTTGTGTCAACAAAATGATT
TCACACATCCAAGATGTGGAGTCCAAAATACAGGAGCATTGAAAAGGTTTGAAACTTCT
TTTGAAGAATGGAGCAGAACTTCTTCCACAAAAGACCTGAAAGAAGATTGGAGTGAAC
ACACCAGTGAAAGAGGTCAAACCAGGAGAAAAGAGAGATGAAAAGTGTCCAGAGTTAAAG
CAGGAAATGGAACATTGCTCTCAGAGGCCATTGCTCTCATTAAAAGTCTAGAACTGAC
CGGGCAGACGCTGAAGAAGCTTTAAAACAACAGAGATCAAGAAAGAATATGATTAACATG
AAAATTGACTCTTGGTCAGTCTGGAACCTCAAGAAGTCCCATTTGGCTGTGCAGAAAGAA
CATGAGGCCTATTTGAGTGATGTTATAGAATTACAATGGCATCTTGAAGATAAAGCTAAT
CAACTACAACATTTTGAAGAACAAAAGACAGAGTTAGAAGAAGCAAATGCAAAGATTCAA
GCAGACATAGACTACATGAATGAACATGGCCCTCTACTGGACTCTAAGCAGAATCAGGAA
CTTCAAGATCTGAAGAACCATTATAAAAAAAAAATGGAGGTAATGGACCTACACAGAAAA
GTTAATGAAGAACTTGAAGAAGCTTTAGAAGCCTGTGAAAATGCCAGATTGAAGGCTCAG
CAAATTAAGAAGAGATTGATAAGGATATTTACCAGGATGAAAAACCATAGAGGCCTAC
AAGAGAGAGATATCAACTTAACAGTCTATTTGATCATTACTCTTCATCAGTGATAAT
GTTAATACTAATATTGAGGAGAAGGAAGAGGAAGTGAAGCAATAAGGGAACAAAG
TCATCAAAAAATGAATTACATTTCTATCAAAAATGCTGGAAGATTTGAGAAGAGTTTAT
GACCAACTAACCTGGAAGCAAAAAAGTCATGAAAATCAGTATCTGGAAGCAGTTAATGAT
TTTTATGCTGCAAAAAAACATGGGATATTGAGCTTTCTGATGTTGCAAAAGATTTTCA
GCTATTTCTTTGGCATGTACAAACTGACGGAAGACAATAAAAACTTGAGATTGATATT
AACAAAATAACAGTAAAAACCAATGAAAGCATACGGAATAATCAAAATACGAATCTGAA
ATAAAATATTTGACAATAATGAAGTTAAAGATGATAAACATCTCAAGAACATCTATAAG
GAGGCTTATCGCATTGGTACTCTTTCCACCTAACCAACACAAGACAGATGAAATGGAA
GATAAAATAGCAGAAGTGAGAAGAAAGTTCAAGGGTAGAGAAGAATTCCTGAAAAACTC
ACTCAAGGTGAAGTGGCTGCTGGAATGGTGCTTCAGAAAAACTATATTCCATTTACGAA
GTCCAGGCACTTGAGCGGAAAGAGCTTATAAAAAATAGAGCAATATGTGCCATGTCACTG
GCAGAAGTACAGGAACCTCTGCTTCAACTAGAAGATGAAGCTGAAAGAATCAGAAGTCTC
GACAAAGAACATTCTGTTAGTAAACGCTCAGCAATTTTTAAAGACCTAGAAGCAACTAA
AGTAAGACAATGATTTTTATGCAAAAAATAATGAATTGAATGAGGAATTAAGCAAAA
GAAGAAGAAAAGAAAAGTTTGTATCAGACACTTGAATATTAAGAACAAATTTATAACT
ATGAGATTTAAAGGGAACATGCACAACTGTGTTTGATCATTATATGCAAGAGAAAAAA
GACTGTGAAGAGAGAATCTTTGAGGAAGATCAGAGATTTAGAGTGCTCCTTGCTGTAAGA
CAAAAACTCTTCAAGATACCAAAAAATAATAGCTGATTCATTGAAGAAAACTGCGT
TTAGCTCAAGAGTATCAACAGCTACAGATTACATTCTTAAAGAAAAGGACAATTATTTT
AATATATATGATAAACAGCTATCACTTGATACTTCAATTAGAGATAAGAAACAGCTCTGT
CAGCTGCAGAGAAGGATGCACAACTGTGGCAGGAGCACTTCAAAGTGGTGGTCTCTTC
AGCCAGATGAGGCTGGCCAACTCCAGACAGACTCTCAGGAGAGTATTGAGAAAATATTA
GCTGTGCAGGAGGAATCTTCAAATTTAATGCAACACATCTTAGGTTTCTTCCAGACTTTG
ACAGATGGCACATGCGAAAACGATGGTTAA

```

**Restriction Sites:** Please inquire

**ACCN:** NM\_198995

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_198995.1, NP_945346.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq Size:</b>           | 3087 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 2490 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 374864                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>            | <u>Q5BJE1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 18q12.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |