

## Product datasheet for **SC309817**

### ARSK (NM\_198150) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ARSK (NM_198150) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARSK
Synonyms:	TSULF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC309817 representing NM\_198150.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCTACTGCTGTGGGTGTGGTGGTGCAGCCTTGGCGTGGCGGTACTGGCCCCGGAGCAGGGGAG
CAGAGGCGGAGAGCAGCCAAAGCGCCCAATGTGGTCTGGTCTGAGCGACTCCTTCGATGGAAGTTA
ACATTTTCATCCAGGAAGTCAGGTAGTGAACTTCCTTTTATCAACTTTATGAAGACACGTGGGACTTCC
TTTCTGAATGCCTACACAACTCTCCAATTTGTTGCCATCACGCGCAGCAATGTGGAGTGGCCTCTTC
ACTCACTTAACAGAATCTTGAATAATTTTAAGGGTCTAGATCCAAATTATACAACATGGATGGATGTC
ATGGAGAGGCATGGCTACCGAACACAGAAATTTGGGAACTGGACTATACTTCAGGACATCACTCCATT
AGTAATCGTGTGGAAGCGTGGACAAGAGATGTTGCTTTCTTACTCAGACAAGAAGCAGGCCCATGGTT
AATCTTATCCGTAACAGGACTAAAGTCAGAGTGATGGAAGGGATTGGCAGAATACAGACAAAGCAGTA
AACTGGTTAAGAAAGGAAGCAATTAATTACACTGAACCATTTGTTATTTACTTGGGATTAATTTACCA
CACCTTACCCTTACCATCTTCTGGAGAAAATTTGGATCTTCAACATTTACACATCTCTTTATTGG
CTTGAAAAAGTGTCTCATGATGCCATCAAAATCCCAAAGTGGTCACCTTTGTGAGAAATGCACCCTGTA
GATTATTACTCTTCTTATACAAAAAAGTGCCTGGAAAGATTTACAAAAAAGAAATTAAGAAATTTAGA
GCATTTTATTATGCTATGTGTGCTGAGACAGATGCCATGCTTGGTGAATTTTGGCCCTTCATCAA
TTAGATCTTCTTCAGAAAATTTGTCATATACTCCTCAGACCATGGAGAGCTGGCCATGGAACATCGA
CAGTTTTATAAAATGAGCATGTACGAGGCTAGTGCACATGTTCCGCTTTTGATGATGGGACCAGGAAT
AAAGCCGGCCTACAAGTATCAAATGTGGTTTCTTGTGGATATTTACCCTACCATGCTTGATATTGCT
GGAATTCCTCTGCCTCAGAACCTGAGTGGATACTTTTGTGGCATTATCATCAGAAACATTTAAGAA
GAACATAAAGTCAAAAACCTGCATCCACCCTGGATTCTGAGTGAATCCATGGATGTAATGTGAATGCC
TCCACCTACATGCTTCGAACTAACCACTGGAAATATATAGCCTATTCGGATGGTGCATCAATATTGCT
CAACTCTTTGATCTTTCCTCGGATCCAGATGAATTAACAATGTTGCTGTAAAATTTCCAGAAATTA
TATTCTTTGGATCAGAAGCTTCAATCCATTATAAACTACCCTAAAGTTTCTGCTTCTGTCCACCAGTAT
AATAAAGAGCAGTTTATCAAGTGGAAACAAAGTATAGGACAGAATTATTCAAACGTTATAGCAAATCTT
AGGTGGCACCAAGACTGGCAGAAGGAACCAAGGAAGTATGAAAATGCAATTGATCAGTGGCTTAAAAAC
CATATGAATCCAAGAGCAGTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

- Restriction Sites:** Sgfl-Mlul
- ACCN:** NM\_198150
- Insert Size:** 1611 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_198150.2</u>
<b>RefSeq Size:</b>	3345 bp
<b>RefSeq ORF:</b>	1611 bp
<b>Locus ID:</b>	153642
<b>UniProt ID:</b>	<u>Q6UWY0</u>
<b>Cytogenetics:</b>	5q15
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>MW:</b>	61.5 kDa
<b>Gene Summary:</b>	Sulfatases (EC 3.1.5.6), such as ARSK, hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans, and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling, and degradation of macromolecules (Sardiello et al., 2005 [PubMed 16174644]).[supplied by OMIM, Mar 2008]