

Product datasheet for **SC309818**

AFAP (AFAP1) (NM_198595) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AFAP (AFAP1) (NM_198595) Human Untagged Clone
Tag:	Tag Free
Symbol:	AFAP1
Synonyms:	AFAP; AFAP-110; AFAP110
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC309818 representing NM_198595.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAAGAGTTAATAGTTGAACCTCGTCTCTTTCTTGAACCTCTGGACCATGAATATCTAACCTCAACT
GTCAGGGAGAAAAAGGCAGTGATAACCAACATTCTGCTAAGAATACAGTCATCCAAGGTTTTGATGTG
AAGGACCATGCTCAGAAGCAGGAGACCGCTAACAGCCTGCCAGCCCTCCTCAGATGCCCTGCCGGAG
ATCCCTCAGCCCTGGCTGCCTCCTGACAGTGGGCTCCACCATTGCCAACATCCTCCCTCCAGAAGGT
TATTATGAGGAAGCTGTGCCGTGAGCCCCGAAAAGCTCCGGAATACATCACATCAAATATGATTCC
GATGCGATGAGCAGCTCTATGAGTCGTATGATGAAGAGGAGGAGGATGGGAAGGGGAAGAAAACCCGG
CACCAGTGGCCCTCCGAGGAGCCTCCATGGACCTGGTCAAGGACGCCAAAATCTGCGCCTTCTGCTG
CGGAAGAAGCGGTTCCGCCAGTGGACCAAGTTGCTCTGCGTCATCAAAGACACCAAATGCTGTGCTAT
AAAAGTTCCAAGGACCAGCAGCCTCAGATGGAACCTCCACTCCAAGGCTGTAACATTACGTACATCCCG
AAAGACAGCAAAAAGAAGCAGCAGCTGAAGATTACTCAGCAGGGCACGGACCCGCTTGTTCGCC
GTCCAGAGCAAGGAACAGGCCGAGCAGTGGCTGAAGGTGATCAAAGAAGCCTACAGTGGTTGTAGTGGC
CCCGTGGATTAGAGTGTCTCCTCCACCAAGCTCCCCGGTGCACAAGGCAGAATGGAGAAGAAAATG
TCTTCAGAGAGACCCAGCTCAGATGGGGAGGGTGTGTGGAAAATGGAATTACCACATGTAATGGAAAG
GAGCAAGTGAAGAGGAAGAAAAGTTCCAATCAGAGGCCAAGGGCACTGTGTCGAAAAGTCACTGGGAAA
AAAATCACAAGATCATCAGTCTGGGAAAGAAAAGCCGTCCACAGACGAGCAGACCTCCTCAGCTGAG
GAAGATGTTCCACCTGCCGCTATCTGAACGTGCTCTCCAACAGCCGCTGGCGAGAGCGCTGGTGCCGA
GTGAAAGATAACAAGCTCATTTTCCACAAGGACAGGACCGACCTGAAGACCCATATTGTGTCTATTCCG
TCCCGTGGCTGCGAGGTGATCCCGGGTTTGGATTCTAAACATCCTCTGACGTTCCGGCTGCTGCGCAAC
GGCCAGGAGTTGCAGTATTGGAGGCATCTTCTTGAAGACATGGGCAGGTGGATTGGGATTTTACTC
GCAGAGACGGGATCGTCCACAGACCCGGAGGCTCTGCACTATGACTACATTGATGTGGAGATGTCTGCA
AGTGTCTATTAGACAGCCAAACAGACCTTCTGTTTCATGAACAGGCGTGTATATCTGCTAACCCATAT
CTAGGGGGCACCTCCAACGGCTATGCCACCCAGCGGGACGGCACTTATTATGACGATGTCCCGTGC
ATCAACGGCTCGCTCAAGGGTAAAAGCCCCCGTGGCGTCTAATGGGGTACAGGAAAAGGGAAGACT
CTGAGCAGTCAGCCAAAGAAAGCGGATCCCGGGCTGTTGTGAAAAGGACGGGTTCAATGCTGCCAG
TACAAGTATGGCAAGAACCGGTAGAAGCAGATGCCAAGCGGCTACAGACCAAAGAGGAGGAGTGTG
AAGAGGAAAGAGGCCCTCGGAAATAGGCTGGCCAGCTCCGCAAGGAAAGAAAAGACCTTCGAGCGCT
ATTGAAGTGAACGCCGGCAGGAAGCCGAGGCGATCCTGGAGGAGAAGCTGAAGCAGCTGGAGGAGGAG
TGCCGGCAGAAGGAGGCGGAGCGTGTGACCTGGAGCTGGAGCTGACGGAGGTAAGGAGAGCCTGAAG
AAAGCGCTGGCGGGCGGAGTCAACCTGGGGCTGGCCATCGAGCCCAAGTCAGGGACATCGAGTCCACAG
TCTCCAGTGTTCGGCACCGGACCTGGAAAATCGCCATCTCCAGCTGTGACACCAGTACACCGAG
GGCCCCGTGCCGGTGAACAGCGCGGCGCTTGAAGAAGAGCCAGGCTGCCCGGGCAGCTCCCCCTGC
CGAGGGCATGTGCTGCGGAAGGCCAAGGAATGGGAATTGAAGAACGGGACCTAG
ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_198595

Insert Size: 2193 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198595.2
RefSeq Size:	7534 bp
RefSeq ORF:	2193 bp
Locus ID:	60312
UniProt ID:	Q8N556
Cytogenetics:	4p16.1
MW:	80.7 kDa
Gene Summary:	<p>The protein encoded by this gene is a Src binding partner. It may represent a potential modulator of actin filament integrity in response to cellular signals, and may function as an adaptor protein by linking Src family members and/or other signaling proteins to actin filaments. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2008]</p> <p>Transcript Variant: This variant (2) contains an alternate in-frame exon in the 3' coding region, compared to variant A. Variant 2 encodes a shorter protein (isoform B), compared to isoform A.</p>