

## Product datasheet for **SC310088**

### **FYB1 (NM\_001465) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	FYB1 (NM_001465) Human Untagged Clone
Tag:	Tag Free
Symbol:	FYB1
Synonyms:	ADAP; FYB; PRO0823; SLAP-130; SLAP130; THC3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC310088 representing NM\_001465.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGAAATATAACACGGGGGCAACCCGACAGAGGATGTCTCAGTCAATAGCCGACCCTTCAGAGTC
ACAGGGCCAAACTCATCTTCAGGAATACAAGCAAGAAAGAACTTATTCAACAACCAAGGAAATGCCAGC
CCTCTGCAGGACCCAGCAATGTACCTAAGTTTGGGTCCCAAAGCCACCTGTGGCAGTCAAACCTTCT
TCTGAGGAAAAGCCTGACAAGGAACCAAGCCCCGTTTCTAAAGCCCACTGGAGCAGGCCAAAGATTC
GGAACACCAGCCAGCTTGACCACCAGAGACCCCGAGGCGAAAGTGGGATTTCTGAAACCTGTAGGCCCC
AAGCCCATCAACTTGCCCAAAGAAGATTCCAAACCTACATTTCCCTGGCCTCCTGGAAACAAGCCATCT
CTTCACAGTGTAAACCAAGACCATGACTTAAAGCCACTAGGCCCGAAATCTGGCCTACTCCTCAAACC
TCAGAAAATGAACAGAAGCAAGCGTTTCCCAAATTGACTGGGGTTAAAGGGAAATTTATGTCAGCATCA
CAAGATCTTGAACCAAGCCCTCTTCCCAAACCCGCCCTTGGCCAGAAGCCGCCCTAAGTACCGAG
AACTCCCATGAAGACGAAAGCCCATGAAGAATGTGTCTCATCAAAGGGTCCCAGCTCCCCTGGGA
GTCAGGTCCAAAAGCGGCCCTTTAAAACAGCAAGGGAAGACTCAGAAAATAAAGACCATGCAGGGGAG
ATTTCAAGTTTGCCTTTCTGGAGTGGTTTTGAAACCTGCTGCGAGCAGGGGAGGCCAGGTCTCTCC
AAAAATGGTGAAGAAAAAAGGAAGATAGGAAGATAGATGCTGCTAAGAACACCTTCAGAGCAAAAATA
AATCAGGAAGAGTTGGCCTCAGGGACTCCTCCTGCCAGTTCCCTAAGGCCCTTCTAAGCTGACAGTG
GGGGGGCCATGGGGCCAAAGTCAGGAAAAGGAAAGGGAGACAAGAATTCAGCCACCCCGAAACAGAAAG
CCATTGCCTCCCTGTTTACCTTGGGTCCACCTCCACAAAACCAACAGACCACCAAATGTTGACCTG
ACGAAATCCACAAAACCTCTTCTGAAACAGTACTAGCAAAGGCCAGACGTCTTACTCAACAACCTCC
CTGCCACCACTCCACATCCCATCCGCGCCAGCCAAACCATTGCCAGCATCTACCCATCACAACCA
CCAGTCCCAAGCCTACCTCCAGAAACATTAACCTCCGTTTGACCTAAAAGCCCTGCAATGAAGAC
AATCAAGATGGTGTACGCACTCTGATGGTGTGGAATCTAGATGAGGAACAAGACAGTGAAGGAGAA
ACATATGAAGACATAGAAGCATCCAAAGAAAGAGAGAAGAAAAGGGAAAAGGAAGAAAAGAGAGTTA
GAGCTGGAGAAAAAGGAACAGAAAGAGAAAAGAAAAGAAACAAGAAATAAAGAAGAAATTTAAACTA
ACAGGCCCTATTCAAGTCATCCATCTTGCAAAAGCTTGTGTGATGTCAAAGGAGGAAAGAATGAACTG
AGCTTCAAGCAAGGAGAGCAAAATGAAATCATCCGCATCACAGACAACCCAGAAGGAAAATGGTTGGGC
AGAACAGCAAGGGGTTTATATGGCTATATTAACAACCTGCTGTAGAGATTGACTATGATCTTTGAAA
CTGAAAAAAGACTCTCTTGGTGGCCCTTCAAGACCTATTGAAGATGACCAAGAAGTATATGATGATGT
GCAGAGCAGGATGATATTAGCAGCCACAGTCAGAGTGGAAAGTGGAGGGATATTCCCTCCACCACCAGAT
GATGACATTTATGATGGGATTGAAGAGGAAGATGCTGATGATGGCTCCACACTACAGTTCAAGAGAAG
AGTAATACGTGGTCTGGGGGATTTTGAAGATGTTAAAGGGAAAAGATGACAGAAAAGAAAAGTATACGA
GAGAAACCTAAAGTCTCTGACTCAGACAATAATGAAGGTTTCATCTTTCCCTGCTCCTCTAAACAATTG
GACATGGGAGATGAAGTTTACGATGATGTGGATACCTCTGATTTCCCTGTTTCATCAGCAGAGATGAGT
CAAGGAACTAATGTTGAAAAAGCTAAGACAGAAGAAAAGGACCTTAAAGAGCTAAAAAGCAGGAAAAA
GAAGAAAAAGACTTCAGGAAAAAATTTAAATATGATGGTGAATTAGAGTCTTATTTCAACTAAAGTT
ACAACCTCCATAACTTCTAAAAAGTGGGAACCCAGAGATCTACAGGTAACCTGGTGAATCTCTAGAA
GTTATACAAACCACAGATGACACAAAAGTTCTCTGCAGAAATGAAGAAGGAAATATGGTTATGTCCTT
CGGAGTTACCTAGCGGACAATGATGGAGAGATCTATGATGATATTGCTGATGGCTGCATCTATGACAAT
GACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001465  
**Insert Size:** 2490 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001465.4</a></u>
<b>RefSeq Size:</b>	4858 bp
<b>RefSeq ORF:</b>	2490 bp
<b>Locus ID:</b>	2533
<b>UniProt ID:</b>	<u><a href="#">O15117</a></u>
<b>Cytogenetics:</b>	5p13.1
<b>Domains:</b>	SH3
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	90.6 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene is an adapter for the FYN protein and LCP2 signaling cascades in T-cells. The encoded protein is involved in platelet activation and controls the expression of interleukin-2. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011]</p> <p>Transcript Variant: This variant (1) differs in the 5' UTR and uses a downstream start codon compared to variant 3. The resulting isoform (1) has a shorter N-terminus compared to isoform 3. Variants 1 and 4 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>