

Product datasheet for **SC310097**

Neuroigin 4 (NLGN4Y) (NM_014893) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Neuroigin 4 (NLGN4Y) (NM_014893) Human Untagged Clone
Tag:	Tag Free
Symbol:	Neuroigin 4
Synonyms:	HNL4Y
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310097 representing NM_014893.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTTGCGTCCCCAGGGACTGCTATGGCTCCCTTTGTTGTTACCTCTGTCTGTGTGATGTTAAACTCC
AATGTTCTTCTGTGGATAACTGCTCTTGCCATCAAGTTCACCCTCATTGACAGCCAAGCACAGTATCCA
GTTGTCAACACAATTATGGTAAAAATCCAGGGCCTAAGAACCATTACCCAGTGAGATCTTGGGTCCA
GTGGAGCAGTACTTAGGGTCCCTATGCCTACCCCCAACTGGAGAGAGCGGTTTCAGCCACCAGAA
TCCCCATCCTCTGGACTGGCATCCGAAATGCTACTCAGTTTTCTGCTGTGTGCCCCAGCACCTGGAT
GAAAGATTCTTATTCATGACATGCTGCCATCTGGTTTACCACCAGTTGGATACTTTGATGACCTAT
GTTCAAGATCAAATGAAGACTGCCTTTACTTAAACATCTATGTGCCATGGAAGATGATATTCATGAA
CAGAACAGTAAGAAGCCTGTTATGGTCTATATCCATGGGGATCTTACATGGAGGAAACCGTAACATG
ATTGATGGCAGCATTTTGGCCAGCTATGGGAACGTCATCGTTATCACCATTAACACCGTCTGGGAATA
CTAGGGTTTTTAAGTACCGGTGACCAGGCAGAAAAGGCAACTATGGGCTCCTGGATCAGATTCAAGCA
CTGAGGTGGATTGAGGAGATGTCGGAGCCTTTGGCGGGACCCCAAGAGAGTACTATCTTTGGCTCG
GGGCTGGGGCCTCCTGTGTGAGCCTGTTGACCCTGTCCCACTACTCAGAAGGTCTCTTCCAGAAGGCC
ATCATTGAGAGCGGCACTGCCCTGTCCAGCTGGGCACTGAACTACCAGCCGGCAAGTACTCGGATA
TTGGCAGACAAGTGGCTGCAACATGCTGGACACCACGGACATGGTAGAATGTCTGAAGAACAAGAAC
TACAAGGAGTCTATCCAGCAGACCATCACCCGGCCACCTACCACATAGCCTTTGGGCGGTGATCGAC
GGCAGCTCATCCAGACGACCCAGATCCTGATGGAGCAAGGCGAGTTCCTCAACTACGACATCATG
CTGGGCGTCAACCAAGGGGAAGGCCTGAAGTTCGTGGACGGCATCGTGGATAACGAGGACGGTGTGACG
CCCAACGACTTTGACTTCTCCGTGTCCTCAACTTCGTGGACAACCTTTACGGCTACCTGAAGGGAAGAC
ACTTTGCCGGGAGACTATCAAGTTCATGTACACAGACTGGGCCGATAAGGAAAACCCGGAGACGCGGCGG
AAAACCTGGTGGCTCTCTTTACTGACCATCAGTGGGTGGCCCCCGCGTGGCCACCGCCGACTGCAC
GCGCAGTACGGCTCCCCACCTACTTCTATGCCTTCTATCATCACTGCCAAAGCGAAATGAAGCCAGC
TGGGCAGATTCGGCCATGGCGATGAAGTCCCCTATGTCTTCGGCATCCCCATGATCGGTCCCACAGAG
CTCTTCAGTTGTAATTTCTCCAAGAACGACGTCATGCTCAGTCCCGTGGTATGACCTACTGGACGAAC
TTGCCAAAACCTGGTATCCAAACCAACGATTCCTCAGGATACCAAGTTCATTCATACAAAACCCAAT
CGCTTTGAAGAAGTGGCCTGGTCCAAGTATAATCCAAAGACCAGCTCTATCTGCATATTGGCTTGAAA
CCCAGAGTGAGAGATCACTACCGGCAACGAAAGTGGCTTCTGGTTGGAATTGGTTCCTCATTTGCAC
AACTTGAACGAGATATCCAGTATGTTTCAACAACCACAAGGTTCTCCACCAGACATGACATCATTT
CCCTATGGCACCCGGCGATCTCCCGCAAGATATGGCCAACCACCAACGCCAGCAATCACTCCTGCC
AACAAATCCAAACACTCTAAGGACCCTCACAAAACAGGGCCCGAGGACACAACCTGTCTCATTGAAACC
AAACGAGATTATCCACCGAATTAAGTGTACCAATTGCCGTGGGGCGTCCGCTCTTCTCAACATC
TTAGCCTTTGGCGCCTGTACTACAAAAAGGACAAGAGGCGCCATGAGACTCACAGGCACCCAGTCCC
CAGAGAAACACCACAAATGATATCACTCACATCCAGAACGAAGAGATCATGTCTCTGCAGATGAAGCAG
CTGGAACACGATCAGAGTGTGAGTCGCTGCAGGCACACGACACGCTGAGGCTCACCTGCCCTCCAGAC
TACACCCTCACGCTGCGCCGGTCCCGGATGACATCCCATTTATGACGCCAAACACCATCACCATGATT
CCAAACACATTGATGGGGATGCAGCCTTTACACACTTTTAAAACCTTTCAGTGGAGGACAAAACAGTACA
AATTTACCCACGGACATTCCACCACTAGAGTATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_014893
Insert Size: 2451 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_014893.4
RefSeq Size:	5409 bp
RefSeq ORF:	2451 bp
Locus ID:	22829
UniProt ID:	Q8NFZ3
Cytogenetics:	Yq11.221
Domains:	COesterase
Protein Families:	Druggable Genome, Transmembrane
MW:	92 kDa
Gene Summary:	<p>This gene encodes a type I membrane protein that belongs to the family of neuroligins, which are cell adhesion molecules present at the postsynaptic side of the synapse, and may be essential for the formation of functional synapses. Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Mar 2011]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (1).</p>