

Product datasheet for **SC310102**

RASA4 (NM_006989) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RASA4 (NM_006989) Human Untagged Clone
Tag:	Tag Free
Symbol:	RASA4
Synonyms:	CAPRI; GAPL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310102 representing NM_006989.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCCAAGCGCAGCTCGCTGTACATCCGCATCGTGGAGGGGAAGAACCTTCCCGCCAAGGACATCACT
GGCAGCAGCGACCCCTACTGCATCGTGAAGGTGGACAATGAGCCCATCATCAGGACAGCCACAGTGTGG
AAGACCCCTGTGCCCTTCTGGGGTGAGGAGTACCAAGTGACCTGCCGCCACCTTCCACGCTGTGGCT
TTCTATGTCATGGATGAGGATGCCCTCAGCCGGGACGACGTTATCGGAAAGGTCTGCCTTACAAGGGAC
ACCATAGCCTCTCACCTAAGGGTTTCAGCGGGTGGGCCACCTGACGGAGGTGACCCCGATGAGGAG
GTGCAGGGCGAGATCCACCTGCGGCTGGAAGTGTGGCCAGGGGCCGGGCTGCCGGCTACGCTGCTCT
GTGCTGGAGGCCAGGGATCTGGCCCAAGGACCGCAATGGCACATCTGACCCCTTCGTCCGAGTGCGC
TACAAGGGCCGGACACGGGAGACCTCGATCGTGAAGAAGTCATGCTACCCACGCTGGAATGAGACGTTT
GAATTTGAGCTGCAGGAGGGGCCATGGAGGCGCTGTGCGTGGAGGCTGGGACTGGGACCTTGTACGC
CGAAACGACTTCTGGGCAAAGTGGTGATTGATGTCCAGAGACTGCGGGTGGTGCAGCAGGAGGAGGGC
TGGTTCCGGCTGCAGCCGACCAAGTCCAAGAGCCGGCGCATGACGAGGGCAACCTGGGTCCTTGACG
CTGGAGGTGCGGCTGCGGGACGAGACGGTGCTGCCCTCCAGCTACTACGACCACTGGTGCACCTGCTG
TGCCACGAGGTCAAGCTGGGATGCAGGGCCAGGGCAGCTGATCCCACTCATCGAGGAGACAACGAGC
ACCGAGTGTGCCAGGACGTGGCCACGAACCTGCTCAAGCTCTTCTGGGGCAGGGGCTGGCCAAGGAC
TTCCTGGACCTGCTCTTCCAGCTGGAGCTGAGTCGCACCAAGTGCAGCAACACCTGTTCCGGAGCAAC
TCTCTGGCCTCAAAGTCCATGGAGTCTTTCTGAAGGTGGCCGGGATGCAGTACCTGCACGGCGTCTG
GGCCCCATCATCAACAAGGTGTTTGAGGAGAAGAAGTACGTGGAGCTGGACCCAGCAAAAGTGAAGTT
AAGGATGTAGGGTGCTCCGGGCTGCACCGCCCGCAGACCGAGGCCGAGGTGCTGGAGCAGAGCGCGCAG
ACGCTGCGCGCCACCTGGGGGCCCTGCTGAGCGCGCTCAGCCGCTCGGTTGCGCGGTGCCCGCGGTG
GTGCGCGCCACCTTCCGCCAGCTCTTCCGGCGCTGCGCGAGCGCTTCCCGCGCGCCAGCACGAGAAT
GTACCGTTTCATCGCCGTACCAAGCTTCTGTGCTGCGCTTCTTCTCTCCCGCCATCATGTGCCCAAG
CTCTTCCACCTGCGGGAGCGCCACGCGGACGCCCGACCAAGCCGACCCCTGCTCTGTTGCCAAGGCA
GTCCAGAACGTGGGCAACATGGACACGCCGGCTTCCAGGGCCAAGGAGGCTTGGATGGAGCCGCTGCAG
CCCACCGTGCGCCAGGGCGTGGCGCAGCTGAAGGACTTCATCACAAGCTCGTGGACATCGAGGAGAAG
GACGAGCTGGACCTGCAGCGACGCTGAGTTTGAGGCGCCACCTGTGAAGGAGGGGCCACTCTTCATC
CACAGGACCAAGGCAAGGGCCCCCTCATGTCTCTCTCTCAAGAAGCTCTACTTCTCCCTCACTACC
GAGGCCCTCAGCTTCGCGAAGACGCCAGCTCCAAGAAAAGCGCCCTCATCAAGTTAGCCAACATCCGG
GCAGCGGAAAAGGTTGAGGAAAAGAGCTTTGGCGGCTCGCACGTCATGCAGGTATCTACACGGACGAC
GCCGGCAGGCCCCAGACTGCCTACCTGCAGTGCAAGTGTGTGAATGAGCTTAACCAAGTGGCTGTCTGCG
CTGCGGAAGGTGAGCATCAACAACACGGGACTGCTGGGCTCCTACCACCTGGCGTCTTCCGTGGGGAC
AAGTGGAGCTGCTGCCACCAAAAGAGAAGACAGGTGAGGGCTGCGATAAGACCCGGTCACGGGTGACC
CTGCAGGAGTGAATGACCTCTTGACCATGACCTTGAAGGCCAGCTCATCTACCGGCACCTGCTGGGC
GTGGAGGCCATGCTGTGGAGAGGCACCGGAGCTGAGCGGGGGCGCAGAGGCAGGCACGGTGCCACG
AGCCCTGGCAAAGTCCCCGAGGACTCATTGGCCCGGCTGCTCCGGGTGCTGCAGGACCTCCGCGAGGCC
CATAGCTCCAGCCCGGGCGGCTCCCAACCTCAGAGCCCAACTGCCTCCTGGAGCTGCAGACGTGA
ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_006989
Insert Size: 2412 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_006989.5</u>
RefSeq Size:	5606 bp
RefSeq ORF:	2412 bp
Locus ID:	10156
UniProt ID:	<u>O43374</u>
Cytogenetics:	7q22.1
Domains:	C2, BTK, PH, RasGAP
Protein Families:	Druggable Genome
MW:	90.5 kDa
Gene Summary:	This gene encodes a member of the GAP1 family of GTPase-activating proteins that suppresses the Ras/mitogen-activated protein kinase pathway in response to Ca(2+). Stimuli that increase intracellular Ca(2+) levels result in the translocation of this protein to the plasma membrane, where it activates Ras GTPase activity. Consequently, Ras is converted from the active GTP-bound state to the inactive GDP-bound state and no longer activates downstream pathways that regulate gene expression, cell growth, and differentiation. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.