

Product datasheet for **SC310110**

PHF15 (JADE2) (NM_015288) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PHF15 (JADE2) (NM_015288) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHF15
Synonyms:	JADE-2; PHF15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310110 representing NM_015288.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAAGAGAAGAGGCGAAAACTCCATCAGCAGTGACAACCTGACACCACTGACAGTCATGCGACA
TCTACATCCGCATCAAGATGCTCCAACTGCCAGCAGCACCAAGTCGGGCTGGCCCCGACAGAACGAA
AAGAAGCCCTCCGAGGTTTTCCGGACAGACTTGATCACAGCCATGAAGATCCCGGACTCATACCAGCTC
AGCCCGGATGACTACTACATCCTGGCAGACCCATGGCGACAGGAATGGGAGAAAGGTGTGCAGGTGCCT
GCCGGGGCAGAGGCCATCCCAGAGCCCGTGGTGAGGATCCTCCACCACTGGAAGCCCCCTGCCAG
GCATCCCCGAGCAGCACCATGCTTGGTGAGGGCTCCCAGCCTGATTGGCCAGGGGCAGCCGCTATGAC
TTGGACGAGATTGATGCCTACTGGCTGGAGCTCATCAACTCGAGCTTAAGGAGATGGAGAGGCCGGAG
CTGGACGAGCTGACATTAGAGCGTGTGCTGGAGGAGCTGGAGACCCTGTGCCACCAGAATATGGCCAGG
GCCATTGAGACGAGGAGGGGCTGGGCATCGAGTACGACGAGGATGTTGTCTGCGACGTGTGTCGCTCT
CCTGAGGGCGAGGATGGCAACGAGATGGTCTTCTGTGACAAGTGAACGCTGTGTGCATCAGGCATGC
TACGGGATCCTCAAGGTGCCACGGGCAGCTGGCTGTGCCGAGCTGTGCCCTGGGTGTCCAGCCAAAG
TGCTGTCTGTGCCCAAGCGAGGAGGAGCCTTGAAGCCCACTAGAAGTGGGACCAAGTGGGTGCATGTC
AGCTGTGCCCTATGGATTCTGAGGTGAGCATCGGCTGCCAGAGAAGATGGAGCCCATCACCAAGATC
TCGCATATCCCAGCCAGCCGCTGGGCTCTGTCTGCAGCCTCTGCAAGGAATGCACAGGCACCTGCATC
CAGTGTTCCATGCCTTCTGCGTCACAGCGTTCATGTGCATGCGCCTTTGACCACGGCCTGGAATG
CGGACTATATTAGCAGACAACGATGAGGTCAAGTTCAAGTCATTCTGCCAGGAGCAGTGACGGGGGC
CCACGTAATGAGCCACATCTGAGCCCACGGAACCCAGCCAGGCTGGCGAGGACCTGGAAGAGGTGACC
CTGCGCAAGCAGCGGCTGCAGCAGCTAGAGGAGGACTTCTACGAGCTGGTGAGCCGGCTGAGGTGGCT
GAGCGGCTGGACCTGGCTGAGGCACTGGTCGACTTCATCTACCAGTACTGGAAGCTGAAGAGGAAAGCC
AATGCCAACCAGCCGCTGCTGACCCCCAAGACCGACGAGGTGGACAACCTGGCCAGCAGGAGCAGGAC
GTCTCTACCGCCGCTGAAGCTCTTACCCATCTGCGGCAGGACCTAGAGAGGGTTAGAAATCTGTGC
TACATGGTGACAAGGCGCGAGAGAACGAAACACGCCATCTGCAAACTCCAGGAGCAGATATTCCACCTG
CAGATGAACTTATTGAACAGGATCTGTGTCGAGGCTGTCCACCTATTCCCCATCGATGGCACCTTC
TTCAACAGCTGGCTGGCAGCTCGGTGCAGATCACAGCAGAGAACATGGCCATGAGCGAGTGGCCACTG
AACATGGGCACCGCGAGGACCCTGCTCCAGGGCTGCTGTGAGAGGAAGTGTGACGAGCAGGAGACA
CTGCTCAGCTTATGCGGGACCCCTCGTGCACCTGGTGACCTGCTAGGAAGGCCCGAGGCCGCACC
CGCTGCCTGCCAAGAAGAAACCAACCACCAACCACCGCAGGACGGGCTGGTTACGGACGACTCCA
GACAAAGCCCCCAAGAAGACCTGGGGCCAGGATGCAGGCAAGTGGCAAGGGGGGTCAAGGGCCACCTACC
AGGAAGCCACACGTCGGACATCTTCTCACTTGCCGTCCAGCCCTGCAGCCGGGGACTGTCCCATCTTA
GCCACCCCTGAAAGCCCCCGCCACTGGCCCTGAGACCCCGACGAGGACGCTCAGTAGCTGTGAC
TCAGATGTCCAAGTGCCTGGCCCTGCAGCAAGCCCTAAGCCTTTGGGCCGGCTCCGGCCACCCCGCGAG
AGCAAGGTAACCCGAGATTGCCGGGTGCCAGGCTGATGCTGGGATGGGACACCTTCAGCTGTGGCT
GAGAGGGCCCAAGGTCAGCCTGCATTTTGACACTGAGACTGATGGTACTTCTCTGATGGGAGATGAGC
GACTCAGATGTAGAGGCCGAGGACGGTGGGGTGCAGCGGGTCCCGGGAGGCAGGGGCAGAGGAGGTG
GTCCGCATGGGCGTACTGGCCTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_015288
Insert Size: 2373 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015288.5</u>
RefSeq Size:	6478 bp
RefSeq ORF:	2373 bp
Locus ID:	23338
UniProt ID:	<u>Q9NQC1</u>
Cytogenetics:	5q31.1
Protein Families:	Druggable Genome, Transcription Factors
MW:	87.5 kDa
Gene Summary:	Component of the HBO1 complex which has a histone H4-specific acetyltransferase activity, a reduced activity toward histone H3 and is responsible for the bulk of histone H4 acetylation in vivo (PubMed:16387653). Acts as a E3 ubiquitin-protein ligase mediating the ubiquitination and subsequent proteasomal degradation of target protein histone demethylase KDM1A (PubMed:25018020). Also acts as a ubiquitin ligase E3 toward itself. Positive regulator of neurogenesis (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) uses an alternate in-frame splice site in the 3' coding region, compared to variant 1. The encoded isoform (3) is shorter than isoform 1.