

Product datasheet for **SC310111**

K Cadherin (CDH6) (NM_004932) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	K Cadherin (CDH6) (NM_004932) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDH6
Synonyms:	CAD6; KCAD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310111 representing NM_004932.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAGAATTACCGTACTTCTTGCTCTTTTGGGTGGGCAGCCCTACCAACTCTCTCAACTCCA
CTATCAAAGAGGACTAGTGGTTTCCCAGCAAAGAAAAGGGCCCTGGAGCTCTCTGGAACAGCAAAAAT
GAGCTGAACCGTTCAAAAAGGAGCTGGATGTGGAATCAGTTCTTTCTCCTGGAGGAATACACAGGATCC
GATTATCAGTATGTGGGCAAGTTACATTACAGACCAGGATAGAGGAGATGGATCACTTAAATATATCCTT
TCAGGAGATGGAGCAGGAGATCTCTTATTATTAATGAAAACACAGGCGACATACAGGCCACCAAGAGG
CTGGACAGGGAAGAAAACCCGTTTACATCCTTCGAGCTCAAGCTATAAACAGAAGGACAGGAGACCC
GTGGAGCCCGAGTCTGAATTCATCATCAAGATCCATGACATCAATGACAATGAACCAATATTCACCAAG
GAGGTTTACACAGCCACTGTCCCTGAAATGTCTGATGTCGGTACATTTGTGTCCAAGTCACTGCGAGC
GATGCAGATGATCCAACATATGGGAACAGTGCTAAAGTTGTCTACAGTATTCTACAGGGACAGCCCTAT
TTTTCAGTTGAATCAGAAAACAGTATTATCAAGACAGCTTTGCTCAACATGGATCGAAGAAAACAGGGAG
CAGTACCAAGTGGTGATTCAAGCCAAGGATATGGGCGGCCAGATGGGAGGATTATCTGGGACCACCACC
GTGAACATCACACTGACTGATGTCAACGACAACCCCTCCCGATTCCCCAGAGTACATACCAGTTTAAA
ACTCCTGAATCTTCTCCACCGGGGACACCAATTGGCAGAATCAAAGCCAGCGACGCTGATGTGGGAGAA
AATGCTGAAATTGAGTACAGCATCACAGACGGTGGGGCTGGATATGTTTGATGTCATACCGACCAG
GAAACCCAGGAAGGGATTATAACTGTCAAAAAGCTCTGGACTTTGAAAAGAAGAAAGTGTATACCCCTT
AAAGTGAAGCCCTCAATCCTTATGTTGAGCCACGATTTCTACTTGGGCCTTTCAAAGATTCAGCC
ACGTTTGAATTTGGTGGAGGATGTAGATGAGCCACCTGTCTTCAGCAAACCTGGCCTACATCTTACAA
ATAAGAGAAGATGCTCAGATAAACACACAATAGGCTCCGTACAGCCCAAGATCCAGATGCTGCCAGG
AATCCTGTCAAGTACTCTGTAGATCGACACACAGATATGGACAGAATATCAACATTGATTCTGGAAAT
GGTTCGATTTTTACATCGAAACTTCTTGACCGAGAAAACACTGCTATGGCACAACATTACAGTATAGCA
ACAGAGATCAATAATCCAAGCAAAGTAGTCGAGTACCTCTATATATAAAGTTCTAGATGTCAATGAC
AACGCCCCAGAATTTGCTGAGTTCTATGAAACTTTTGTCTGTGAAAAAGCAAAGGCAGATCAGTTGATT
CAGACCTGCATGCTGTTGACAAGGATGACCCTTATAGTGGACACCAATTTTCGTTTTCTTGGCCCT
GAAGCAGCCAGTGGCTCAAACCTTACCATTCAAGACAACAAGACAACACGGCGGGAATCTTAACCTCGG
AAAAATGGCTATAATAGACACGAGATGAGCACCTATCTCTTGCCTGTGGTCATTTACAGACAACGACTAC
CCAGTTCAAAGCAGCACTGGGACAGTACTGTCCGGTCTGTGCATGTGACCACCAGGGAACATGCAA
TCCTGCCATGCGGAGGCGCTCATCCACCCACGGGACTGAGCACGGGGCTCTGGTTGCCATCCTTCTG
TGCATCGTGATCCTACTAGTGACAGTGGTGTGTTTGCAGCTCTGAGGCGGCAGCGAAAAAAGAGCCT
TTGATCATTTTCAAAGAGGACATCAGAGATAACATTGTCAGTTACAACGACGAAGGTGGTGGAGAGGAG
GACACCCAGGCTTTTGATATCGGCACCCTGAGGAATCCTGAAGCCATAGAGGACAACAAATTACGAAGG
GACATTGTGCCGAAGCCCTTTTCTACCCGACGGACTCCAACAGCTCGCGACAACACCGATGTCAGA
GATTTTCATTAACCAAAGGTTAAAGGAAAATGACACGGACCCCACTGCCCGCCATACGACTCCTTGGCC
ACTTACGCCTATGAAGGCACTGGCTCCGTGGCGGATTCCCTGAGCTCGCTGGAGTCACTGACCACGGAT
GCAGATCAAGACTATGATTACCTTAGTGACTGGGGACCTCGATTCAAAAAGCTTGCAGATATGTATGGA
GGAGTGGACAGTGACAAAGACTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_004932

Insert Size: 2373 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_004932.3</u>
RefSeq Size:	8571 bp
RefSeq ORF:	2373 bp
Locus ID:	1004
UniProt ID:	<u>P55285</u>
Cytogenetics:	5p13.3
Domains:	Cadherin_C_term, CA
Protein Families:	Transmembrane
MW:	88.3 kDa
Gene Summary:	<p>This gene encodes a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. The encoded protein is a type II cadherin and may play a role in kidney development as well as endometrium and placenta formation. Decreased expression of this gene may be associated with tumor growth and metastasis. [provided by RefSeq, May 2011]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>