

## Product datasheet for **SC310139**

### ZNF280C (NM\_017666) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF280C (NM_017666) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF280C
Synonyms:	SUHW3; ZNF633; ZPET
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC310139 representing NM\_017666.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGATGACGACAAAACCTTTTCAACCAAAAAACATTTCAAAAATGGCAGAACTCTTTATGGAATGTGAA
GAAGAAGAGCTAGAGCCATGGCAGAAAGATAGAAGAACTCAGGATGAGGATGACGATGAACTGATC
TTTGTTGGAGAGATATCAAGTTCAAACCAGCCATTTCAAATATTTTGAACAGAGGCCACTCCAGCTCA
TCTTCAAAGGAATAAAGAGTGAGCCACACAGTCCAGGTATTCCTGAAATATTCAGGACTGCAAGTCAA
CGCTGCAGAGACCACCATCAAATCCAGTGGTGCCTCGCTAGATTTTCATCTTGTATCTAAATCTTCA
CAAAGCTCTGTTACTGTTGAGAAATGCGTCTAAACCTGATTTTACAAGAATTACAAGTTGGATCGGAT
AATTCTTCAATTTACTGTTGACTCGACCCAGGAATCACTACCACCATCCAAGACATACCAGCAATT
TTTAGAGAAGGTATGAAAAACTTTCATATGTGTTGAAACATCCTTCTACTTCTAAAGTAAACAGTGTT
ACTCAAAAAAACCAAAGACCAGTGAAGATGTTCTCAGATAAAATCCCTCCACTTCATTGCCTTTAATT
GGCTCTCTCCAGTGACATCCTCCAAGTTATGCTGTCAAAGGTACAAATACCTCATCTCCATATGAT
GCTGGAGCAGATTACCTAAGAGCTTGCCAAAGTGCAATGTTCAAGTTCAATCTTTGGATCCTTTGAAA
TACCACATGAAGCATTGTTGTCAGACATGATAACTAAATTTTGGGAGTAATTGTTAAATCAGAACGT
CCATGTGATGAAGACAAGACTGATTCAGAGACAGGAAAGTTGATCATGTTAGTCAATGAGTTTTATTAT
GGAAGGCATGAAGGAGTCACTGAGAAAGAGCCAAAGACTTACACAACCTTTAAATGCTTCAGTTGCTCG
AAAGTTCTTAAAAATAATATTAGGTTTATGAACCACATGAAACATCACTTGGAACTTGAGAAGCAGAAC
AATGAAAGCTGGGAAAACCACACCCTGCCAGCACTGTTACCGCAATATCCCACACCTTCCAAGTGC
CAGTGCCACATTGAGAGTACACACACTCCCATGAGTTTTCTACTATTTGCAAAAATCTGTGAATATCA
TTTGAACAGAGCATATCTTTTACAACATATGAAGGACACCCATAAAACCTGGTAAATGCCATATGTT
TGCCAGGTTTGCCAGTTTAGATCATCAACATTTTCTGATGTAGAAGCTCATTTTAGAGCAGCCCATGAA
AACACTAAGAACTTGCTATGTCCATTTTGCCTCAAAGTTAGTAAAATGGCAACCCCTACATGAATCAT
TACATGAAGCATCAGAAAAAGGAGTTCATCGTTGCCAAAATGCAGACTACAATTTTGGACAGCAAG
GAGAAAGCTGAACATAAGGCGCAGCATCGTACATTTATAAAGCTAAAGAAGTAGAAGGATTGCCTCCT
GGAGCAAAAGTTACTATTCGAGCTTCACTTGGACCTCTCCAGTCAAAATACCAACTGCACCTTTCGGT
TGCGCTCCAGGCACTTCTTTTCTCAGGTCACACCTCCGACATCTCAAAATACAAGTCTAGAAATCCT
AGAAAATCTAATGCCAGTAGATCTAAGACAAGTAAGCTTCATGCAACTACATCCACTGCAAGTAAAGTT
AATACAAGTAAGCAAGGGGACGTATAGCTAAGTCAAAGCAAAACCTCTTACAAGCAAAAGCGACAG
CGCAACAGAAAAATAAAATGAGCCTTGCTTTGAAGAACATAAGGTGTCGTGGGGAATTCACAAGTGC
ATTGAATGTCATTCAAAATAAAAGATTTTGAAGCCACTTTTCTATATACATCCACTGCAGTTTTTGC
AAGTACAATACTAACTGTAACAAAGCCTTTGTAATCATATGATGAGCTCTCATAGCAACCATCCAGGT
AAACGGTTTTGTATTTCAAGAAGCATTAGGAACTCTCAGGGGCATTACTCTAGTGTGCCTTAAATGT
GATTTCTAGCTGATTCTCCGGCTTAGATCGTATGGCTAAACACTTAAGTCAACGTAAGTCAACTACT
TGCCAAGTTATAATAGAGAATGTTTCAAAGTACCTCAACTTCTGAACCCACTACTGGCTGCTCATTG
AAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:** □

**ACCN:** NM\_017666

**Insert Size:** 2214 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_017666.4</a>
<b>RefSeq Size:</b>	4695 bp
<b>RefSeq ORF:</b>	2214 bp
<b>Locus ID:</b>	55609
<b>UniProt ID:</b>	<a href="#">Q8ND82</a>
<b>Cytogenetics:</b>	Xq26.1
<b>Domains:</b>	zf-C2H2
<b>MW:</b>	83.1 kDa
<b>Gene Summary:</b>	This gene encodes a member of the zinc finger domain-containing protein family. This family member contains multiple Cys2-His2(C2H2)-type zinc finger domains, the most common type of zinc finger domain that self-folds to form a beta-beta-alpha structure that binds a zinc ion. [provided by RefSeq, Aug 2011]