

Product datasheet for **SC310152**

LRRFIP2 (NM_006309) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRFIP2 (NM_006309) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRFIP2
Synonyms:	HUFI-2
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_006309, the custom clone sequence may differ by one or more nucleotides

```

ATGGGGACTCCTGCTTCTGGAAGGAAAAGAACACCTGTGAAAGACCGATTTTCTGCAGAA
GATGAAGCTTTGAGTAACATTGCCAGAGAGGCAGAGGCAAGGCTGGCAGCAAAACGGGCT
GCCCGGGCAGAAGCAAGAGATATACGCATGAGAGAAGTGGAAACGACAACAAAAAGAGTAC
TCTCTTCATTCTTTGATCGGAAGTGGGGACAGATTCAGAAGTGGCTGGAAGATTGCGAA
AGGGCCAGGTATCCCACCGGTCCAGTCACCATCGTCCTTATCTGGGAGTTGAGGATGCA
TTGTCCATTGCAAGTGTGGCAGTCACAGGTATGATATGTTCAAGGATAGATCATCAAGA
CTTTCATCATTAAATCATTCTTACAGTCACTCTCATGGAATGAAGAAGAGGTCTTCTGAT
TCTCATAAAGACCTACTGAGTGGCCTCTACTTTGACCAGAGAACTATAGCAGTCTTAGA
CATAGCAAACCCACCTCTGCCTACTACACTCGGCAGTCTTCTTCCCTGTACAGTGACCCT
CTGGCAACATATAAGAGTGACAGGGCCTCTCTACTGCAAATTCTGGTCTGCTGAGAAGT
GCCAGTCTGGCATCATTGTACAATGGTGGATTATATAACCCTTATGGTCTCGAACTCCA
TCTGAATGCAGTTATTATTCATCAAGAATAAGTTCAGCCCGAAGCAGTCCAGGGTTTACC
AACGATGACACTGCAAGCATTGTGTCTTCTGATCGTGCCAGTCTGGACGAAGGGAGAGT
GTGGTATCTGCCGCTGATTATTTTCAGTCGCTCCAATCGTAGGGGAAGTGTGTCTCTGAG
GTGGATGATATCAGTATCCCAGATTTGTCCAGCTTGGATGAAAAATCTGCAAAACAGTAT
GCTGAAAATTATACAAGACCTTCATCTCGAAATTTCTGCCTCAGCAACAACCCCTTAAGT
GGAAACTCATCCAGACGAGGAAGTGGGGACACCAGCAGCTTAATAGATCCAGACACTTCA
TTAAGTGAATTGCGGGATATCTATGACCTTAAGGACCAGATACAGGATGTAGAAGGGAGA
TACATGCAGGGGCTTAAAGAACTAAAGGAATCTTTGTCTGAAGTGAAGAAAAATACAAG
AAAGCCATGGTTTCCAATGCACAGTTAGACAATGAGAAGAACAATTTGATCTACCAAGTA
GACACACTCAAGGATGTTATTGAAGAGCAGGAGGAACAGATGGCAGAATTTTATAGAGAA
AATGAAGAAAAATCAAAGGAGTTAGAAAAGCAGAAACATATGTGTAGTGTGCTGCAGCAT
AAGATGGAAGAACTTAAAGAAAGGCTGCGGCAAAGAGATGAGCTTATTGAGGAGAAGCAG
CGCATGCAGCAGAAAAATAGACACCATGACAAAAGAGGTGTTTGACCTCCAGGAGACACTT
CTTTGAAAGATAAAAAAATGGGGCCCTAGAGAAACAGAAAGAATACATTGCCTGCCTT
AGGAATGAGCGAGATATGCTCAGAGAGGAGCTGGCTGACCTGCAGGAGACAGTGAAGACG
GGAGAGAAACATGGCTTAGTTATAATCCCCGATGGCACTCCCAATGGTGTGTCAGTCAT
GAACCAGTGGCTGGAGCCATCACTGTTGTGTCTCAGGAAGCTGCTCAGGCTTGGAGTCA
GCAGGAGAAGGGCCATTAGATGTAAGGCTACGAAAACCTGCTGGAGAGAAGGAAGAACTA
CTGTACAGATTAGAAAACCTGAAGCTTCAGTTAGAGGAGGAACGACAGAAATGCTCCAGG
AATGATGGCACAGTGGGTGACCTGGCAGGACTGCAGAATGGCTCAGACTTGCAGTTCATC
GAAATGCAGAGAGATGCCAATAGACAAATTAGCGAATACAAATTTAAGCTTTCAAAGCA
GAACAGGATATAACTACCTTGGAGCAAAGTATTAGCCGGCTTGGGGACAGGTTCTGAGA
TATAAACTGCTGCTGAGAATGCTGAGAAAAGTTGAAGATGAATTGAAAGCAGAAAAACGG
AAGCTACAACGAGAGTTACGAACAGCACTGGACAAGATTGAGGAGATGGAGATGACCAAC
AGCCACCTGGCCAAGCGGCTGGAGAAGATGAAGGCAATAGGACAGCACTTCTGGCCAG
CAGTAG

```

Restriction Sites: Please inquire

ACCN: NM_006309

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_006309.1</u> , <u>NP_006300.1</u>
RefSeq Size:	2412 bp
RefSeq ORF:	2166 bp
Locus ID:	9209
UniProt ID:	<u>Q9Y608</u>
Cytogenetics:	3p22.2
Gene Summary:	<p>The protein encoded by this gene, along with MYD88, binds to the cytosolic tail of toll-like receptor 4 (TLR4), which results in activation of nuclear factor kappa B signaling. The ubiquitin-like protein FAT10 prevents the interaction of the encoded protein and TLR4, thereby inactivating the nuclear factor kappa B signaling pathway. In addition, this protein can downregulate the NLRP3 inflammasome by recruiting the caspase-1 inhibitor Flightless-I to the inflammasome complex. [provided by RefSeq, Jan 2017]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a, also known as 1).</p>