

Product datasheet for **SC310157**

LRCH3 (NM_032773) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRCH3 (NM_032773) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRCH3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310157 representing NM_032773.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCGGCCCGGGCTTGGTCGCTGTGGCAGCGGCTCCGAGTACTCTGGCACGGTAGCGTCGGGAGGT
AACCTCCCTGGTGTTCCTGCTGCGGCCAAGCTCCGGGCGAGGCCCTGGTTTTGGCCGGGCTCGTGGAGC
CGCTCTCTCGATCGAGCCCTGGAGGAGCGCGGCTCACTGGGGTCTGAGCCTGAGCGGCCGAAACTG
AGGGAGTTTCCCGGGGAGCGGCCAACACGACCTGACGACACCCCGGGCGGACCTGTCGCGAAAT
CGCCTTTCAGAAATTCCTATAGAAGCATGTCACCTTGTCTCTGGAAAACTCAACTTGTACAAAAT
TGTATTCGTTATATTCCAGAGGCAATTTAAACCTACAAGCTCTAACATTCTTAAATATTAGTCGGAAC
CAACTGTCAACATTGCCGTACACTTGTGTAATTTGCCATTGAAAGTCTAATTGCTAGTAATAACAAA
TTGGTGTCACTCCAGAAGAAATGGACACCTTAGACATTTGATGGAACCTGATGTGAGCTGCAATGAA
ATTCAAACTATACCTCCCAAATGGTAACCTGGAGGCCCTGAGAGACCTAATGTAAGAAGAAATCAC
CTAGTACATTTGCCTGAAGAGCTGGCGGAGTTGCCTTTGATACGGTTAGACTTCTCATGCAATAAAAT
ACCACAATCCCTGTTTGTATCGGAACCTCAGGCACCTACAGACGATCACCTAGATAACAAATCCACTA
CAATCACCTCCTGCACAGATATGTATAAAAGGCAAAGTCCACATATTTAAATACCTGAACATACAAGCT
TGTAAGATTGCTCCAGATCTGCCGATTATGATAGGAGACCGTTGGGTTTTGGCTCCTGCCATGAAGAA
CTGTACTCAAGTCGCCCTTATGGAGCCCTTGATTCAGGCTTCAATAGTGTGGACAGTGGTGATAAGAGA
TGGTCAGGGAATGAACCTACAGATGAATTTTTCAGATCTGCCTCTTCGAGTAGCAGAGATTACTAAAGAA
CAAAGACTACGAAGAGAAAGCCAGTACCAAGAGAACCGCGGCAGTTTGGTAGTAACAAACGGCGGAGTG
GAACATGATCTGGATCAGATTGACTACATAGACAGCTGCACCGCAGAGGAAGAGGAGGCCGAGGTGAGA
CAGCCCAAGGGACCAGACCAGACGCTTAGTTACAGTTTATGGCGTATATTGAACAGCGCGAATC
TCTCATGAGGTTTACCAGTAAAGCCAGTAGCCATTAGGGAGTTTCAAAAAACAGAAGATATGAGAAGA
TATTTACATCAAAACAGGGTTCAGCTGAGCCATCTTCCCTCCTGCACTATCAGCAAGTACAATCAG
CTGTCACACACAGACCTGGAACCTCATCAGAGAAGGGAGCAGTTAGTAGAGCGCACTCGGAGAGAGGCT
CAGCTTGTGCCCTGCAGTATGAGGAGGAGAAAAAAGGACCAAGCAGATCCAGAGAGATGCTGTCTG
GACTTTGTCAAACAAAAAGCATCACAAGTCCAAAAACAGCACCCGCTCCTAGATGGCGTAGATGGT
GAGTGGCCCTTCCCATCCAGAAGGTCTCAGCACACTGATGATAGTGCCTTGTGCATGTGCTGTCAGGG
TTGAATCAAGTGGGCTGTGCTGCTACCTGCCTCATTCTTCTGCCTTACGCCTCTTAAGAGTATGAC
AGACCTAATGCTCTATTAAGTTCACCTGCAACAGAAACAGTTCATCATTCCCTGCATATCTTTTCT
GCTGCTATCCAGAGAAATCAGCCTCAGCGCCCTGAAAGCTTCTTTTCCGAGCAGGTGTCAGGGCAGAA
ACCAACAAAGGTCATGCTTACCCCTTCTCCATCTGCTGCACCTACCACTGATTCTACAGATCCATA
ACAGGACAGAAATCAAGACAGAGAGAAGAAGAGCTGGAATTAATAGACCAACTGCGTAAACATATTGAG
TACCGGTTGAAAGTGTCTTACCTTGTGATCTCGGAGCAGCTCTAACTGACGGTGTGTTCTTTGCCAT
TTGGCCAATCATGTGCGACCTCGATCTGTCCCAAGCATTATGTTCCCTCACCAGCTGTAGTAAGTGA

```

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: SgfI-MluI

ACCN: NM_032773

Insert Size: 2139 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032773.3
RefSeq Size:	2325 bp
RefSeq ORF:	2139 bp
Locus ID:	84859
UniProt ID:	Q96II8
Cytogenetics:	3q29
Domains:	LRR
Protein Families:	Secreted Protein
MW:	78.8 kDa
Gene Summary:	<p>As part of the DISP complex, may regulate the association of septins with actin and thereby regulate the actin cytoskeleton.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 3' UTR and coding sequence compared to variant 1. The resulting isoform (2) has a shorter and distinct C-terminus compared to isoform 1.</p>