

## Product datasheet for **SC310163**

### **ApoER2 (LRP8) (NM\_017522) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	ApoER2 (LRP8) (NM_017522) Human Untagged Clone
Tag:	Tag Free
Symbol:	ApoER2
Synonyms:	APOER2; HSZ75190; LRP-8; MCI1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)



**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** The ORF of this clone has been fully sequenced and found to be a perfect match to NM\_017522.3.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_017522.3](#), [NP\\_059992.3](#)

**RefSeq Size:** 3739 bp

**RefSeq ORF:** 2103 bp

**Locus ID:** 7804

**UniProt ID:** [Q14114](#)

**Cytogenetics:** 1p32.3

**Domains:** ldl\_recept\_b, EGF\_CA, ldl\_recept\_a, EGF

**Protein Families:** Druggable Genome, Secreted Protein, Transmembrane

**Gene Summary:**

This gene encodes a member of the low density lipoprotein receptor (LDLR) family. Low density lipoprotein receptors are cell surface proteins that play roles in both signal transduction and receptor-mediated endocytosis of specific ligands for lysosomal degradation. The encoded protein plays a critical role in the migration of neurons during development by mediating Reelin signaling, and also functions as a receptor for the cholesterol transport protein apolipoprotein E. Expression of this gene may be a marker for major depressive disorder. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2011]

Transcript Variant: This variant (3) lacks three exons in the coding region but maintains the reading frame, compared to variant 1. The encoded isoform (3) is shorter than isoform 1.

Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.