

## Product datasheet for **SC310172**

### **TBC1D23 (NM\_018309) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	TBC1D23 (NM_018309) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBC1D23
Synonyms:	NS4ATP1; PCH11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC310172 representing NM\_018309.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGAAGGAGAAGATGTGCCCGCGCTGCCAACGTGAGCGGGCAGCGCTGGGAAAAAGATCTTGAA
GAAGCTCTGGAAGCAGGAGGTTGTGATCTTGAAACGTTGAGAAAATAAATCAAGGAAGACCGCTGCCT
GCTGATCTGAGGGCCAAAGTTTGAAGATTGCTCTGAATGTTGCAGGAAAAGGTGATAGTTTGGCATCA
TGGGATGGTATTTTAGACTTGCCAGAACAGAACACTATTCACAAAAGATTGCCTGCAGTTTATTGACCAG
CTTTCAGTGCCAGAGGAGAAGGCAGCAGAATTACTTTTGGATATTGAATCTGTAATTACCTTTTATTGT
AAATCACGTAACATTAATATAGCACATCCCTTAGCTGGATACATCTACTGAAACCATTGGTGCATCTT
CAACTGCCACGCAGGATTTATACAACGCTTTTATGCCATAATGAATAAGTACATCCCAGGGATTGT
TCCCAGAAAGGGAGACCATTTTCATCTCTCAGGTTGCTCATCCAATACCATGAGCCTGAGCTTTGTTCT
TATCTTGATACAAAGAAAATTACTCCAGACTCCTATGCACTCACTGGCTTGAAGTCTTTTTCATGT
TACTGTTCCACTGAAGTCACTCAGGCAATATGGGATGGATATCTACAACAAGCAGATCCATTTTTTATT
TATTTCTTAATGTTAATTATCCTTGTTAATGCAAAAGAAGTTATTTTAAACAAGAGTCAGACAGCAAA
GAAGAAGTTATCAAGTCTTGAAAAACTCCATCCAGTCTGAATATAGAAGATATAGAAGACCTTTTC
TCTCTGGCTCAGTATTATTGCAGCAAAACCCGGCTTCTTTTAGGAAGGATAATCACCATCTCTTGGT
AGTACTTTGTTGGGAATTAAGGATGATGATGCAGATCTGAGTCAGGCTCTTTGTCTGGCCATCTCCGTG
TCAGAGATCCTTCAAGCGAATCAGCTACAAGGGGAAGGAGTCCGGTCTTTGTGGTGGATTGCCGTCT
GCAGAACAATAAATGCTGGGCATTTATCAACTGCTTCCACTTAGATTCAGACCTGATGCTTCAGAA
CCATCTGAGTTTGACAGTCAGTAAAATCCTTGCTGGAAGCACAGAAGCAGTCCATTGAGTCTGGCTCC
ATAGCTGGTGGGAGCACCTCTGTTTTATGGGCAGTGGCAGGGAGGAAGAAGACATGTATATGAACATG
GTCTGGCACACTTTTTACAGAAAAACAAGAATATGTGAGTATTGCCAGTGGAGGATTTATGGCACTG
CAGCAGCACCTGGCAGACATTAATGTGGATGGACCAGAAAATGGATATGGCCATTGGATTGCTAGTACC
TCAGGCTCAAGGAGCAGTATCAATTCTGTTGATGGTGAATCTCCTAATGGCTCAAGTATAGAGGAATG
AAATCACTAGTAAATAAAATGACTGTGGCTTTGAAGACAAAATCCGTTAATGTCAGGGAAAAAGTTATC
AGTTTTATAGAGAATACATCAACTCCTGTGGATCGGCATGTGAGCAGCAGTGACAGAGTGGCAAGCCT
TACCGTGGCGTAAAGCCTGTTTTCAGCATTGGGGATGAAGAAGAATACGACACAGATGAAATTGACAGT
TCTTCAATGTCAGATGATGATAGAAAAGAGTTGTAACATTGACTTGGATAAAACAAACCAGATGTC
AAACATCATTTTCTTGTAAAGAAGTAAAGAAAGTGGACACATGTTTCCAGTCACTGTTGGTTACT
GCAACACATATGTACTGTTTAAAGGAGATTGTTTACGGAAAGGATTGGCTTATATACAGTCTCGACAA
GGCTGAATTCTGTAGTTAAAATTACATCCAAAAAAAACATCCTGAACTCATTACCTTCAAGTATGGA
AATAGCAGTGTTCAGGAATAGAAATCTTGCAATCGAAAGGATTTGATTCCAATGCAGGGGATGCA
ACTAAAGCCATAAAACAGCAGATCATGAAAGTTTTGGATGCTTTGGAAAGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**

**ACCN:** NM\_018309

**Insert Size:** 2055 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_018309.4</a>
<b>RefSeq Size:</b>	3841 bp
<b>RefSeq ORF:</b>	2055 bp
<b>Locus ID:</b>	55773
<b>UniProt ID:</b>	<a href="#">Q9NUY8</a>
<b>Cytogenetics:</b>	3q12.1-q12.2
<b>Domains:</b>	TBC
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	76.5 kDa
<b>Gene Summary:</b>	<p>Putative Rab GTPase-activating protein which plays a role in vesicular trafficking (PubMed:28823707). Involved in endosome-to-Golgi trafficking. Acts as a bridging protein by binding simultaneously to golgins, including GOLGA1 and GOLGA4, located at the trans-Golgi, and to the WASH complex, located on endosome-derived vesicles (PubMed:29084197, PubMed:29426865). Together with WDR11 complex facilitates the golgin-mediated capture of vesicles generated using AP-1 (PubMed:29426865). Plays a role in brain development, including in cortical neuron positioning (By similarity). May also be important for neurite outgrowth, possibly through its involvement in membrane trafficking and cargo delivery, 2 processes that are essential for axonal and dendritic growth (By similarity). May act as a general inhibitor of innate immunity signaling, strongly inhibiting multiple TLR and dectin/CLEC7A-signaling pathways. Does not alter initial activation events, but instead affects maintenance of inflammatory gene expression several hours after bacterial lipopolysaccharide (LPS) challenge (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) lacks an exon in the 3' coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (2) is shorter than isoform 1.</p>