

## Product datasheet for **SC310184**

### AGAP8 (AGAP4) (NM\_133446) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	AGAP8 (AGAP4) (NM_133446) Human Untagged Clone
Tag:	Tag Free
Symbol:	AGAP4
Synonyms:	AGAP-4; AGAP-8; AGAP8; CTGLF1; CTGLF5; MRIP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC310184 representing NM\_133446.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGAACATACTGACCTGTCGTGTGCACCCTAGCGTCAGCCTCGAGTTTGACCAGCAGCAGGGGTCCG
GTGTGTCCCTCTGAATCTGAGATCTATGAGGCAGGAGCTGGGGACAGGATGGCAGGAGCGCCCATGGCT
GCTGCTGTACAGCCTGCTGAGGTGACTGTTGAAGTTGGTGAGGACCTCCACATGCACCACGTTTCGTGAC
CGGGAGATGCCTGAAGCTTTGGAGTTTAAACCTTCTGCCAATCCAGAGGCAAGCACAATATTCCAGAGG
AACTCTCAAACAGATGTTGTAGAAATAAGAAGAAGCAACTGTACAAACCATGTATCTACTGTGCGTTTC
AGTCAACAATACAGCTTGTGTTGACAATATTCCTTGATGACAGCACAGCCATCCAGCATTATCTTACA
ATGACAATAATATCTGTGACCTGGAGATACCTCATCATATCACACAAAGAGATGCAGATAGATCTTTG
AGCATACCTGATGAACAGTTACTCATTTGCGGTTTCCACCGTGCACATTATGAAGAAAAGAAATGGA
GGTGGGAGTTAAATAACTATTCCTCCTCCATTCCATCGACTCCCAGCACCAGCCAGGAGGACCCCTCAG
TTCAGTGTTCCTCCACTGCCAACACACCCACCCCGTTTGAAGCGGTCCATGCGCTGGTCCAACCTG
TTTACATCTGAGAAAGGGAGTGACCCAGACAAAGAGAGGAAAGCCCGGAGAATCATGCTGACACCATC
GGGAGCGGCAGAGCCATCCCATTAACAGGGCATGCTCTTAAAGCGAAGTGGGAAATGGCTGAAGACA
TGGAAAAAGAAATACGTCACCCTGTGTTCCAATGGCGTGCTCACCTATTATTCAAGCTTAGGTGATTAT
ATGAAGAATATTCATAAAAAAGAGATTGACCTTCAGACATCTACCATCAAAGTCCAGGAAAGTGGCCA
TCCCTAGCCACATCGGCCTGCACACCCATCTCCACCTCTAAAAGCAATGGCCTATCCAAGGACATGGAC
ACCGGGCTGGGTGACTCCATATGCTTCAGCCCCAGTATCTCCAGCACCACCAGCCCCAAGCTCAACCCG
CCCCCTCTCCTCATGCCAATAAAAAGAAACACCTAAAGAAGAAAAGCACCAACAACCTTTATGATTGTG
TCTGCCACTGGCAAACGTGGCACTTTGAAGCCACGAGTATGAGGAGCGGGATGCCTGGTCCAAGCC
ATCCAGAGCCAGATCCTGGCCAGCTGCGATCATGCGAGAGCAGTAAAAGCAAGTCCCAGCTGACCAGC
CAGAGCAAGGCCATGGCCCTGCAGTCGATCCTAAAACATGCGTGGGAACGCCCACTGTGTGGACTGTGAG
ACCCAGAATCCTAAGTGGCCAGTTTGAACCTGGGAGTCTCATGTGATTGAATGCTCAGGAATCCAC
CGCAGTCTTGGCACCCGCTTTCCCGTGTGCGATCTCTGGAGCTGGATGACTGGCCAGTTGAGCTCAGG
AAGTTATGTCATCTATTGGCAATGACCTAGCCAACAGCATCTGGGAAGGGAGCAGCCAGGGGCAGACA
AAACCCTCAGAAAAGTCCACGAGGGAAGAGAAGGAACGGTGGATCCGTTCCAAATATGAGGAGAAGCTC
TTTCTGGCCCCACTACCCTGCAGTGTGCTGCTGGCCAGCAGCTGCTGCGGGCCACCCTGATGAG
GACCTGCAGACAGCCATCCTGCTGCTGGCAGATGGCTCCCGTGGAGGTTGAACGAGACCTGTGGGGAG
GGAGACGGCTGCACGGCGCTCCATCTGGCTGCCGAAGGGAAATGTGGTCTGGCGCAGCTCCTGATC
TGGTACGGGGTGGACGTATGGCCCAGATGCCACGGGAACACAGCGCTGACCTACGCCCGCAGGCC
TCCAGCCAGGAGTGCATCAACGTGCTTCTGCACTACGGCTGCCCCGACAAGTGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:** □

**ACCN:** NM\_133446

**Insert Size:** 1992 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_133446.3</a>
<b>RefSeq Size:</b>	2451 bp
<b>RefSeq ORF:</b>	1992 bp
<b>Locus ID:</b>	119016
<b>UniProt ID:</b>	<a href="#">Q96P64</a>
<b>Cytogenetics:</b>	10q11.22
<b>Domains:</b>	ArfGap, PH, ANK
<b>MW:</b>	73.1 kDa
<b>Gene Summary:</b>	Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) lacks an exon in the 5' coding region, compared to variant 1. The encoded isoform (2) is shorter, compared to isoform 1.