

## Product datasheet for **SC310190**

### **RUFY2 (NM\_017987) Human Untagged Clone**

#### **Product data:**

|                           |                                        |
|---------------------------|----------------------------------------|
| Product Type:             | Expression Plasmids                    |
| Product Name:             | RUFY2 (NM_017987) Human Untagged Clone |
| Tag:                      | Tag Free                               |
| Symbol:                   | RUFY2                                  |
| Synonyms:                 | RABIP4R; ZFYVE13                       |
| Mammalian Cell Selection: | Neomycin                               |
| Vector:                   | pCMV6-Entry (PS100001)                 |
| E. coli Selection:        | Kanamycin (25 ug/mL)                   |



[View online »](#)

**Fully Sequenced ORF:** >SC310190 representing NM\_017987.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGACTTTTCAGGTTTGGGGTGGCGGAGGGAAGACGCATCTCAGGTTTGGCCTGGGTGCCTGACGCC
GAGGGCGGGAGACGGGAATGCTGACGCGCAGGAGCCTGGCTACAAAAGACCCACAGCTGTAGAGAGA
GCAAACTTGTTAAACATGGCTAACTGAGTATCAAAGGACTCATTGAATCTGCTCTGAGCTTTGGCCGC
ACTTTGGATTCTGACTATCCCCCTTGCAAGCAATCTTTGTTGTTATGGAACATTGCCTGAAACACGGT
CTTAAAGTAAGAAAATCATTTTTGTGTTACAACAAAACCATCTGGGGCCCTTTGGAAGTGGTGGAGAAG
CTGTACCCCAAGCAGAGGAAATAGGAGCTAGTGTCCGGGATCTACCTGGTCTGAAGACCCCTCTGGGT
CGAGCAAGAGCGTGGCTTCGATTAGCCCTCATGCAAAAAAATGGCCGATTACTTACGTTGCTTAATT
ATTCAGAGGGATCTCTTGAGTGAGTTTATGAGTATCACGCACTAATGATGGAAGAAGAAGGAGCAGTA
ATTGTTGGGCTGCTGGTGGCCTGAATGTGATCGATGCTAATCTGTGTGAAGGAGAGGATTTAGAC
TCACAAGTTGGAGTGATTGATTTTCTATGTATTTAAAGAATGAAGAAGATATTGAAATAAAGAAAGG
AATGTTCAAATTGCTGCCATATTAGACCAAAAGAAATTATGTTGAAGAATTAATAGACAACCTGAACAGC
ACAGTCAGCAGCCTCCATTCAAGAGTTGATTCTATTAGAAAAGTCAAACTAAGCTGATTGAAGAGTTA
GCAATAGCAAAGAATAACATCATTAACTCCAGGAAGAAAATCATCAATTACGAAGTGAATAAATTG
ATTTTAAATGAAAACACAGCAGCAGCTAGAGGTTACCAAGTAGATGTGGAAGTGAAGCTTCAACATAT
AAGCATTTCTCGTCAGGGGCTAGATGAAATGTACAATGAAGCCAGAAGGCAGCTTCGAGATGAATCTCAG
TTACGACAGGATGTAGAGAATGAGCTAGCAGTACAAGTTAGTATGAAGCATGAGATTGAAGTGGCATG
AAGTTGCTGGAGAAAGATATCCATGAGAAACAAGATACTCTGATAGGCCCTTCGACAACAAGTGAAGGAA
GTTAAAGCAATTAACATAGAGATGTATCAAAAGTTGCAGGTTCTGAAGATGGCTTGAAAGAAAAAAT
GAAATAATTGCCGACTAGAAGAAAAAACCAATAAAATTACTGCAGCCATGAGGCAGCTGGAACAAAGA
TTGCAGCAAGCAGAGAAGGCGCAATGGAAGCTGAAGATGAGGATGAGAAATATCTACAAGAATGTCTC
AGTAAATCTGATAGTCTGCAGAAACAAATCTCCAAAAGGAGAAACAGCTGGTGAAGTGAAGTGAAGTGA
TTGAAGATTGAGAAGGAATGGAGGCAGACTTTGCAGGAAGATCTTCAAAAGGAGAAAGATGCCTTATCT
CATCTTAGAAATGAGACTCAACAAATCATTAGTCTTAAAAAGAGTTCCTTAACCTCCAGGATGAAAT
CAGCAGTTGAAAAAATATATCATGAACAAGAGCAAGCTCTTCAAGAACTCGGCAACAAGCTTAGCGAA
TCAAACTTAAATTTGAAGACATAAAAGAAGCCAAAGCATTGCAGGAGTGGTTGGCTGAAAGAC
AAAGAAGCAACACATTGTAACTTTGTGAAAAGGAATTCTCACTCTTAAGAGAAAGCACCAGTGTAGA
AATTGTGGGAAATTTCTGTAATGCCTGCTCTGACAACGAAGTACCTTTGCCTTCTTCAACAAAACCA
GTACGGGTTTGTGATTCTGTGTCATGCACTGCTCATTGAGAGATGCTCATCTAAGTTCGCTTGA
ACGGTACGGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-MluI

**ACCN:** NM\_017987

**Insert Size:** 1926 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_017987.4](#)

**RefSeq Size:** 4520 bp

**RefSeq ORF:** 1926 bp

**Locus ID:** 55680

**UniProt ID:** [Q8WXA3](#)

**Cytogenetics:** 10q21.3

**Domains:** FYVE, RUN

**MW:** 74.1 kDa