

Product datasheet for **SC310191**

Rho GTPase activating protein 24 (ARHGAP24) (NM_031305) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rho GTPase activating protein 24 (ARHGAP24) (NM_031305) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP24
Synonyms:	FILGAP; p73; p73RhoGAP; RC-GAP72; RCGAP72
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310191 representing NM_031305.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCTGAAGACCGGAATTCGGGGGTGCCCGCTGGTGCCTAGCCTCAACTCCTTTCATCCCTAAA
ACTACATACAGAAGAATCAAACGGTGTGTTAGTTTTTCGAAAGGCATTTTTGGACAGAACTGGAGGAT
ACTGTTTCGTTATGAGAAGAGATATGGGAACCGTCTGGCTCCGATGTTGGTGGAGCAGTGCCTGGACTTT
ATCCGACAAAGGGGGCTGAAAGAAGAGGGTCTCTTTCGACTGCCAGGCCAGGCTAATCTTGTAAAGGAG
CTCCAAGATGCCTTGTACTGTGGGGAGAAGCCATCATTGACAGCAACACAGATGTACACACGGTGGCA
TCACTTCTTAAGCTGTACCTCCGAGAATTCCAGAACCAGTTATTCCTTATGCGAAGTATGAAGATTTT
TTGTCATGTGCCAACTGCTCAGCAAGGAAGAGGAAGCAGGTGTTAAGGAATTAGCAAAGCAGGTGAAG
AGTTTGCCAGTGGTAAATTACAACCTCCTCAAGTATATTTGCAGATTCTTGGATGAAGTACAGTCTAC
TCGGGAGTTAACAAAATGAGTGTGCAGAACTTGGCAACGGTCTTTGGTCTAATATCCTGCGCCCAAAA
GTGAAGATCCTTTGACTATCATGGAGGGCACTGTGGTGGTCCAGCAGTTGATGTCAGTGTGATTAGC
AAACATGATTGCCTCTTCCCAAAGATGCAGAACTACAAAGCAAGCCCAAGATGGAGTGAGCAACAAAC
AATGAAATTCAGAAGAAAGCCACCATGGGGCAGTTACAGAACAAGGAGAAACAATAACACCAAGGACAGC
CCTAGTAGGCAGTGTCTCTGGGACAAGTCTGAGTACCCCAAGAGAAGCAGCATGAACAATGGATCCCCC
ACAGCTCTATCAGGCAGCAAAACCAACAGCCCAAGAAGCAGTGTTCACAAGCTAGATGTGTCTAGAAGC
CCCCCTCATGGTCAAAAAGAACCAGCCTTTAATAAGGGTAGTGGGATAGTTACCAATGGGTCCTTC
AGCAGCAGTAATGCAGAAGGTCTTGAGAAAACCAACACCCCAATGGGAGCCTACAGGCCAGAAGG
AGCTCTTCACTGAAGGTATCTGGTACCAAAATGGGCACGCACAGTGTACAGAATGGAACGGTGCATG
GGCATTTTGAACAGCGACACACTCGGGAACCCCAAAATGTTGAAACATGAGCTGGCTGCCAAATGGC
TATGTGACCCTGAGGGATAACAAGCAGAAAGAACAAGCTGGAGAGTTAGGCCAGCACAACAGACTGTCC
ACCTATGATAATGTCCATCAACAGTTCTCCATGATGAACCTTGTGACAAGCAGAGCATTGACAGTGT
ACCTGGTCCACTTCTCTGTGAAATCTCCCTCCCTGAGAACTCCAACCTCTGTGCTCTTCTACCACC
ACCTGCCAGAGCAAGACTTTTTTGGGGGAACTTTGAGGACCCTGTTTTGGATGGGCCCCCGCAGGAC
GACCTTTCCACCCAGGGACTATGAAAGCAAAAGTGACCACAGGAGTGTGGGAGTTCGAAGTAGTCGT
GCCACCAGTAGCAGTGACAACAGTGAGACATTTGTGGCAACAGCAGCAGCAACCACAGTGCCTGCAC
AGTTTAGTTTCCAGCCTGAAACAGGAAATGACCAACAGAAGATAGAGTATGAGTCCAGGATAAAGAGC
TTAGAACAGCGAAACTTGACTTTGAAACAGAAATGATGAGCCTCCATGATGAACTGGATCAGGAGAGG
AAAAAGTTCACAATGATAGAAATAAAAATGCGAAATGCCGAGCGAGCAAAAGAAGATGCCGAGAAAAGA
AATGACATGCTACAGAAAGAAATGGAGCAGTTTTTTTCCAGTGTGGAGAACTGACAGTGAACCCAGG
AGAACCGAGAGAGGAAACACAATATGGATTCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_031305

Insert Size: 1968 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031305.2
RefSeq Size:	4574 bp
RefSeq ORF:	1968 bp
Locus ID:	83478
UniProt ID:	Q8N264
Cytogenetics:	4q21.23-q21.3
Domains:	RhoGAP
MW:	73.4 kDa
Gene Summary:	<p>This gene encodes a Rho-GTPase activating protein, which is specific for the small GTPase family member Rac. Binding of the encoded protein by filamin A targets it to sites of membrane protrusion, where it antagonizes Rac. This results in suppression of lamellae formation and promotion of retraction to regulate cell polarity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (2) has a shorter and distinct N-terminus compared to isoform 1.</p>