

Product datasheet for **SC310233**

CCDC148 (NM_138803) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC148 (NM_138803) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC148
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310233 representing NM_138803.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTGTGCAGCTTCTGCTTCTCCAGATAATCTGGTATTCCATATGAAAAATGAGATGAGAAACATCAAG
TACAAACCACTAGACTATCAACAATTGCGTGCATTAAGTGAAGCAAAGAAATGGCTTCTGCCTCTGCA
AAGCTAAAGATCAGAAAAGCAATGTTGACTTCAAAGTTATCCAAAGAACAACACTAATAAAACAACAC
AAGCAAGTGTGGTGGCAGGAATACCAGAGGCTGAATGAAGTCAGATGTAAATGGAATCTGAAATAAAA
TCCCTTCTCAATGAAGAGAACATTGGAATGAATGTCTTTGTGACCTTACAAATTTTGTGCAAGAGCTA
TCAGAACAACAGTGCACATATCTTAAAAATGTAATAATCCTATTCAGCAGCTGAGAGCAGATCTAAAA
TACAGACAGCATCACACTTTGCAGCATTACACCCACATATTGAGTTAACTCCATGAAAGTATTGGAA
GAGGTTGACTTTGTGAAGAAACAATTGAAACTGTCTTTGAAAGACTTAGGCTGGAGCAACAGAGAATA
GAAATGATCTTTCAGACTGGAGTATAAAGATTCTAGATCATTCTTTGGAAGAAAAGACTAACCCGCTT
AGTGAAGTCCCATTTGAATTAGAAAGTTTGAATGCCCATACCCTGATTTGAAATCTTCAATTCTCAGT
GAGTTCTATAAGTTTACACAGAAATATCAGAAGAAGCTTCAAGACTTTAATCTGCAGTTGGAAGACATA
TACAGAACTGTCAACTAAGTGAAGAAGACCACTGGATTACCAGGCTATTTGGATCAGTACCCTGGA
GATCTCTTTGGAAGAAGGACTCTGTATCTGGACATGTTACAAAGATATTTTCTCACAATCTAGGCAT
GATTTGGTTGAACACGAGAAATATTGTGACCAATATCGCTTGTCTATAGAGCAGCAAAATATCCTGATA
TCAAAATGGAATAAAAAAAGAAAGACTTTATACAGAAGGCTGTGCTGACACTCACTGAGGCTTGTGCA
ACACATGAAATGGAGAGCATGTTGGCTAAGGACAAGAAAAAGCAACAGGAATTGTGTCTGATCTGAAA
GCCAAGGTTCTCAATGGCGAGCCCAAGAAGAGGTAGCAAGACTGGAATGGAATTTCTGCCAGA
AGAAGAGAAAAGGAAGAGGAGAGAAAGAGAACTGTGGAAGAAGAAGGAATTGTTGCAAAGAGCAGAGAAG
AAAAAAAATAAAAAAATACTGGGCCAAGAAAAACAGAAGTGGCAAGAAATGGAATGAGAGATCTT
CAGCGTCTAGAAGAACTGAAGAAATTAATCGCTGAACAGTCACTAAAAGACAGAGAAAGGTTAAATAT
AGACAAGAATTATTGAAAGGCGTTTGTGAGAGAAAAGGAAGTGGCCCTTCAAGAAGCTCATGAAGAC
AAAGAGAGAGCACGGCGCTAGAAGCCCTTAGGAAACAGGTTGCTGTTGTTGCTCAATTTGATCCTGTT
AGAATGATGTCAGATACAATGGCATCAAAGCTAGAATGGGCATTGAAATGAAGAAGAATTTATCTT
CAAAAGCCACTCTTCACATTGAATACATACAATGAACAGCAGATAATTTCTGACCCTAGACTTCGCTTC
GAGTTAGCACTTCGAGAAGCTGGACTTCATAGAACACTTTATGCTAAAGAGATATTACAAAAATTAGT
CCTCAAAACCTCCAAGAAAGGACATGGAGTCTACTGTATTTAAATAAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_138803

Insert Size: 1776 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_138803.3](#)

RefSeq Size: 2846 bp

RefSeq ORF: 1776 bp

Locus ID: 130940

UniProt ID: [Q8NFR7](#)

Cytogenetics: 2q24.1

MW: 71.1 kDa