

Product datasheet for **SC310250**

MAP11 (NM_018275) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAP11 (NM_018275) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAP11
Synonyms:	C7orf43; MCPH25; TRAPPC14
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_018275, the custom clone sequence may differ by one or more nucleotides <pre> ATGGAGTCCCAGTGCGACTACTCGATGTACTTCCCGGCCGTGCCGCTGCCGCCGCGCGCG GAGCTGGCAGGGGATCCGGGCCGGTACCGGGCGCTGCCCGGCCGCAACCATCTGTACTTG GGGAGACTGTCCGTTTTCTGCTGGTGTGCGCTGCCGGGGCGGTGCGGGGTCCGGCACC GGGGCGGCCCGGGCTTGGGCTCCAGAGGAGCCTGGGCAGAACTGGCAACCGCCCTGGCC GCCCTGGCCTCGGTCAGCGCCGGAGGCGGGATGCCGGGGGCGGCGCGCCGCGACCAAG GATTCGAGAGCCCCAGGGGGAGGGGATCCTGGGGGTGGGGTTTTGTTCCGAGGCTGCAGC CCCCTTCTCACCCACGGCCCGGGCCCTGCTACCTCAGGGGGAGCGACCACGCTGCCTGTG GAGGAACCGATTGTGTCCACAGATGAGGTCATCTTCCCACTACCGTTTCACTGGATAGA CTGCCCCCAGGGACACCTAAGGCCAAGATTGTAGTGACTGTGTGAAGCGGGAGATTGAG GCACCAGAGGTCAAGATCAAGGCTACCTGCGATTGCTGCAGACCCGATCTCCTGGGGAG ACATTCGGGGGCGAGCAGAGCGCTTTCAAGGCCAAGTGAGCACGCTGCTGACTCTGCTG CCCCCTCCGTTCTGAGATGCCGGCAGTTCACTGTGGTGGAAAACACTTGACCGTGCTC AAGGTGTGAACAGCTCCTCTCAGGAGGAAATCTCCATCTGGGATATCCGAATCCTCCCC AACTTCAACGCCAGTTATCTACCTGTGATGCCCGATGGCTCTGTGCTGCTGGTGGACAAT GTCTGTACACAGTCTGGGGAAGTCTCCATGGGCTCCTTCTGCCGCTACCCGGGACCTCT GGCTGCTTCCCCTGCCGCTGAATGCCCTGGAGGAACACAACCTTCTGTTTCAGCTGAGA GGGGGTGAGCAGCCCCCTCCAGGGGCCAAGGAGGGCCTGGAAGTCCCCTGATTGCTGTG GTTCACTGGTCTACCCCAAAGCTGCCCTTCACTCAGAGCATCTACACCCACTACCGCCTG CCCAGTGTCCGTTGGACCGCCCGTGTGTTGTGATGACCGCTTCTGTAAGTCCCCTGTT CGGACCTACGAGCGTTTCACTGTACCTACACGCTGCTTAACAATCTCCAAGACTTCCTT GCTGTGAGGCTCGTGTGGACCCAGAGCATGCACAGGCTGGAAAGCAGCTGTGTGAGGAG GAGCGCCGGGCCATGCAGGCTGCCCTGGACTCCGTCGTCTGCCACACGCCCTCAACAAC CTTGGCTTTTCCCGGAAGGGCAGCGCGCTCACCTCAGTGTGGCCTTCCAGGCTCTGAGG ACGGGGCTCTTCGAGCTAAGCCAGCACATGAAACTGAAGCTGCAGTTACCGCCAGCGTG TCCACCCCTCCACCCAGGCCCCGCCCTCTCCCGCAAGAGCAGCCCAGCAGCCCTGCT GTCCGGGACTTGGTGGAGAGGCATCAGGCTAGCCTGGGCCGCTCCAGTCTTCTCCAC CAGCAGCCTTCCCGAAGCCACCTCATGAGGTCGGGCAGTGTGATGGAGCGCAGAGCCATC ACGCCCCCTGTGGCCTCTCCTGTTGGCCGCCCTCTACCTGCCCGGACAAGGCTGTG TTGTCTCTGGACAAGATTGCCAAGCGCAGTGCAAGTCTGGTGGTGGAAACCGTCAAG TAG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_018275
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_018275.3, NP_060745.3</u>
RefSeq Size:	2563 bp
RefSeq ORF:	1743 bp
Locus ID:	55262
UniProt ID:	<u>Q8WVR3</u>
Cytogenetics:	7q22.1
Gene Summary:	Microtubule associated protein wich seems to play a role in spindle dynamics and cell division (PubMed:30715179, PubMed:30447097). Modulates YAP1 activity as transcriptional regulator (PubMed:30447097).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).