

## Product datasheet for SC310352

### GABRR1 (NM\_002042) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	GABRR1 (NM_002042) Human Untagged Clone
Tag:	Tag Free
Symbol:	GABRR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310352 representing NM_002042. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTGGCTGTCCAAATATGAGATTTGGCATCTTTCTTTTGTGGTGGGGATGGGTTTTGGCCACTGAA
AGCAGAAATGCACTGGCCCGGAAGAGAAGTCCACGAGATGTCTAAGAAAGGCAGGCCCAAAGACAAAAGA
CGAGAAGTACATGAAGATGCCACAAGCAAGTCAGCCAAATCTGAGACGAAGTCTGACATCACAAA
TCGCCTCTGACAAAGTCAGAACAGCTTCTGAGGATAGATGACCATGATTTAGCATGAGGCGCTGGCTTT
GGAGGCCCTGCCATTCTGTTGGTGTGGATGTGCAGGTGGAGAGTTTGGATAGCATCTCAGAGGTTGAC
ATGGACTTTACGATGACCTCTACCTGAGGCACTACTGGAAGGACGAGAGGCTGTCTTTTCCAAGCACCC
AACAACTCAGCATGACGTTTGACGGCCGGCTGGTCAAGAAGATCTGGTCCCTGACATGTTTTTTCGTG
CACTCCAAACGCTCCTTCATCCACGACACCACCACAGACAACGTCATGTTGCGGGTCCAGCCTGATGGG
AAAGTGCTCTATAGTCTCAGGGTTACAGTAACTGCAATGTGCAACATGGACTTCAGCCGATTTCCCTTG
GACACACAACGCTGCTCTTTGAAATTGAAAGCTATGCCTATACAGAAGATGACCTCATGCTGTACTGG
AAAAAGGGCAATGACTCCTTAAAGACAGATGAACGGATCTCACTCTCCAGTTCTCATTAGGAATTC
CACACCACCACAACTGGCTTTCTACAGCAGCACAGGCTGGTACAACCGTCTCTACATTAATTTACAG
TTGCGTCGCCACATCTTCTTCTTCTGCTCCAACTTATTTCCCGCTACCCTGATGGTCATGCTGTCC
TGGGTGCTCTTCTGGATCGACCCGAGCCGTGCCTGCCAGAGTCCCCTTAGGTATCACAACGGTGCTG
ACCATGTCCACCATCATACGGGCGTGAATGCCTCCATGCCGCGGTCTCCTACATCAAGCCGTGGAC
ATCTACCTCTGGGTCAGCTTTGTGTTCTGTTCTCTCGGTGCTGGAGTATGCGGCCGTCAACTACCTG
ACCACTGTGCAGGAGAGGAAGGAACAGAAGCTGCGGGAGAAGCTTCCCTGCACCAGCGGATTACCTCCG
CCCCGCACTGCGATGCTGGACGGCAACTACAGTGTGGGGAGGTGAATGACCTGGACAACATACGCCA
GAGAATGGAGAGAAGCCCGACAGGATGATGGTGCAGCTGACCTGGCCTCAGAGAGGAGCTCCACAG
AGGAAAAGTCAGAGAAGCAGCTATGTGAGCATGAGAATCGACACCACGCCATTGATAAATACTCCAGG
ATCATCTTTCCAGCAGCATACATTTTATTCAATTAATACTGGTCTATTTTCTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



<b>Restriction Sites:</b>	Sgfl-Mlul
<b>Plasmid Map:</b>	□
<b>ACCN:</b>	NM_002042
<b>Insert Size:</b>	1440 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_002042.4</a>
<b>RefSeq Size:</b>	3174 bp
<b>RefSeq ORF:</b>	1440 bp
<b>Locus ID:</b>	2569
<b>UniProt ID:</b>	<a href="#">P24046</a>
<b>Cytogenetics:</b>	6q15
<b>Protein Families:</b>	Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane
<b>Protein Pathways:</b>	Neuroactive ligand-receptor interaction
<b>MW:</b>	55.9 kDa
<b>Gene Summary:</b>	<p>GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA receptors, which are ligand-gated chloride channels. GABRR1 is a member of the rho subunit family. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>