

## Product datasheet for **SC310355**

### UPF3B (NM\_023010) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	UPF3B (NM_023010) Human Untagged Clone
Tag:	Tag Free
Symbol:	UPF3B
Synonyms:	HUPF3B; MRX62; MRX82; MRXS14; RENT3B; UPF3BP1; UPF3BP2; UPF3BP3; Upf3p-X; UPF3X
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC310355 representing NM\_023010.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAGGAAGAGAAGGAGCACAGGCCTAAGGAGAAGCGAGTAACCTGTAAACCCCGCCGGGGCCACA
GGCAGCGGTGGTGGGACCTCGGGGACAGCTCCAAGGGGAAGATAAGCAGGATCGCAACAAGGAGAAG
AAAGAAGCGCTGAGCAAGGTGGTAATTCGAAGATTACCTCCCACTTTGACCAAGGAGCAGCTTCAGGAA
CATCTTCAACCTATGCCTGAGCATGATTATTTTGTAGTTTTTTCTAATGATACGAGTTTGTATCCTCAT
ATGTATGCCAGAGCATACATCAACTTTAAAAACCAAGAGGACATTATTTTGTTCAGGGATCGCTTTGAT
GGTTATGTATTCCTTGACAATAAAGGTCAGGAATATCCCGCTATAGTAGAATTTGCACCTTTTAAAAA
GCTGCAAAAAAGAAGACTAAGAAAAGAGATACCAAGTCGGGACTATCGATGATGATCCAGAATATAGA
AAGTTTTTGGAAAGTTATGCCACAGATAATGAGAAAATGACATCTACTCCAGAGACACTGCTAGAGGAA
ATAGAAGCAAAAAATAGAGAATTAATAGCTAAAAAGACAACCCCACTTTTGTAGCTTCTGAAAAACAAG
CAGAGAATGAGAGAAGAAAAGAGAGAAGAAAGGAGGAGGAGAGAAATAGAAAGAAAAAGACAAAGAGAA
GAAGAGAGGAGGAAATGGAAAGAAGAAGAGAAACGAAAAAGGAAAGATATAGAAAAGCTAAAGAAGATA
GACAGAATTCCAGAAAGGACAAATTAAGGATGAACCAAGATTAAGCTGCTCAAGAAGCCAGAAAAA
GGAGATGAAAAAGAATTGGACAAAAGAGAAAAAGCCAAGAAATTGGACAAAGAGAATCTCAGTGATGAA
AGAGCCAGTGGGCAAGTTGTACATTGCCCAAGCGTTCTGATAGCGAACTTAAAGATGAAAAACCAAG
AGACCTGAAGATGAGAGCGGCAGAGACTATAGGGAGAGGGAACGGGAATATGAACGAGATCAGGAGCGC
ATCTTCGAGAAAGAGAGAGGCTGAAGCGGCAAGAAGAAGAGCGCGTAGGCAGAGGAGCGCTATGAG
AAAGAGAAGACTTTTAAGAGAAAAGAAGAAGAAATGAAAAAGAGAAAGACACTTCGGGATAAAGGA
AAGAAGGCTGAAAGTACAGAATCAATAGGCAGCTCAGAAAAACTGAAAAGAAAGAAGTGGTCAAG
AGAGATCGAATAAGAAACAAGGATCGTCCAGCGATGCAGCTTTACCAACCAGGAGCTCGAAGCCGAAT
CGACTCTGTCCCCTGATGACAGCACCAAGTCTGGAGATTCAGCAGCAGAAAGGAAGCAGGAAAGTGGT
ATTAGCCATAGAAAAAGAAGGAGGAGAGGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_023010

**Insert Size:** 1413 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_023010.3</u>
<b>RefSeq Size:</b>	2365 bp
<b>RefSeq ORF:</b>	1413 bp
<b>Locus ID:</b>	65109
<b>UniProt ID:</b>	<u>Q9BZ17</u>
<b>Cytogenetics:</b>	Xq24
<b>Domains:</b>	Smg4_UPF3
<b>MW:</b>	56.2 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein that is part of a post-splicing multiprotein complex involved in both mRNA nuclear export and mRNA surveillance. The encoded protein is one of two functional homologs to yeast Upf3p. mRNA surveillance detects exported mRNAs with truncated open reading frames and initiates nonsense-mediated mRNA decay (NMD). When translation ends upstream from the last exon-exon junction, this triggers NMD to degrade mRNAs containing premature stop codons. This protein binds to the mRNA and remains bound after nuclear export, acting as a nucleocytoplasmic shuttling protein. It forms with Y14 a complex that binds specifically 20 nt upstream of exon-exon junctions. This gene is located on the long arm of chromosome X. Two splice variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) lacks an alternate coding segment, compared to variant 1, and encodes the shorter protein (isoform 2).</p>