

Product datasheet for SC310405

JNK1 (MAPK8) (NM_139047) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	JNK1 (MAPK8) (NM_139047) Human Untagged Clone
Tag:	Tag Free
Symbol:	JNK1
Synonyms:	JNK; JNK-46; JNK1; JNK1A2; JNK21B1/2; PRKM8; SAPK1; SAPK1c
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310405 representing NM_139047. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGAGCAGAAGCAAGCGTGACAACAATTTTATAGTGTAGAGATTGGAGATTCTACATTCACAGTCCTG
 AAACGATATCAGAATTTAAACCTATAGGCTCAGGAGCTCAAGGAATAGTATGCGCAGCTTATGATGCC
 ATTCTTGAAAGAAATGTTGCAATCAAGAAGCTAAGCCGACCATTTTCAAGTATCAGACTCATGCCAAGCGG
 GCCTACAGAGAGCTAGTTCTTATGAAATGTGTTAATCACAAAAATATAATTGGCCTTTTGAATGTTTTC
 ACACCACAGAAATCCCTAGAAGAATTTCAAGATGTTTACATAGTCATGGAGCTCATGGATGCAAAATCTT
 TGCCAAGTGATTGAGTGGAGCTAGATCATGAAAGAATGTCCTACCTTCTCTATCAGATGCTGTGTGGA
 ATCAAGCACCTTCATTCTGCTGGAATTATTCATCGGGACTTAAAGCCAGTAATATAGTAGTAAATCT
 GATTGCACTTTGAAGATCTTGACTTCGGTCTGGCCAGGACTGCAGGAACGAGTTTTATGATGACGCTT
 TATGTAGTGACTCGCTACTACAGAGCACCCGAGGTCATCCTTGGCATGGGCTACAAGGAAACGTTGAC
 ATTTGGTCAGTTGGGTGCATCATGGGAGAAATGATCAAAGGTGGTGTGTTTGTCCAGGTACAGATCAT
 ATTGATCAGTGAATAAAGTTATTGAACAGCTTGAACACCATGTCCTGAATTCATGAAGAACTGCAA
 CCAACAGTAAGGACTTACGTTGAAAACAGACCTAAATATGCTGGATATAGCTTTGAGAACTCTTCCCT
 GATGTCCTTTTCCAGCTGACTCAGAACACAACAACTTAAAGCCAGTCAGGCAAGGGATTTGTTATCC
 AAAATGCTGGTAATAGATGCATCTAAAAGGATCTCTGTAGATGAAGCTCTCCAACACCCGTACATCAAT
 GTCTGGTATGATCCTTCTGAAGCAGAAGCTCCACCACCAAGATCCCTGACAAGCAGTTAGATGAAAGG
 GAACACACAATAGAAGAGTGGAAGAATTGATATATAAGGAAGTTATGGACTTGGAGGAGAGAACCAAG
 AATGGAGTTATACGGGGCAGCCCTCTCCTTTAGGTGCAGCAGTGATCAATGGCTCTCAGCATCCATCA
 TCATCGTCGTCTGCAATGATGTGCTTCAATGTCAACAGATCCGACTTTGGCCTCTGATACAGACAGC
 AGTCTAGAAGCAGCAGCTGGGCCTCTGGGCTGCTGTAG**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_139047
Insert Size:	1284 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_139047.1</u>
RefSeq Size:	1412 bp
RefSeq ORF:	1284 bp
Locus ID:	5599
Cytogenetics:	10q11.22
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase
Protein Pathways:	Adipocytokine signaling pathway, Colorectal cancer, Epithelial cell signaling in Helicobacter pylori infection, ErbB signaling pathway, Fc epsilon RI signaling pathway, Focal adhesion, GnRH signaling pathway, Insulin signaling pathway, MAPK signaling pathway, Neurotrophin signaling pathway, NOD-like receptor signaling pathway, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway, Type II diabetes mellitus, Wnt signaling pathway
MW:	48.1 kDa

Gene Summary:

The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. This kinase is activated by various cell stimuli, and targets specific transcription factors, and thus mediates immediate-early gene expression in response to cell stimuli. The activation of this kinase by tumor-necrosis factor alpha (TNF-alpha) is found to be required for TNF-alpha induced apoptosis. This kinase is also involved in UV radiation induced apoptosis, which is thought to be related to cytochrom c-mediated cell death pathway. Studies of the mouse counterpart of this gene suggested that this kinase play a key role in T cell proliferation, apoptosis and differentiation. Several alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Apr 2016]

Transcript Variant: This variant (JNK1-b2) encodes the longer of the two JNK1 beta isoforms (JNK1 beta2). The JNK1-b2 variant differs from the JNK1-a2 variant in the use of an alternate internal coding exon of the same length. Thus, JNK1 beta2 isoform is the same length as JNK1 alpha2 isoform, with a few aa difference in an internal protein segment.